



**EVOLUTION<sup>®</sup>**  
Flexible, Modular, Column-Based System

WE BRING MORE TO THE TABLE

# As your needs evolve, the Evolution® System adapts. Beautifully.

With its sleek lines and durable construction, Kewaunee's Evolution® is the perfect example of form and function coming together in one modular workplace system.

Evolution is a streamlined, column-based system that was designed with your changing needs in mind. It offers infinite height adjustability for a wide range of components. And almost limitless reconfiguration possibilities. Think of Evolution as the starting point for a complete redesign of your work environment, without costly relocation of walls and other structures.

Ideal for healthcare, research and development, light manufacturing and assembly, high-tech, testing operations, and scientific and technical labs, the Evolution system is an example of form, function, and flexibility at its best. New technician on staff? Evolution adapts. OSHA/ADA requirements changing? Evolution can meet most any workplace challenge.



## What makes Evolution so special?

Evolution's unmatched range of features makes it the most flexible modular system available today:

- Workstations available in single and double-sided configurations
- Advanced laminates, cured maple, and Kemresin® worksurfaces are available in standard widths of 36" to 96"; and standard depths of 24", 30", and 36"
- Corner units create a variety of "U," "L," and "T" configurations; available in radius, square, and as "sit-at" corners
- Bridge units convert the opening between two free-standing workstations into a continuous worksurface
- Metal shelves can be mounted above or below the worksurface
- Fixed and mobile floor shelving are also available
- Lockable metal drawers feature smooth, ball-bearing glides
- Lockable steel overhead storage units provide 15" of height and depth
- Powered backsplashes are UL Listed and switched with built-in circuit breaker
- Dual-channel utility raceways promote orderly, segregated management of electrical wiring, air and gas lines, and cabling
- Optional overhead and undershelf task lights



# There's an Evolution<sup>®</sup> System for every environment.

## Standard Duty Evolution

- Rated for 1,000 lbs.
- Worksurface rated for 700 lbs.
- Shelving rated for 150 lbs., with optional 500 lbs.
- Power, data and services can be placed in columns
- Infinitely height adjustable
- Variety of optional casters available

## Wide Column Evolution

- Rated for 1,750 lbs.
- Worksurface rated for 1,000 lbs.
- Shelving rated for 150 lbs., with optional 500 lbs.
- Wide channel design for easy placement of significant power, data, and service elements in columns
- Infinitely height adjustable
- Variety of optional casters available



## Heavy Duty Evolution

- Rated at 2,500 lbs.
- Shelving rated for 150 lbs., with optional 500 lbs.
- Power, data and services can be placed in columns
- Incrementally height adjustable
- Variety of optional casters available

## Evolution ADJUSTAbench<sup>®</sup>

- Height-adjustable workstation for multi-shift service
- Adjust the worksurface at the push of a button to any height within the 28" to 44" design range
- Smoothly and quietly assumes a worker's precise comfort setting
- Can be programmed to "remember" up to three individual settings for quick, pre-set repositioning from shift to shift, or person to person
- Standard 36", 48", 60", or 72" widths
- Offers a wide range of accessories and options, including ESD protection

## The Evolution line-up



Standard Duty Evolution



Wide Column Evolution



Heavy Duty Evolution



ADJUSTAbench<sup>®</sup>

# Customize your workspace with Evolution<sup>®</sup> accessories.

To accommodate your requirement for flexibility, Kewaunee offers a complete line of lighting and computer accessories designed specifically to complement the Evolution - and the way you work.



**EVOLUTION<sup>®</sup>**  
Flexible, Modular, Column-Based System

## Lighting



LED Under-Shelf Task Light



Fluorescent Under-Shelf Task Light



Fluorescent Overhead Task Light

## Computer Accessories



LCD Monitor Arm/Worksurface Mounted



LCD Monitor Arm/Column Mounted



CPU Holder



Keyboard/Mouse Tray


# Making Evolution<sup>®</sup> work for you.

At Kewaunee, flexibility is a specialty. Our expert engineering services, advanced facilities, and computerized manufacturing processes deliver a product that effectively meets even the most specific need on a cost- and time-efficient basis.

Whether it's as simple as non-standard workstation dimensions or as complex as specialized laboratory and utility requirements, Kewaunee can meet your needs.




## Evolution Steel + Accent Color Palette



## Laminate Worksurface Palette



## Kemresin<sup>®</sup> Worksurface Palette



## Phenolic Worksurface Palette



## About KEWAUNEE Technical Furniture Group

Headquartered in Statesville, North Carolina, the Kewaunee Technical Furniture Group specializes in intelligent, highly flexible furniture solutions for high-tech, healthcare, industrial, institutional, and manufacturing environments.

The Technical Furniture Group offers a broad choice of systems, including Sturdilite® pedestal-based modular design, and Evolution® column-based modular design. The Kewaunee line also includes BASIKBench®, SturdiKwik®, and K4B stand-alone benches. In addition to this wide range of products, Kewaunee provides design and engineering services to ensure the best match of furniture to application.

The Technical Furniture Group is a division of Kewaunee Scientific Corporation, since 1906 a recognized leader in the design and manufacture of laboratory systems, technical and computer furniture, research and development workstations, industrial benches, enclosures, and more.



WE BRING MORE TO THE TABLE™

Technical Furniture Group  
2700 West Front Street  
Statesville, North Carolina 28677  
800-824-6626  
[MoreToTheTable.com](http://MoreToTheTable.com)