



Miscellaneous Accessories

Writing Boards & Tackboards




Black Chalkboard
F-5336-00 18¹/₄"H x 24¹/₄"L
F-5337-00 24¹/₄"H x 36¹/₄"L
F-5338-00 36¹/₄"H x 48¹/₄"L

Surface: washable Duracite
Frame Material: anodized aluminum
Trough Material: anodized aluminum




White Magnetic Markerboard
F-5342-00 18¹/₄"H x 24¹/₄"L
F-5343-00 24¹/₄"H x 36¹/₄"L
F-5344-00 36¹/₄"H x 48¹/₄"L

Surface: 1/4" white magnetic
Frame Material: grey mylar vinyl



Green Chalkboard
F-5333-00 18¹/₄"H x 24¹/₄"L
F-5340-00 24¹/₄"H x 36¹/₄"L
F-5341-00 36¹/₄"H x 48¹/₄"L


Surface: washable Duracite
Frame Material: anodized aluminum
Trough Material: anodized aluminum



Cork Tackboard
F-5418-24 18"H x 24"L
F-5424-36 24"H x 36"L
F-5436-48 36"H x 48"L

Surface: Resilient natural cork
Frame Material: anodized aluminum


Kemresin Balance Table



30" High x 24" Deep x 35" Long
K7-302435-BK Black
K7-302435-GR Grey
K7-302435-PT Putty
K7-302435-SL Slate

Work Top: 2" thick Kemresin
Work Top Size: 2"T x 24"D x 35"L
Pedestal Legs: 2" thick Kemresin
Total Weight: 300 pounds

Marvel Refrigerators & Freezers



L-110 Refrigerator/Freezer
Size: 34¹/₂" H x 24" D x 23⁷/₈" L
 6.1 cubic feet
Door Swing: Right Hand
 (User changeable with swing reversing kit)
Shelves: (3) Adjustable
Service: Single-phase,
 60 hertz, 120 volt AC

L-105 Refrigerator
 Same as above except
 without freezer compartment.


L-115 Freezer
 Same as above except
 entire unit is a freezer.

Gillis Stainless Steel Cart



F-5349-00
Overall Size: 33³/₈"H x 19"D x 31"L
Shelf Size: (3) 18"D x 27"L
 (1" high lip on 3 sides of
 top & middle shelf)
Capacity: 400 pounds
Casters: (4) Black rubber
 4" Diameter – Swivel

Lund Key Cases



F-5332-00
Height: 11"
Depth: 2¹/₂"
Length: 17⁵/₁₆"
Key Capacity: 90 Keys

F-5332-01
Height: 24³/₄"
Depth: 2¹/₂"
Length: 17⁵/₁₆"
Key Capacity: 200 Keys