

**BASIK**  
BENCH



K E W A U N E E

## VERSATILE, FUNCTIONAL, ADJUSTABLE



BASIKBench accessory uprights are designed to make height adjustments quick and easy and its above-worksurface adaptability enables you to configure your work environment for increased productivity.



Choose from a full range of accessories, including adjustable metal or wire shelves, powered back-splashes, task lighting, hanging drawers/cabinets and component handling systems. Several hanging drawers/cabinets are available including 4" high personal drawers, 7" high box drawers or 14" high file drawers and cabinets.

Engineered with our "C-Frame" design for strength, stability, and under-worksurface mobility, BASIKBench gives you the option to economically design your own workstation. BASIKBench is constructed from heavy-gauge rectangular steel and includes a 29" deep front radius edge laminate work surface or a square edge maple work surface. The BASIKBench is available in widths of 48", 60" or 72" and heights of 30" or 36".

The BASIKBench can be configured with a variety of components suitable for test and measurement labs, production/assembly lines, medical and healthcare applications, clean rooms, instrumentation labs, computer and CAD operations, and technical support areas.

ESD work surfaces are supplied with a complete grounding kit. This kit contains two grounding terminuses, one working as the single workstation grounding point and the second as a continuity test terminus. These terminuses are designed for quick and easy installation and will provide continuity with the static dissipative laminate and wrist straps with grounding blocks.

Call

**1-800-824-6626**

or Fax

**1-800-624-5536**

to contact your Kewaunee Sales Representative  
for more information about the  
BASIKBench by Kewaunee.



2700 West Front Street  
Statesville NC 28677  
1-800-824-6626

Email: [custsvc@kewaunee-tpg.com](mailto:custsvc@kewaunee-tpg.com)