


# 6" Partition Modules

## Features:

- 6" Partition Modules can be used as Island Modules, or as a Partition System
- 6" Partition Modules provide support for a wide combination of accessories, including:
  - adjustable worksurface frames
  - adjustable tables
  - hanging cabinetry — upper and lower
  - floor mounted cabinetry
  - adjustable shelving
  - lighting
  - electrical and plumbing raceway
  - service drops
  - tack boards and marker boards
- 6" Partition Modules provide 3½" of service space and are built to receive standard pipe hangers for ease of installing service piping and electrical conduit.
- System design allows shelving to be suspended both above and below the worksurface.
- Easy and positive structural alignment and attachment to adjacent modules.
- Available with easy to remove access panels. May be ordered without access panels for use with floor mounted base cabinets.

## SEFA 10 | Class 5 – Configuration 2 – Category 2



# 6" Partition Modules

## Specifications:


6" Partition Modules are 37" high, with Upper Carriers, 84" high. Both are 6" deep and available in lengths shown.

Modules are available with or without Access Panels. Order without when used with floor mounted base cabinets.

Each module is furnished with hardware for attaching adjacent modules.

**Ceiling Support and Module Ledge must be ordered separately.**


### 6" Deep Partition Module – 37" High



Length	with Access Panels	without Access Panels
24"	<b>F08M3724-00</b>	<b>F08-3724-00</b>
30"	<b>F08M3730-00</b>	<b>F08-3730-00</b>
36"	<b>F08M3736-00</b>	<b>F08-3736-00</b>
42"	<b>F08M3742-00</b>	<b>F08-3742-00</b>
48"	<b>F08M3748-00 *</b>	<b>F08-3748-00</b>
54"	<b>F08M3754-00 *</b>	<b>F08-3754-00</b>
60"	<b>F08M3760-00 *</b>	<b>F08-3760-00</b>
66"	<b>F08M3766-00 *</b>	<b>F08-3766-00</b>
72"	<b>F08M3772-00 *</b>	<b>F08-3772-00</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

### 6" Deep Partition Module with Upper Carrier – 84" High



Length	with Access Panels	without Access Panels
24"	<b>F38M3724-00</b>	<b>F38-3724-00</b>
30"	<b>F38M3730-00</b>	<b>F38-3730-00</b>
36"	<b>F38M3736-00</b>	<b>F38-3736-00</b>
42"	<b>F38M3742-00</b>	<b>F38-3742-00</b>
48"	<b>F38M3748-00 *</b>	<b>F38-3748-00</b>
54"	<b>F38M3754-00 *</b>	<b>F38-3754-00</b>
60"	<b>F38M3760-00 *</b>	<b>F38-3760-00</b>
66"	<b>F38M3766-00 *</b>	<b>F38-3766-00</b>
72"	<b>F38M3772-00 *</b>	<b>F38-3772-00</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

# 6" Partition Modules

**Specifications:**


Service Riser Partition Modules are 47" high, with Upper Carrier, 84" high. Both are 6" deep and in lengths shown.

Modules are available with or without Access Panels. Order without when used with floor mounted base cabinets.

Each module is furnished with hardware for attaching adjacent modules.

**Ceiling Support and Module Ledge must be ordered separately.**


**6" Deep Service Riser Partition Module – 47" High**



Length	with Access Panels	without Access Panels
24"	<b>F08M4724-00</b>	<b>F08-4724-00</b>
30"	<b>F08M4730-00</b>	<b>F08-4730-00</b>
36"	<b>F08M4736-00</b>	<b>F08-4736-00</b>
42"	<b>F08M4742-00</b>	<b>F08-4742-00</b>
48"	<b>F08M4748-00 *</b>	<b>F08-4748-00</b>
54"	<b>F08M4754-00 *</b>	<b>F08-4754-00</b>
60"	<b>F08M4760-00 *</b>	<b>F08-4760-00</b>
66"	<b>F08M4766-00 *</b>	<b>F08-4766-00</b>
72"	<b>F08M4772-00 *</b>	<b>F08-4772-00</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

**6" Deep Service Riser Partition Module with Upper Carrier – 84" High**



Length	with Access Panels	without Access Panels
24"	<b>F38M4724-00</b>	<b>F38-4724-00</b>
30"	<b>F38M4730-00</b>	<b>F38-4730-00</b>
36"	<b>F38M4736-00</b>	<b>F38-4736-00</b>
42"	<b>F38M4742-00</b>	<b>F38-4742-00</b>
48"	<b>F38M4748-00 *</b>	<b>F38-4748-00</b>
54"	<b>F38M4754-00 *</b>	<b>F38-4754-00</b>
60"	<b>F38M4760-00 *</b>	<b>F38-4760-00</b>
66"	<b>F38M4766-00 *</b>	<b>F38-4766-00</b>
72"	<b>F38M4772-00 *</b>	<b>F38-4772-00</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

# 6" Partition Modules

## Specifications:

Partition Modules are 37" high, with Upper Carriers, 84" high. Both are 6" deep and available in lengths shown.


Modules are available with or without Access Panels. Order without when used with floor mounted base cabinets.

24" and 30" long modules have two single-gang electric box cutouts, all other sizes have three.

Each module is furnished with hardware for attaching adjacent modules.

**Ceiling Support and Module Ledge must be ordered separately.**


## 6" Deep Standard Partition Module with Electric Cutouts – 37" High



Length	with Access Panels	without Access Panels
24"	<b>F08M3724-E0</b>	<b>F08-3724-E0</b>
30"	<b>F08M3730-E0</b>	<b>F08-3730-E0</b>
36"	<b>F08M3736-E0</b>	<b>F08-3736-E0</b>
42"	<b>F08M3742-E0</b>	<b>F08-3742-E0</b>
48"	<b>F08M3748-E0 *</b>	<b>F08-3748-E0</b>
54"	<b>F08M3754-E0 *</b>	<b>F08-3754-E0</b>
60"	<b>F08M3760-E0 *</b>	<b>F08-3760-E0</b>
66"	<b>F08M3766-E0 *</b>	<b>F08-3766-E0</b>
72"	<b>F08M3772-E0 *</b>	<b>F08-3772-E0</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

## 6" Deep Standard Partition Module with Electric Cutouts & Upper Carrier – 84" High



Length	with Access Panels	without Access Panels
24"	<b>F38M3724-E0</b>	<b>F38-3724-E0</b>
30"	<b>F38M3730-E0</b>	<b>F38-3730-E0</b>
36"	<b>F38M3736-E0</b>	<b>F38-3736-E0</b>
42"	<b>F38M3742-E0</b>	<b>F38-3742-E0</b>
48"	<b>F38M3748-E0 *</b>	<b>F38-3748-E0</b>
54"	<b>F38M3754-E0 *</b>	<b>F38-3754-E0</b>
60"	<b>F38M3760-E0 *</b>	<b>F38-3760-E0</b>
66"	<b>F38M3766-E0 *</b>	<b>F38-3766-E0</b>
72"	<b>F38M3772-E0 *</b>	<b>F38-3772-E0</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

# 6" Partition Modules

## Specifications:

Service Riser Partition Modules are 47" high, with Upper Carrier, 84" high. Both are 12" deep and in lengths shown.


Modules are available with or without Access Panels. Order without when used with floor mounted base cabinets.

24" and 30" long modules have two single-gang electric box cutouts, all other sizes have three.

Each module is furnished with hardware for attaching adjacent modules.

**Ceiling Support and Module Ledge must be ordered separately.**


### 6" Deep Service Riser Partition Module with Electric Cutouts – 47" High



Length	with Access Panels	without Access Panels
24"	<b>F08M4724-E0</b>	<b>F08-4724-E0</b>
30"	<b>F08M4730-E0</b>	<b>F08-4730-E0</b>
36"	<b>F08M4736-E0</b>	<b>F08-4736-E0</b>
42"	<b>F08M4742-E0</b>	<b>F08-4742-E0</b>
48"	<b>F08M4748-E0 *</b>	<b>F08-4748-E0</b>
54"	<b>F08M4754-E0 *</b>	<b>F08-4754-E0</b>
60"	<b>F08M4760-E0 *</b>	<b>F08-4760-E0</b>
66"	<b>F08M4766-E0 *</b>	<b>F08-4766-E0</b>
72"	<b>F08M4772-E0 *</b>	<b>F08-4772-E0</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

### 6" Deep Service Riser Partition Module with Electric Cutouts & Upper Carrier – 84" High



Length	with Access Panels	without Access Panels
24"	<b>F38M4724-E0</b>	<b>F38-4724-E0</b>
30"	<b>F38M4730-E0</b>	<b>F38-4730-E0</b>
36"	<b>F38M4736-E0</b>	<b>F38-4736-E0</b>
42"	<b>F38M4742-E0</b>	<b>F38-4742-E0</b>
48"	<b>F38M4748-E0 *</b>	<b>F38-4748-E0</b>
54"	<b>F38M4754-E0 *</b>	<b>F38-4754-E0</b>
60"	<b>F38M4760-E0 *</b>	<b>F38-4760-E0</b>
66"	<b>F38M4766-E0 *</b>	<b>F38-4766-E0</b>
72"	<b>F38M4772-E0 *</b>	<b>F38-4772-E0</b>

\* Modules 48" long and over are provided with two half-length Access Panels.

# 6" Partition Modules

## Specifications:


Partition Modules with Paneled Upper Carrier are the same as standard Partition Modules; except, they are furnished with easy to remove panels to close off the area between the upper carrier uprights. They may be used as a service drop, to hide piping or tubing, or as a room divider.

When ordered with electric cutouts, 24" and 30" long modules have two single-gang electric box cutouts, all other sizes have three.

Each module furnished with hardware for attaching adjacent modules.

**Ceiling Support must be ordered separately.  
Module Ledge is not required.**


### 6" Deep Standard Partition Module with Paneled Upper Carrier – 84" High



Length	with Lower Access Panels	without Lower Access Panels
24"	<b>F39M3724-00</b>	<b>F39-3724-00</b>
30"	<b>F39M3730-00</b>	<b>F39-3730-00</b>
36"	<b>F39M3736-00</b>	<b>F39-3736-00</b>
42"	<b>F39M3742-00</b>	<b>F39-3742-00</b>
48"	<b>F39M3748-00 *</b>	<b>F39-3748-00 *</b>
54"	<b>F39M3754-00 *</b>	<b>F39-3754-00 *</b>
60"	<b>F39M3760-00 *</b>	<b>F39-3760-00 *</b>
66"	<b>F39M3766-00 *</b>	<b>F39-3766-00 *</b>
72"	<b>F39M3772-00 *</b>	<b>F39-3772-00 *</b>

\* Modules 48" long and over are provided with half-length removable panels.

### 6" Deep Standard Partition Module with Electric Cutouts & Paneled Upper Carrier – 84" High



Length	with Lower Access Panels	without Lower Access Panels
24"	<b>F39M3724-E0</b>	<b>F39-3724-E0</b>
30"	<b>F39M3730-E0</b>	<b>F39-3730-E0</b>
36"	<b>F39M3736-E0</b>	<b>F39-3736-E0</b>
42"	<b>F39M3742-E0</b>	<b>F39-3742-E0</b>
48"	<b>F39M3748-E0 *</b>	<b>F39-3748-E0 *</b>
54"	<b>F39M3754-E0 *</b>	<b>F39-3754-E0 *</b>
60"	<b>F39M3760-E0 *</b>	<b>F39-3760-E0 *</b>
66"	<b>F39M3766-E0 *</b>	<b>F39-3766-E0 *</b>
72"	<b>F39M3772-E0 *</b>	<b>F39-3772-E0 *</b>

\* Modules 48" long and over are provided with half-length removable panels.

# 6" Partition Modules

## Specifications:


Service Riser Partition Modules with Paneled Upper Carrier are the same as standard Service Riser Partition Modules; except, they are furnished with easy to remove panels to close off the area between the upper carrier uprights. They may be used as a service drop, to hide piping or tubing, or as a room divider.

When ordered with electric cutouts, 24" and 30" long modules have two single-gang electric box cutouts, all other sizes have three.

Each module furnished with hardware for attaching to adjacent modules.

**Ceiling Support must be ordered separately.  
Module Ledge is not required.**


### 6" Deep Service Riser Partition Module with Paneled Upper Carrier – 84" High



Length	with Lower Access Panels	without Lower Access Panels
24"	<b>F39M4724-00</b>	<b>F39-4724-00</b>
30"	<b>F39M4730-00</b>	<b>F39-4730-00</b>
36"	<b>F39M4736-00</b>	<b>F39-4736-00</b>
42"	<b>F39M4742-00</b>	<b>F39-4742-00</b>
48"	<b>F39M4748-00 *</b>	<b>F39-4748-00 *</b>
54"	<b>F39M4754-00 *</b>	<b>F39-4754-00 *</b>
60"	<b>F39M4760-00 *</b>	<b>F39-4760-00 *</b>
66"	<b>F39M4766-00 *</b>	<b>F39-4766-00 *</b>
72"	<b>F39M4772-00 *</b>	<b>F39-4772-00 *</b>

\* Modules 48" long and over are provided with half-length removable panels.

### 6" Deep Service Riser Partition Module with Electric Cutouts & Paneled Upper Carrier – 84" High



Length	with Lower Access Panels	without Lower Access Panels
24"	<b>F39M4724-E0</b>	<b>F39-4724-E0</b>
30"	<b>F39M4730-E0</b>	<b>F39-4730-E0</b>
36"	<b>F39M4736-E0</b>	<b>F39-4736-E0</b>
42"	<b>F39M4742-E0</b>	<b>F39-4742-E0</b>
48"	<b>F39M4748-E0</b>	<b>F39-4748-E0</b>
54"	<b>F39M4754-E0 *</b>	<b>F39-4754-E0 *</b>
60"	<b>F39M4760-E0 *</b>	<b>F39-4760-E0 *</b>
66"	<b>F39M4766-E0 *</b>	<b>F39-4766-E0 *</b>
72"	<b>F39M4772-E0 *</b>	<b>F39-4772-E0 *</b>

\* Modules 48" long and over are provided with half-length removable panels.

# 6" Partition Modules


## Specifications:

Partition Module Ceiling Supports consist of two cold rolled steel tube assemblies that telescope into the Partition Module Upper Carrier. They are designed to be fastened to channels (provided by others) attached to either the drop ceiling grid or to the roof supports.

Partition Module Ceiling Supports are required when the module is not adjacent to a fixed wall, or supported by a wing wall or fixed furniture.


Ceiling Partitions are available to finish off the space between the Ceiling Supports above the Partition Module to form a continuous floor-to-ceiling divider. They are available in a number of heights and in either glass or with a solid steel panel. Solid partitions may be finished to match the Partition Module or in a contrasting color.

### 6" Deep Partition Module Ceiling Support




Length	24" High	36" High	48" High	60" High
24"	<b>F48-2424-00</b>	<b>F48-3624-00</b>	<b>F48-4824-00</b>	<b>F48-6024-00</b>
30"	<b>F48-2430-00</b>	<b>F48-3630-00</b>	<b>F48-4830-00</b>	<b>F48-6030-00</b>
36"	<b>F48-2436-00</b>	<b>F48-3636-00</b>	<b>F48-4836-00</b>	<b>F48-6036-00</b>
42"	<b>F48-2442-00</b>	<b>F48-3642-00</b>	<b>F48-4842-00</b>	<b>F48-6042-00</b>
48"	<b>F48-2448-00</b>	<b>F48-3648-00</b>	<b>F48-4848-00</b>	<b>F48-6048-00</b>
54"	<b>F48-2454-00</b>	<b>F48-3654-00</b>	<b>F48-4854-00</b>	<b>F48-6054-00</b>
60"	<b>F48-2460-00</b>	<b>F48-3660-00</b>	<b>F48-4860-00</b>	<b>F48-6060-00</b>
66"	<b>F48-2466-00</b>	<b>F48-3666-00</b>	<b>F48-4866-00</b>	<b>F48-6066-00</b>
72"	<b>F48-2472-00</b>	<b>F48-3672-00</b>	<b>F48-4872-00</b>	<b>F48-6072-00</b>

### 6" Deep Partition Module Glass Ceiling Partition



Length	12" High	18" High	24" High	30" High
24"	<b>F49G1224-00</b>	<b>F49G1824-00</b>	<b>F49G2424-00</b>	<b>F49G3024-00</b>
30"	<b>F49G1230-00</b>	<b>F49G1830-00</b>	<b>F49G2430-00</b>	<b>F49G3030-00</b>
36"	<b>F49G1236-00</b>	<b>F49G1836-00</b>	<b>F49G2436-00</b>	<b>F49G3036-00</b>
42"	<b>F49G1242-00</b>	<b>F49G1842-00</b>	<b>F49G2442-00</b>	<b>F49G3042-00</b>
48"	<b>F49G1248-00</b>	<b>F49G1848-00</b>	<b>F49G2448-00</b>	<b>F49G3048-00</b>
54"	<b>F49G1254-00</b>	<b>F49G1854-00</b>	<b>F49G2454-00</b>	<b>F49G3054-00</b>
60"	<b>F49G1260-00</b>	<b>F49G1860-00</b>	<b>F49G2460-00</b>	<b>F49G3060-00</b>
66"	<b>F49G1266-00</b>	<b>F49G1866-00</b>	<b>F49G2466-00</b>	<b>F49G3066-00</b>
72"	<b>F49G1272-00</b>	<b>F49G1872-00</b>	<b>F49G2472-00</b>	<b>F49G3072-00</b>

### 6" Deep Partition Module Solid Ceiling Partition



Length	12" High	18" High	24" High	30" High
24"	<b>F49M1224-00</b>	<b>F49M1824-00</b>	<b>F49M2424-00</b>	<b>F49M3024-00</b>
30"	<b>F49M1230-00</b>	<b>F49M1830-00</b>	<b>F49M2430-00</b>	<b>F49M3030-00</b>
36"	<b>F49M1236-00</b>	<b>F49M1836-00</b>	<b>F49M2436-00</b>	<b>F49M3036-00</b>
42"	<b>F49M1242-00</b>	<b>F49M1842-00</b>	<b>F49M2442-00</b>	<b>F49M3042-00</b>
48"	<b>F49M1248-00</b>	<b>F49M1848-00</b>	<b>F49M2448-00</b>	<b>F49M3048-00</b>
54"	<b>F49M1254-00</b>	<b>F49M1854-00</b>	<b>F49M2454-00</b>	<b>F49M3054-00</b>
60"	<b>F49M1260-00</b>	<b>F49M1860-00</b>	<b>F49M2460-00</b>	<b>F49M3060-00</b>
66"	<b>F49M1266-00</b>	<b>F49M1866-00</b>	<b>F49M2466-00</b>	<b>F49M3066-00</b>
72"	<b>F49M1272-00</b>	<b>F49M1872-00</b>	<b>F49M2472-00</b>	<b>F49M3072-00</b>




# 6" Corner Column Modules

**Specifications:**


Corner Column Modules allow multiple configurations of Partition Modules to form a variety of assemblies. The Column Modules are available in three heights to match the Partition Modules. Each module is furnished with removable access panels for two sides and attachment hardware.

**Module Ledge, when required must be ordered separately.**

**6" Deep Partition Module Corner Column**



Height	Part Number
37"	<b>F40M3706-00</b>
47"	<b>F40M4706-00</b>
84"	<b>F40M8406-00</b>



**Possible Alpha System Assembly configurations using 6" Corner Column Modules**


# 6" Partition Module End Filler

**Specifications:**


Partition Module End Fillers are designed to close off the end of Partition Modules to offer a more finished appearance. A single filler is used to finish Standard Partition Modules. Two fillers are used to finish Partition Modules with Paneled Upper Carriers.

Partition Module End Fillers are provided with all required installation hardware.

**6" Deep Partition Module End Fillers**



Height	Lower Filler
35"	<b>F00M0636-00</b>
45"	<b>F00M0646-00</b>



Height	Upper Filler
38"	<b>F00M0638-00</b>
48"	<b>F00M0648-00</b>