

Free-Standing Adjustable Height Table


Specifications:

Free-Standing Adjustable Height Tables are available for either 24" deep or 30" deep worksurface tops in lengths as shown. Telescoping tubular steel legs allow the table to be height adjustable from 30" to 37" high on one inch increments. The table is furnished with a cold rolled steel frame and four adjustable legs with glides. Adjustable Height Tables will accommodate Suspended Cabinets.

Worksurface Top must be ordered separately.

Worksurface Tops are 1" thick. Front and side edges are treated as shown on page 69. Rear edge is finished square.

Free-Standing Adjustable Height Table Frame



22" Deep Table Frame

Length	Part Number
34"	FC50M362234-00
40"	FC50M362240-00
46"	FC50M362246-00
52"	FC50M362252-00
58"	FC50M362258-00
64"	FC50M362264-00
70"	FC50M362270-00

28" Deep Table Frame


Length	Part Number
34"	FC50M362834-00
40"	FC50M362840-00
46"	FC50M362846-00
52"	FC50M362852-00
58"	FC50M362858-00
64"	FC50M362864-00
70"	FC50M362870-00

Optional Rear Hanging Rail

Length	Part Number
34"	FC89M020234-00
40"	FC89M020240-00
46"	FC89M020246-00
52"	FC89M020252-00
58"	FC89M020258-00
64"	FC89M020264-00
70"	FC89M020270-00

Required for 28" deep Table Frames only when mounting suspended cabinets.

Free-Standing Adjustable Height Table Top



24" Deep Worksurface Top

Length	Kemresin®	Phenolic Resin	High Pressure Laminate
36"	FC96R012436-00	FC96T012436-00	FC96P012436-00
42"	FC96R012442-00	FC96T012442-00	FC96P012442-00
48"	FC96R012448-00	FC96T012448-00	FC96P012448-00
54"	FC96R012454-00	FC96T012454-00	FC96P012454-00
60"	FC96R012460-00	FC96T012460-00	FC96P012460-00
66"	FC96R012466-00	FC96T012466-00	FC96P012466-00
72"	FC96R012472-00	FC96T012472-00	FC96P012472-00

30" Deep Worksurface Top

Length	Kemresin®	Phenolic Resin	High Pressure Laminate
36"	FC96R013036-00	FC96T013036-00	FC96P013036-00
42"	FC96R013042-00	FC96T013042-00	FC96P013042-00
48"	FC96R013048-00	FC96T013048-00	FC96P013048-00
54"	FC96R013054-00	FC96T013054-00	FC96P013054-00
60"	FC96R013060-00	FC96T013060-00	FC96P013060-00
66"	FC96R013066-00	FC96T013066-00	FC96P013066-00
72"	FC96R013072-00	FC96T013072-00	FC96P013072-00