
Research Collection

Steel Laboratory Furniture











Contents

Steel Furniture Styles & Construction	2
Catalog Number Explanation	4
Steel Cabinet Options	6
Sitting Height Cabinets	. 11
ADA Height Cabinets	. 14
Standing Height Cabinets	. 22
Mobile Cabinets	.34
Apron Assemblies	.38
Wall Mounted Storage Cabinets	46
Full Height Storage Cabinets	. 53
Fillers, Scribes & Panels	. 58
Sloping Tops & Enclosures	. 66
Fume Hoods & Ventilated Enclosures	. 68
Work Tops	. 70
Center Bench Shelving	. 72
Wall Shelving	. 74
Reagent Racks & Shelving	. 76
Pegboards	. 78
Miscellaneous Accessories	. 79
General Plumbing Information	82
Service Fittings	. 83
Drain Accessories	. 91
Sinks & Troughs	. 92
Electrical Fittings	. 96
Seating	. 98
Technical Information	100
Alphabetical Index	106
Numerical Index	108
Typical Installations	112

Beautifully finished, smooth, clean, easy to maintain surfaces. Rugged, strong, heavy gauge steel, full-frame, four sided, flush front construction. That's Kewaunee's Steel Laboratory Furniture.



Kewaunee Steel Cabinet Construction

Kewaunee's fully-welded, heavy-gauge, steel cabinets are designed to SEFA 8 performance standards and built to rigorous benchmarks. The Kewaunee classic "unit" design allows doors

and drawers to be reconfigured in place. The configurable design provides a rich offering of hardware and aesthetic options.

DRAWER OPTIONS




Standard Drawers feature one-piece drawer body construction with full-extension, side mounted, soft-close, slides. (model F-8400-00)




Optional Concealed Slide Drawers feature integral concealed, full-extension, soft-close slides





- Four sided, flush front, frame construction
- Solid finished back with optional removable back
- Full-depth adjustable shelf
- Leveling bolts accessible through cabinet top



DOOR CATCH



Positive Catch
Catch furnished with 5-knuckle hinges only

DOOR HINGE OPTIONS



Inset 5-Knuckle



Overlay Semi-concealed



Overlay 5-Knuckle



Overlay Concealed 170° Swing



Kewaunee Steel Cabinet Styles



Inset - Square Edge Doors and drawer heads are flush with face of

Doors and drawer heads are flush with face of cabinet. Available with Positive Catch and a wide selection of pulls.



Overlay - Square Edge Doors and drawer heads overlay face of cabinet.

Doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.



Overlay - Wood on Steel

Oak or maple veneer doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls. (Flush Pull 9 not available on wood)




Overlay - Bevel Edge Doors and drawer heads, with beveled top edge,

Doors and drawer heads, with beveled top edge, overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.




Steel Cabinet Information

Catalog Number Explanation




Door & Drawer Head Style - Hardware - Material Selections



^{*} Stainless Steel Cabinets available in D&D Style 1 only

Steel Cabinet Information

Drawer and Cupboard Size Coding



E40C352215L	15"L
E40C352218L	18"L
E40C352221L	21"L
E40C352224L	24"L

Base cabinets are designed on a "unit" system to allow maximum flexibility in the selection of drawer and/or cupboard door sizes.

Standing height cabinets (35" high) have an opening 60 "units" high. Sitting height cabinets (29" high) have an opening 48 "units" high. ADA height cabinets (32" high) have an opening 54 "units" high. This means that any combination of drawer and cupboard door heights equaling 60, 54 or 48 units can be used in a base cabinet.

Mobile standing height cabinets (33" high) have an opening 56 "units" high. Mobile sitting height cabinets (261/4" high) have an opening 42 "units" high.

The adjacent tables indicate the "unit" height and inside clearance of the Inset style base cabinets shown in this catalog. Simply match the unit number shown on the cabinet drawers and cupboard to the appropriate number in the tables.

	Drawer	Inside Dimensions		
"Units" High	Head Height	Clear/ Height	Length (Front-to-Back)	
8*	3¾"	23/8"	18"	
9*	41/4"	23/8"	18"	
10*	43⁄4"	3 ⁷ /16"	18"	
11	5 ³ ⁄16"	3 ¹⁵ /16"	18"	
12	5 ¹¹ /16"	4 ⁷ /16"	18"	
14	65/8"	53/8"	18"	
15	71/8"	5 ⁷ /8"	18"	
16	75⁄8"	63/8"	18"	
20	9 %16"	8 ⁵ ⁄16"	18"	
24	11 ¹ /2"	10 ¹ / ₄ "	18"	
3*	1 3⁄16"	† Pull	Out Shelf	

^{*} Not available with Option Z Concealed Slide Drawer † Pull Out Shelf mounted on full extension slides only

Drawer Size Codes

Cupboard Size Codes

		Inside Dimensions		
"Units" High	Door Height	Height (Top-to-Bottom)	Depth (Front-to-Back)	
30	14 ⁹ /16"	13 ⁹ /16"	201/2"	
36	17 ⁷ /16"	16 ⁷ /16"	20½"	
45	21 ¹³ / ₁₆ "	20 ¹³ / ₁₆ "	201/2"	
48	23 ³ /16"	22 ³ /16"	20½"	
54	26 ¹ / ₁₆ "	25 ¹ /16"	20½"	
60	29"	28"	20½"	

Drawer Width I.D. (side-to-side)

Base	Number of Drawers Across			
Cabinet Length	1	2	3	
12"	83/8"			
15"	11 ³ ⁄8"			
18"	14 ³ ⁄8"			
21"	17 ³ ⁄8"			
24"	203/8"	83/8"		
30"	26 ³ / ₈ "	11 ³ ⁄8"		
36"	32 ³ /8"	14 ³ ⁄8"	83/8"	
42"	38 ³ ⁄8"	17 ³ ⁄8"		
48"	443/8"	203/8"	12 ³ ⁄8"	

Door Opening Clearance Width I.D. (side-to-side)

Base	Number of Drawers Across				Number of Drawers Acros		
Cabinet Length	1	2	3				
12"	91/2"						
15"	12½"						
18"	15½"						
21"	18½"						
24"	21½"	91/2"					
30"	27½"	12½"					
36"	33½"	15½"	91⁄2"				
42"	39½"	18½"					
48"	45½"	21½"	13½"				

Cupboard and Drawer sizes listed in these charts are for Inset style cabinets only. Sizes for Overlay styles are generally slightly larger. Contact your Kewaunee Dealer for more information.

Cabinet Options - (Add to End of Catalog Number)

B – Base Cabinet Caster Option

B11 – 4" Swivel Caster w/ brakes & with counter balance

B21 – 4" Swivel Caster w/o brakes & with counter balance

B31 - 4" Swivel Caster w/ brakes at front & counter balance

B41 – 3" Swivel Caster w/ brakes & with counter balance

B51 - 3" Swivel Caster w/o brakes & with counter balance

B61 – 3" Swivel Caster w/ brakes at front & counter balance

B71 - Fixed rear - Swivel front - with counter balance

B81 – Fixed rear – Glides front – with counter balance

C – Toe space depth increased by ⁷/8" (for B00 & G07)

F – Removable Back (not available on Mobile Cabinets)

H - Label Holders '

*J - Locks w/ Lock Panels J=Disk Locks† J2=Pin Locks†

*K - Locks w/o Lock Panels K=Disk Locks† K2=Pin Locks†

O - Shelf Lips

- P Pullout Trays & Pullout Pans
- P1 One Pullout Tray at bottom of cupboard
- P2 Two Pullout Trays center & bottom
- P3 One Pullout Pan at bottom of cupboard
- P4 Two Pullout Pans center & bottom
- I 1/4" Laminated Safety Glass
- R 1/8" Tempered Safety Glass T - Heavy Duty Full Extension Slide
- U Toe Space on Base Cabinet Ends
- U1 Toe Space on Left End
- U2 Toe Space on Right End
- U3 Toe Space on Both Ends
- Z Concealed Slide Drawer
- Z1 Cold Rolled Steel Drawer Sides
- **Z2** Stainless Steel Drawer Sides

Furnished on all cupboards and drawers on cabinet + Lock keying must be specified.

Option B



Steel Cabinet Options

Base Cabinet Caster Option

Base Cabinet Caster Option adds casters to the base cabinet. Eight caster styles are available and all are furnished with counter balance weight for cabinets that may be top heavy when loaded and/or prone to tipping when opened. (**Note:** Casters ship loose)

Swivel Casters

adds four (4) swivel casters and counter balance weight for stability.

B11 – *4" Swivel Casters with brakes

B21 - *4" Swivel Casters without brakes

B31 - *4" Swivel Casters with brakes at front only

B41 – 3" Swivel Casters with brakes **B51** – 3" Swivel Casters without brakes

B61 – 3" Swivel Casters with brakes at front only

* Adds 1" to height of cabinet

Fixed and Swivel Casters

adds two (2) fixed 4" casters at the rear, two locking swivel casters at the front, and counter balance weight for stability.

B71 – Fixed rear casters, locking front swivel casters

Fixed Casters and Glides

adds two (2) fixed 4" casters at the rear, two (2) adjustable glides at the front, and counter balance weight for stability


B81 - Fixed rear casters, front adjustable glides.

To specify Casters, add **B**__ to the end of the cabinet number.

Example: **E40C352236-0101-B71** (Available on base cabinets only.)

B11 – B21 – B31 Caster Option





B81 Caster Option

Toe Space Depth Increase by 7/8"

Option **C** adds 7/8" to the depth of a base cabinet toe space. This option is only used on base cabinets without doors or drawers (B00 & G07) to

allow the face of the cabinet to be pulled forward to align with the face of the doors and drawers on adjacent cabinets.

To specify increased toe space depth add **C** to the end of the

Removable Back

cabinet number.

Example: **B00C352236-0000-C**

Option C

Option F

All cabinets are furnished with a solid back for improved cleanliness. A removable back is available when required for access to the pipe space behind the cabinet.

Note: Removable Backs are not availabe on Mobile Cabinets.

To specify Removable Back, add **F** to the end of the cabinet number.

Example: **D30C352236-0101-F**



Label Holders Option H

Option **H** adds the F-0260-00 aluminum Label Holder to each drawer and door on the cabinet. Overall size is 3" x 1⁵/₈" and holds a 2¹/₂" x 1¹/₈" card. Double doors

recieve only one Label Holder per

Shipped loose for field installation.

To specify Label Holders, add **H** to the end of the cabinet number.

Example: **E40C352236-0101-H** (Added to all cupboards and drawers.)



One-half Actual Size

Research Collection

Steel Cabinet Options

Shelf Lips Option 0

Option **O** adds a 3/4" high steel lip to the front edge of all steel shelves in the cabinet.

To specify Shelf Lips, add **O** to the end of the cabinet number.

Example: **W00C361336-0101-O** (Added to all shelves in cabinet.)

Pullout Trays & Pans

- P1 One Pullout Tray at bottom of cupboard
- P2 Two Pullout Trays center & bottom
- P3 One Pullout Pan at bottom of cupboard
- P4 Two Pullout Pans center & bottom

Pullout Trays and Pans are fabricated from 18 gauge cold rolled steel and are mounted on full-extension slides. Both have 23/8" high sides and back. Trays are open at the front while pans are enclosed on all four sides. Each is available in two configurations: one tray or pan mounted at bottom of cupboard; or

two trays or pans mounted at bottom and in center of cupboard. When pullout trays or pans are added, the adjustable shelf is omitted.


Pullout Trays and Pans are available for base cabinets only.

To specify a Pullout Tray, add **P1** or **P2** to the end of the cabinet number.

Example: **E40C352236-0101-P2** (Available for base cabinets only.)

To specify a Pullout Pan, add **P3** or **P4** to the end of the cabinet number.

Example: **E40C352236-0101-P3** (Available for base cabinets only.)



Option P

Option I

Option R

Option T

Pullout Pan

Laminated Safety Glass Tempered Safety Glass

Option I replaces the standard float glass in plate glass or glazed door wall, counter mounted and full height cabinets with 1/4" Laminated Safety Glass.

To specify Laminated Saftey Glass, add I to the end of the cabinet number. Example: W20C351336-0101-I

Option **R** replaces the standard 1/8" float glass in glazed door wall, counter mounted and full height cabinets with 1/8" Tempered Safety Glass.

To specify Tempered Glass, add **R** to the end of the cabinet number. Example: **W20C351336-0101-R**

Full Extension Heavy Duty Slides

Full extension ball bearing drawer slide with cushioned in and outstops and positive drawer hold-in detent. Loading capacity of 200 lbs. Replaces standard 100 lbs full-

extension slides on all drawers in cabinet when **T** option is specified.

To specify Heavy Duty Slides, add **T** to the end of the cabinet number.

Example: **E40C352236-1011-T**



Toe Space on Cabinet Ends

When cabinets are used at the end of center benches and in some other special applications, they may require a toe space on one or

both ends of the cabinets. These toe spaces match the toe space at the front of the cabinets and my be added to the left, right, or both ends.

Example: **E40C352236-0101-U2** (Available for base cabinets only.)

Option U

U1 – Toe Space on Left EndU2 – Toe Space on Right End

U3 – Toe Space on Both Ends

To specify End Toe Spaces, add **U**_ to the end of the cabinet number.

SF-05/17-**7**



Steel Cabinet Options

Locks with Lock Panels Locks without Lock Panels

Option J Option K

All cabinets may be ordered with locks by specifying the K option. In addition, base cabinets may be specified with the J option to add lock panels. Lock panels are

installed between drawers and/or drawers and cupboards to provide additional security when the locks on a cabinet are keyed different.

Locks Types

Cabinet locks are available in two types. Disk Tumbler which offers medium security, and Pin Tumbler which offer a higher level of security. Pin Tumbler also allows more keying arrangements. Both are fabricated with die cast sleeves and cylinders with a satin nickel face and zinc plated steel cams. Both locks use the same sleeve and cam assembly allowing the lock cylinder to be removed and switched after the lock has been installed. All locks are furnished with two keys unless otherwise specified.



Lock Panel

Disk Tumbler Locks

Security: Medium

Keying: Different, Alike, Alike in

Sets, Masterkeyed

Key Changes: Up to 2000 Master Key Changes: Up to 5

J - Diek Locke with Lock Pa

J – Disk Locks with Lock PanelsK – Disk Locks w/o Lock Panels

Pin Tumbler Locks

Security: Medium-High

Keying: Different, Alike, Alike in Sets, Masterkeyed, Grandmasterkeyed, Sub-masterkeyed, Mason Keyed

Key Changes: Up to 2000 Master Key Changes: Up to 52

J2 – Pin Locks with Lock Panels K2 – Pin Locks w/o Lock Panels

Lock schedule with keying arrangement must be furnished when locks are specified. See page 9 & 10 for keying arrangement diagrams, and page 104 & 105 for sample lock schedules.

To specify Locks without Lock Panels, add **K** or **K2** to the end of the cabinet number.

Example: **E40C352236-0101K2** (Added to all cupboards and drawers.)

To specify Locks with Lock Panels, add **J** or **J2** to the end of the cabinet number

Example: **E40C352236-0101J** (Added to all cupboards and drawers.)



Disk Tumbler Lock (Medium Security)



Pin Tumbler Lock (Higher Security)

Concealed Slide Drawers

Overlay cabinets may be

steel drawer sides or with stainless steel drawer sides.

Concealed Slide Drawers are

upgraded to a Concealed Slide
Drawer with integral, concealed,
full-extension,soft-close, slides.
Concealed Slide Drawers are
available with either powder-coated

To specify Concealed Slide Drawers with Powder-coated Steel sides, add **Z1** to the end of the cabinet number.

Example: **E40C352236-1101Z1** (Added to all drawers.)

To specify Concealed Slide Drawers with Stainless Steel Sides, add **Z2** to the end of the cabinet number.

available in Overlay cabinets only and

are not available for all drawer sizes.

Example: **E40C352236-1101Z2** (Added to all drawers.)

Option Z

Research Collection

Lock Keying Options


Keyed Different

Disk Tumbler Locks: Available **Pin Tumbler Locks**: Available

All locks are keyed different from each other with individual keys for each lock. No key will open more than one lock.

Individual Locks

Individual Keys



Keyed Alike

Disk Tumbler Locks: Available **Pin Tumbler Locks**: Available

All locks are keyed alike. One key will open all locks.


Individual Locks

Single Key

Keyed Alike in Sets

Disk Tumbler Locks: Available Pin Tumbler Locks: Available


Locks are grouped into sets, with as many locks as needed in each set. All locks in a set can be opened by a single key, but that key can open no other locks. No key will open all locks.



Master Keyed

Disk Tumbler Locks: Available **Pin Tumbler Locks**: Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Muptiple sets are available, but no key can open all locks in all sets.






Lock Keying Options


Grand Master Keyed

Disk Tumbler Locks: N/A
Pin Tumbler Locks: Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, with Grand master key that can open all locks in all sets.



Sub-master Keyed & Mason Keyed



Sub-master Keyed

All locks are keyed different from each other with individual keys for each lock. Sub-master key provided that can open all locks in sub-set. Multiple sub-sets are available, with master key that can open all locks in all sub-sets. Multiple sets are available, with grand master key that can open all locks in all sets and sub-sets.

Mason Keyed

Locks are keyed different from each other with individual keys for each lock. Each individual key can also open all common locks in set. Master key provided that can open all locks in set. Multiple sets are available with grand master key that can open all locks in all sets.



Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Open and Cupboard Cabinets are furnished with one adjustable shelf.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Open Cabinets



B00C292212-00 0 12"L B00C292215-00_0 15"L B00C292218-00_0 18"L 21"L B00C292221-00_0 B00C292224-00 0 24"L


B00C292230-00 0 30"L B00C292236-00_0 36"L B00C292242-00_0 42"L B00C292248-00 0 48"L Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a C to the end of the part number.

Example: **B00C292236-0000-C**

Cupboard Cabinets







C03C292215L 15"L C03C292218L 18"L C03C292221L 21"L C03C292224L 24"L C03C292215-15"L C03C292218-___ 18"L C03C292221-21"L 24"L C03C292224C03C292230-30"L C03C292236-36"L C03C292242-42"L C03C292248-48"L

Drawer Cabinets







D00C292215-15"L 18"L D00C292218-D00C292221-21"L D00C292224-24"L

D00C292230-30"L D00C292236-36"L D00C292242-42"L D00C292248-48"L

D26C292218-18"L 1 Pull Out Shelf

12u





— _{12u}	— _{12u}
— _{12u}	— _{12u}
— _{12u}	—
— _{12u}	— _{12u}


D30C292215-15"L D30C292218-18"L 21"L D30C292221-D30C292224-24"L

30"L D30C292230-D30C292236-36"L D30C292242-42"L D30C292248-48"L

D31C292230-30"L D31C292236-36"L D31C292242-42"L D31C292248-48"L

D32C292230-30"L D32C292236-36"L D32C292242-42"L D32C292248-48"L





Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Drawer Cabinets



D81C292218-____ 18"L **D81C292221-**____ 21"L **D81C292224-**___ 24"L

2 File Drawer with hangers




D84C292218-____ 18"L **D84C292221-**____ 21"L **D84C292224-**___ 24"L


Pull Out Shelf, 1 File Drawer with hangers



1 File Drawer with hangers

Cupboard and Drawer Cabinets






E05C292230L 30"L E05C292230-30"L E05C292236-E05C292236L 36"L 36"L E05C292242-42"L E05C292242L 42"L E05C292248L 48"L E05C292248-48"L







E40C292215L_ E40C292215-___ 15"L 15"L E41C292230-30"L E40C292218L____ E40C292218-___ 18"L E41C292236-18"L 36"L E40C292221-___ E41C292242-21"L E40C292221L 21"L 42"L 24"L E40C292224-24"L E41C292248-E40C292224L 48"L




Specifications:


Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Cupboard and Drawer Cabinets





E45C292230L	30"L	E45C292230	30"L
E45C292236L	36"L	E45C292236	36"l
E45C292242L	42"L	E45C292242	42"l
E45C292248L	48"L	E45C292248	48"l

Sink Cabinets








I I	ĺ
	36u

G00C292215L	15"L	G00C292215	15"L	G00C292230	30"L
G00C292218L	18"L	G00C292218	18"L	G00C292236	36"L
G00C292221L	21"L	G00C292221	21"L	G00C292242	42"L
G00C292224L	24"L	G00C292224	24"L	G00C292248	48"L

Diagonal Corner Cabinets







Sink Support Kits



Part Number	Fits
K-MSS-15	15" Cabinet
K-MSS-18	18" Cabinet
K-MSS-24	24" Cabinet
K-MSS-30	30" Cabinet
K-MSS-36	36" Cabinet
K-MSS-42	42" Cabinet
K-MSS-48	48" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support undermount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: K-MSS-36-D





Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.


Open and Cupboard Cabinets are furnished with one adjustable shelf.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Open Cabinets



B00C322212-00_0 12"L B00C322215-00_0 15"L B00C322218-00_0 18"L B00C322221-00_0 21"L B00C322224-00_0 24"L



 B00C322230-00_0
 30"L

 B00C322236-00_0
 36"L

 B00C322242-00_0
 42"L

 B00C322248-00_0
 48"L


Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.

Example: G00C322236-0000-C

Cupboard Cabinets



C03C322215L 15"L C03C322218L 18"L C03C322221L 21"L C03C322224L 24"L




15"L

18"L

21"L

24"L

C03C322215-____ C03C322218-___ C03C322221-___ C03C322224-___



C03C322230-_____ 30"L C03C322236-____ 36"L C03C322242-____ 42"L C03C322248- 48"L



Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


Drawer Cabinets



D00C322215	15"L
D00C322218	18"L
D00C322221	21"L
D00C322224	24"L




D00C322230	30"L
D00C322236-	36"L
D00C322242	42"L
D00C322248	48"L



D02C322230	30"L
D02C322236-	36"L
D02C322242	42"L
D02C322248-	48"I





_	— 9u
_	— 15u
_	— 15u
_	—

D30C322230	30"L
D30C322236	36"L
D30C322242	42"L
D30C322248	48"L

— 9u	— 9u
_	— _{15u}
_	— _{15u}
_	— 15u

D31C322236- 42 D31C322248- 48
--

— 9u	— 9u
— 15u	—
	— 15u
— 15u	— 15u

D32C322230-	30"L
D32C322236	36"L
D32C322242	42"L
D32C322248-	48"L

— 9u	— 9u
—	—
_	— 15u
_	— 15u
	100

D36C322230	30"L
D36C322236	36"L
D36C322242	42"L
D36C322248	48"L

_	9u
_	9u
_	12u
_	12u
_	12u

050C322215	15"L
050C322218	18"L
050C322221	21"L
050C322224	24"L

— 9u	— 9u
	— 9u
_	— _{12u}
_	—
_	— _{12u}


D51C322230	30"L
D51C322236	36"L
D51C322242	42"L
D51C322248	48"L

	— 9u	— 9u
	— 9u	— 9u
	— _{12u}	— _{12u}
	— _{12u}	—
	— _{12u}	— _{12u}
ı		

D52C322230	30"L
D52C322236	36"L
D52C322242	42"L
D52C322248	48"L



D83C322215	15"l
D83C322218	18"l
D83C322221	21"l
D83C322224	24"l



D84C322215	15"L
D84C322218-	18"L
D84C322221-	21"L
D84C322224-	24"L







Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Cupboard and Drawer Cabinets






E05C322230L	30"L
E05C322236L	36"L
E05C322242L	42"L
E05C322248L	48"L

E05C322230	30"L
E05C322236	36"L
E05C322242	42"L
E05C322248-	48"L








— 9u
'
45u

E40C322215L	15"L
E40C322218L	18"L
E40C322221L	21"L
E40C322224L	24"L

E40C322215-	15"L
E40C322218-	18"L
E40C322221	21"L
E40C322224	24"L

E41C322230	30"L
E41C322236	36"L
E41C322242	42"L
E41C322248-	48"l



	— 9u
— _{15u}	
— 15u	:======
— 15u	45u

— 9u	— 9u
	— _{15u}
======	— 15u
45u	

— 9u	— 9u
- Ju	. 34
15u	11 1
	.======
15u	
— _{15u}	45u
1301	1 430

E42C322230L	30"L
E42C322236L	36"L
E42C322242L	42"L
E42C322248L	48"L

E42C322248- 481	E42C322230 E42C322236 E42C322242 E42C322248	30"L 36"L 42"L 48"L
------------------------	--	------------------------------

E43C322230L	30"L
E43C322236L	36"L
E43C322242L	42"L
E43C322248L	48"L

E43C322230	30"l
E43C322236	36"l
E43C322242	42"l
F43C322248-	48"I

ı	— 9u
I 	
	15u
54u	—

— 9u — 15u	I
	=======
	54u

	9u 	=	- 9u
=====	==:	:===:	
	15u		45u

E45C322230L	30"L
E45C322236L	36"L
E45C322242L	42"L
E45C322248L	48"L

E45C322230	30"L
E45C322236	36"L
E45C322242	42"L
E45C322248	48"L

E56C322230	30"L
E56C322236	36"L
E56C322242	42"L
E56C322248-	48"L

Research Collection

Steel ADA Height Base Cabinets




Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Cupboard and Drawer Cabinets



	9u
— 9u	lı
— _{12u}	
— _{12u}	
— _{12u}	45u

9u
— 9u
— _{12u}
— _{12u}
— _{12u}


— 9u	— 9u
— 9u	
— _{12u}	
— _{12u}	
— _{12u}	45u


E62C322230L	30"L
E62C322236L	36"L
E62C322242L	42"L
E62C322248L	48"L

E62C322248 - 48"L

E63C322242L 42 E63C322248L 48

E63C322230	30"L
E63C322236	36"L
E63C322242	42"L
E63C322248-	48"L






E65C322230L	30"L
E65C322236L	36"L
E65C322242L	42"L
E65C322248L	48"L

E65C322230	30"L
E65C322236	36"L
E65C322242	42"L
E65C322248	48"L







Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Blind Corner Cabinets









G12C322248- & G13C322248-

Acid Storage Cabinets







20" Deep G21C322018L_____ 18"L **G21C322024L**____ 24"L

20" Deep G21C320218-____ 18"L **G21C320224-___** 24"L

20" Deep G21C322030- 30"L G21C322036- 36"L G21C322042- 42"L G21C322048- 48"L

22" Deep G21C322218L 18"L **G21C322224L** 24"L

22" Deep G21C322218-____ 18"L **G21C322224-___** 24"L

 22" Deep

 G21C322230 30"L

 G21C322236 36"L

 G21C322242 42"L

 G21C322248 48"L


Note: G21 Acid Storage Cabinets are availble with five-knuckle hinges only

Diagonal Corner Cabinet









Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


Sink Cabinets



G00C322212L	12"L
G00C322215L	15"L
G00C322218L	18"L
G00C322221L	21"L
G00C322224L	24"L



G00C322212	12"I
G00C322215	15"I
G00C322218	18"l
G00C322221	21"l
G00C322224	24"l



G00C322242 42	5"L 2"L 3"L
---------------	-------------------

Sink Support Kits



Part Number	Fits
K-MSS-12	12" Cabinet
K-MSS-15	15" Cabinet
K-MSS-18	18" Cabinet
K-MSS-24	24" Cabinet
K-MSS-30	30" Cabinet
K-MSS-36	36" Cabinet
K-MSS-42	42" Cabinet
K-MSS-48	48" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support undermount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a -D to the end of the part number.


Example: K-MSS-36-D

Solvent Storage Cabinets




20" Deep G63C322018L_100 18"L **G63C322024L_100** 24"L

22" Deep		
G63C322218L 100	18"L	
G63C322224L 100	24"	



20" Deep	
G63C322018100	18
G63C322024100	24

22" Deep	
G63C322218- 100	18"L
G63C322224- ⁻ 100	24"L



20" Deep		
G63C322030- G63C322036- G63C322048-	100	30"L 36"L 48"L
22" Deep		

	-1-	
G63C322230-	100	30"L
G63C322236-	100	36"L
G63C322248-	100	48"L

G63 Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number. Example: **G63C322230-0100-B7** — Mobile, Inset style cabinet, 22" deep by 30" long.

Bin Storage Cabinet




G23C322218- 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.



Height Base Cabinets




Specifications:


ADA Height Fume Hood Cabinets are 32" high and either 20" or 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Fume Hood Base Cabinets







-0 -500p	
G08C322012L	12"L
G08C322015L	15"L
G08C322018L	18"L
G08C322021L	21"L
G08C322024L	24"L
22" Deep	
G08C322212L	12"L

G08C322215L

G08C322218L

G08C322221L

G08C322224L

G08C322012	12"L
G08C322015	15"L
G08C322018	18"L
G08C322021	21"L
G08C322024	24"L
22" Deep	
G08C322212	12"L
G08C322215-	15"L
G08C322218	18"L
G08C322221	21"L
G08C322224	24"L

· · · · ·	
G08C322030-	30"L
G08C322036	36"L
G08C322042	42"L
G08C322048	48"L

22" Deep	
G08C322230-	30"L
G08C322236-	36"L
G08C322242-	42"L
G08C322248-	48"L


Fume Hood Solvent Storage Base Cabinets


15"L


18"L

21"L

___ 24"L







G68C322024L 100 24"L G68C322024- 100

24"L

22" Deep G68C322224L 100 24"L G68C322224- 100


22" Deep 24"L G68C322030-_100 30"L G68C322036-_100 G68C322048-_100 48"L 22" Deep


G68C322230-_100 G68C322236-_100 G68C322248-_100

G68 Solvent Storage Cabinets are available as mobile. 36"L To add casters to the cabinets above, add -B7 to the end of the part number.

30"L Example: **G68C322018-0100-B7** 36"L — Mobile, Inset style cabinet, 22" 48"L deep by 18" long.

Fume Hood Vacuum Storage Base Cabinets





20" Deep

	17	
		-
þ	193	[
	淵	
	睴	
		1
ľ		45u
	\perp	45

G35C322024L 24"L G35C322024-24"L

22" Deep 22" Deep 24"L 24"L G35C322224L G35C322224-

20" Deep	
G35C322030	30"L
G35C322036	36"L
G35C322048	48"L
22" Deep	
G35C322230	30"L
G35C322236	36"L
G35C322248	48"L

Fume Hood Acid Storage Base Cabinets



20" Deep	
G80C322018L	24"L
G80C322024L	24"L


22" Deep G80C322218L 18"L **G80C322224L** 24"L



20" Deep	
G80C322018	18"l
G80C322024	24"l

 18"L
 G80C322218-____
 18"L

 24"L
 G80C322224-___
 24"L




20 Deep	
G80C322030-	30"L
G80C322036	36"L
G80C322048	48"L
22" Deep	
G80C322230-	30"L
G80C322236	36"L
G80C322248-	48"L


20" Deen

Note: G80 Acid Storage Cabinets are available with five-knuckle hinges only

Wheelchair Accessible Sink Cabinet with Doors




The Wheelchair Accessible Base Cabinet with Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and with the toespace rails attached to the doors so it will match adjacent cabinets when the doors are closed. The cabinet incorporates a 2³/₄" top rail to give a minimum of 29¹/₄" of clearance from the floor to the bottom of the rail.



G03C322236	36"L
G03C322242	42"L
G03C322248	48"L

Wheelchair Accessible Base Cabinet without Doors




The Wheelchair Accessible Base Cabinet without Doors is designed to be used as a sink unit, but may be use anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and without doors for easy accessibility.

Wheelchair Accessible Base Cabinets without Doors may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space, add a "C" to the end of the part number.


Example: G07C322236-0000-C

The cabinet incorporates a 4½" top rail to give a minimun of 27½" of clearance from the floor to the bottom of the rail and a removable sloping knee space panel to enclose piping and drain lines.



G07C322236-00_0 36"L **G07C322242-00_0** 42"L **G07C322248-00_0** 48"L






Specifications:


Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.

Open and Cupboard Cabinets are furnished with one adjustable shelf.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Open Cabinets





B00C352212-00_0 B00C352215-00_0 B00C352218-00_0B00C352221-00_0

12"L B00C3: 15"L B00C3: 18"L B00C3: 21"L B00C3: 24"L


B00C352230-00_0 30"L B00C352236-00_0 36"L B00C352242-00_0 42"L B00C352248-00_0 48"L Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.


Example: **G00C352236-0000-C**

Cupboard Cabinets

B00C352224-00 0







C03C352212L 12"L C03C352215L 15"L C03C352218L 18"L C03C352221L 21"L C03C352224L 24"L

 C03C352212-______
 12"L


 C03C352215-______
 15"L

 C03C352218-______
 21"L

 C03C352224-______
 24"L

60u

Drawer Cabinets





	— 20u
	— 20u
—	— 20u

]
5"L 8"L 21"L 24"L	D00C352230- D00C352236- D00C352242- D00C352248-	30"L 36"L 42"L 48"L


D02C352230	30"L
D02C352236	36"L
D02C352242	42"L
D02C352248	48"L

10u 10u 20u 20u

D08C352215-	15"L
D08C352218-	18"L
D08C352221-	21"L
D08C352224-	24"

_		
	_	— _{10u}
I	_	— _{10u}
	_	— 20u
	_	— 20u

D08C352230	30"L
D08C352236	36"L
D08C352242	42"L
D08C352248-	48"L




Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Drawer Cabinets




12"L
15"L
18"L
21"L
24"L

 D30C352230-____
 30"L

 D30C352236-____
 36"L

 D30C352242-____
 42"L

 D30C352248-____
 48"L



— 15u	— 15u
— 15u	— _{15u}
— 15u	— _{15u}
— 15u	— _{15u}

	— _{15u}
— _{15u}	— 15u
	— 15u
	— 15u

D31C352230	30"L
D31C352236	36"L
D31C352242	42"L
D31C352248	48"L

D32C352242- 42 D32C352248- 48
--

D36C352230-	30"L
D36C352236	36"L
D36C352242	42"L
D36C352248	48"L

_	12u
_	12u

— _{12u}	— _{12u}
_	— _{12u}
_	— _{12u}
_	—
_	— _{12u}

— _{12u}
120
— _{12u}
— _{12u}
— _{12u}
— _{12u}

D50C352215	15"L
D50C352218	18"L
D50C352221	21"L
D50C352224	24"L

D51C352230	30"L
D51C352236	36"L
D51C352242	42"L
D51C352248	48"L

D52C352230	30"L
D52C352236-	36"L
D52C352242	42"L
D52C352248	48"L

-	_	12u
-	_	
		24u
-		
		24u

_	9u
	12u
	12u
	24u

D81C352218	18"L
D81C352221	21"L
D81C352224	24"L

D84C352218-____ 18"L **D84C352221-**____ 21"L **D84C352224-**___ 24"L

2 File Drawers with hangers

Pull Out Shelf 1 File Drawer with hangers






Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.

Cupboards have one adjustable shelf.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection.
See page 4 for full explanation.
Actual product may vary slightly from catalog illustration due to style and hardware selection.

Cupboard and Drawer Cabinets




E05C352230L	30"L
E05C352236L	36"L
E05C352242L	42"L
E05C352248L	48"L

E05C352230	30"L
E05C352236	36"L
E05C352242	42"L
E05C352248-	48"L



	į
∥'	
	45u




_	
======	=======
	45u

	E40C352212L E40C352215L E40C352218L E40C352221L E40C352224L	12"L 15"L 18"L 21"L 24"L
--	---	--------------------------------------



E40C352248- 48"		E40C352230 E40C352236 E40C352242 E40C352248	30" 36" 42" 48"
------------------------	--	--	--------------------------

E41C352230-	30"L
E41C352236	36"L
E41C352242	42"L
E41C352248	48"L



_	
— 15u	I
— _{15u}	:======
— _{15u}	45u

—	— 15u
	— 15u
======	— 15u
45u	— 15u


15u	— 15u — 15u	
[15u[] 45u	— 15u — 15u	45u


E42C352230L	30"L
E42C352236L	36"L
E42C352242L	42"L
E42C352248L	48"L

E42U352248- 48	E42C352230 E42C352236 E42C352242 E42C352248	30"L 36"L 42"L 48"L
-----------------------	--	------------------------------

		30"l 36"l 42"l 48"l
--	--	------------------------------

E43C352230-	30"
E43C352236	36"
E43C352242	42"
F43C352248-	48"





15u	15u
l 45u	45u

E45C352230L	30"L
E45C352236L	36"L
E45C352242L	42"L
E45C352248L	48"L

	J
E45C352230	30"L
E45C352236-	36"L
E45C352242-	42"L
E45C352248-	48"L

30"L
36"L
42"L
48"L




Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.


Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Cupboard and Drawer Cabinets



E60C352212L 12"L E60C352215L 15"L E60C352218L 18"L E60C352221L 21"L E60C352224L 24"L







E62C352230L 30"L E62C352236L 36"L E62C352242L 42"L E62C352248L 48"L



E60C352212	12"L
E60C352215	15"L
E60C352218	18"L
E60C352221	21"L
E60C352224	24"L



E62C352230-30"L E62C352236-36"L E62C352242-42"L 48"L E62C352248-




•12u

E63C352230L 30"L E63C352236L 36"L E63C352242L 42"L E63C352248L 48"L

— _{12u}	
	======
1	48u
	<u> </u>


E61C352230-30"L E61C352236-36"L E61C352242-42"L E61C352248-48"L



E63C352230-30"L 36"L E63C352236-E63C352242-42"L E63C352248-48"L




E65C352230L 30"L E65C352236L 36"L E65C352242L 42"L E65C352248L 48"L



E65C352230-30"L E65C352236-36"L E65C352242-42"L E65C352248-48"L





Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Sink cabinets are provided with lowered back for plumbing access except where noted. Sink support kits are available when needed to support under-mount sinks.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


Sink Cabinets




G00C352212L	12"L
G00C352215L	15"L
G00C352218L	18"L
G00C352221L	21"L
G00C352224L	24"L



G00C352260-____ 60"L



G00C352212	12"L
G00C352215	15"L
G00C352218	18"L
G00C352221	21"L
G00C352224	24"L




G03C352230-	30"L
G03C352236-	36"L
G03C352242	42"L
G03C352248	48"L

G03 Cabinets are constructed without a bottom and with the toespace rail attached to the doors so it will match adjacent cabinets when the doors are closed.

Furnished with full height back



G00C352230	30"L
G00C352236-	36"l
G00C352242	42"L
G00C352248-	48"I
G00C332240	40 L



G04C352252	52"L
G04C352258	58"L
G04C352260	60"L

Center Bench End Sink Cabinet Furnished with full height back

Sink Support Kits



Fits 12" Cabinet 15" Cabinet 18" Cabinet
21" Cabinet
24" Cabinet 30" Cabinet
36" Cabinet
42" Cabinet
48" Cabinet
52" Cabinet 58" Cabinet
60" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support undermount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a -**D** to the end of the part number.

Example: K-MSS-36-D


Specifications

Acid Storage Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with 1/8" thick corrosion resistant polyethylene sheets.

The cabinet doors are provided with louvers and are latched using a nylon roller catch.

An optional 1" deep polyethylene pan is available to contain spills.

Each cabinet is furnished with a PVC vent pipe for venting. (Requires a 2" hole in the counter top above.)




Optional Polyethylene Pan

Seamless 1" deep polyethylene pan designed to contain spills.

Part Number	Fits 20" Deep
F-8440-15	15" Cabinet
F-8440-18	18" Cabinet
F-8440-21	21" Cabinet
F-8440-24	24" Cabinet
F-8440-30	30" Cabinet
F-8440-36	36" Cabinet
F-8440-42	42" Cabinet
F-8440-48	48" Cabinet
Part Number F-8441-15	Fits 22" Deep
Part Number F-8441-15	Fits 22" Deep 15" Cabinet
Part Number F-8441-15 F-8441-18	Fits 22" Deep 15" Cabinet 18" Cabinet
Part Number F-8441-15 F-8441-18 F-8441-21	Fits 22" Deep 15" Cabinet 18" Cabinet 21" Cabinet
Part Number F-8441-15 F-8441-18	Fits 22" Deep 15" Cabinet 18" Cabinet
Part Number F-8441-15 F-8441-18 F-8441-21 F-8441-24	Fits 22" Deep 15" Cabinet 18" Cabinet 21" Cabinet 24" Cabinet
Part Number F-8441-15 F-8441-18 F-8441-21 F-8441-24 F-8441-30	Fits 22" Deep 15" Cabinet 18" Cabinet 21" Cabinet 24" Cabinet 30" Cabinet

Acid Storage Cabinet



Note: G21 Acid Storage Cabinets are available with five-knuckle hinges only



Specifications


Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.

A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks and spills. The cabinet is provided with a full depth adjustable shelf and two diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18 gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped

with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.


All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.




Vent Locations



Solvent Storage Cabinets







20" Deep		Capacity		
G63C352024L	100	30"I	12 a	al

20" Deep	Capacity
G63C352024100	30"L 12 ga

G63C352030100	30"L	12 gal
G63C352036100		
G63C352048100	48"L	20 gal

Capacity

20" Deep


22" Deep		Capacity	
G63C352224L_10	00	30"L	12 gal

22" Deep	Capacity
G63C352224100	30"L 12 gal

22" Deep	Capacity
G63C352230- 100	30"L 12 gal
G63C352236100	36"L 15 gal
G63C352248- 100	48"I 20 dal

G63 Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number. Example: **G63C352036-0100-B7** — Mobile, Inset style cabinet, 20" deep by 36" long.




Specifications:


Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Blind Corner Cabinets have one adjustable shelf.

Note: Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Blind Corner Cabinets










G12C352248- & G13C352248-


Diagonal Corner Cabinets







Provided with Lazy Susan



Waste Storage Cabinet



G20C352215-G20C352218-18"L G20C352221-21"L G20C352224-21"L

Cabinet furnished with 13" deep, cold rolled steel bin, mounted on full extention drawer slides for ease of access. Includes 6" x 5" (opening) stainless steel waste chute

Bin Storage Cabinet




G23C352218-18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.



Fume Hood Base Cabinets




Specifications


Standing Height Steel Fume Hood Base Cabinets are 35" high and 22" deep; ADA Height Steel Fume Hood Base Cabinets are 32" high and 22" deep; lengths as shown.

Cabinet Style Option:


Specify cabinet door and hardware style when ordering by replacing blanks in last four digits of catalog number with style numbers.

Standing Height Fume Hood Base Cabinets









20" Deep	
G08C352012L	12"L
G08C352015L	15"L
G08C352018L	18"L
G08C352021L	21"L
G08C352024L	24"L

20" Deep	
G08C352012-	12"L
G08C352015	15"L
G08C352018	18"L
G08C352021	21"L
G08C352024	24"L

20" Deep	
30"L 36"L 42"L	
48"L	

20" Deep	
G08C352060	60"L

22" Deep	
G08C352212L	12"L
G08C352215L	15"L
G08C352218L	18"L
G08C352221L	21"L
G08C352224L	24"L

22" Deep	
G08C352212-	12"L
G08C352215	15"L
G08C352218	18"L
G08C352221	21"L
G08C352224	24"L

22" Deep	
G08C352230	30"L
G08C352236	36"L
G08C352242	42"L
G08C352248	48"L

22" Deep	
G08C352260	60"L


Acid Storage Fume Hood Base Cabinets

Specifications

Acid Storage Fume Hood Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with a molded one piece linear low density polyethylene tub with coved corners

and a 1" lip at the bottom front. The cabinet doors are lined with 1/8" sheet polyethylene and the doors are latched using a nylon roller catch.


Each cabinet is furnished with a PVC vent pipe for venting to the fume hood above. (Requires a 2" hole in the hood worksurface.)




Optional Phenolic Resin Removable Half-Depth Shelf

Part No.	
019623	Shelf for 18"cabinet
019624	Shelf for 24"cabinet
019625	Shelf for 30"cabinet
019626	Shelf for 36"cabinet
019627	Shelf for 48"cabinet


Standing Height Acid Storage Fume Hood Cabinets




Vent Location
On Center
115/8" from Top of Cabinet



Vent Location
On Center
115/8" from Top of Cabinet



Vent Location 30" - 75/8" Right of Center 36" & 48" - 5" Right of Center 115/8" from Top of Cabinet



Vent Location On Center 11⁵/8" from Top of Cabinet

20" Deep		
G80C352018L G80C352024L	18"L 24"L	G

20" Deep	
G80C352018	18"L
G80C352024	24"L

Zu Deeh	
G80C352030	30"L
G80C352036	36"L
G80C352048	48"L

20" Deep	
G80C352060	60"L

22" Deep		
G80C352218L	18"L	G8
G80C352224L	24"L	G8

22 Deep	
G80C352218 G80C352224	18"L 24"L
G60C352224	24 I

22" Deep		
G80C352230 G80C352236 G80C352248	30"L 36"L 48"L	

22" Deep	
G80C352260	60"L

Note: G80 Acid Storage Cabinets are available with five-knuckle hinges only



Solvent Storage Fume Hood Base Cabinets

Specifications


Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.

A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks and spills. The cabinet is provided with a full depth adjustable shelf and two diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18

gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.

All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.



Vent Locations



Fume Hood Solvent Storage Base Cabinets



20" Deep Capacity **G68C352024L 100** 30"L 12 gal

22" Deep Capacity

G68C352224L 100 30"L 12 gal

A FLAMMABE

A WALMARE

- DO - A MALMARE

- A MALMARE

20" Deep Capacity **G68C352024- 100** 30"L 12 gal

22" Deep Capacity **G68C352224- 100** 30"L 12 gal

= - A FLAMMARI - A MFLAMARI - A MFLAMARI

 20" Deep
 Capacity

 G68C352030-_100
 30"L 12 gal

 G68C352036-_100
 36"L 15 gal

 G68C352048- 100
 48"L 20 gal

 22" Deep
 Capacity

 G68C352230-_100
 30"L 12 gal

 G68C352236-_100
 36"L 15 gal

 G68C352248- 100
 48"L 20 gal

G68 Fume Hood Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G68C352224-1100-B7** — Mobile, Overlay style cabinet, 22" deep by 24" long.


Vacuum Pump Storage Base Cabinets

Specifications


Vacuum Pump Storage Fume Hood Cabinets are designed without a bottom to allow vacuum pumps and other equipment to be rolled in or out of the cabinets. The interior is lined with an easy to clean sound absorbant material. Each cabinet is furnished with


a 120 VAC, 20 amp, duplex receptacle mounted on the inside cabinet back and a pilot lighted toggle switch mounted in the top front rail. (Wiring is not included.) Cabinets are also furnished with a 11/2" diameter PVC vent pipe for venting or access to the


hood above. (Requires a 2" hole in the hood worksurface.) The toespace rail is attached to the door to allow total access to the cabinet. Cabinet inside clearance at the floor is $14^{1}/_{2}$ " or $16^{1}/_{2}$ " front-to-back, 25" high, and 3" less than the overall cabinet length.



Fume Hood Vacuum Pump Storage Storage Cabinets









20" Deep		20" Deep	
G35C352024L	24"L	G35C352024-	24"L

20" Deep	
G35C352030	30"L
G35C352036	36"L
G35C352048	48"l
	G35C352030 G35C352036

22" Deep		22" Deep	
G35C352224L	24"L	G35C352224	24"L

_	22" Deep	
-	G35C352230 G35C352236 G35C352248	30"L 36"L 48"L

Vacuum Pump Carts

Vacuum Pump Carts are designed for use with Kewaunee's Vacuum Pump cabinets shown on this page. The cart consists of a 1" deep pan, fabricated from powder coated or stainless steel and mounted on 4" swivel casters.



Carts for 20" Deep Cabinets

Part Number	Length
VPCC061414-00_0	14"
VPCC061420-00_0	20"
VPCC061426-00_0	26"
VPCC061432-00_0	32"
VPCC061438-00_0	38"
VPCC061444-00_0	44"

Carts for 22" Deep Cabinets

Part Number	Length
VPCC061614-00_0	14"
VPCC061620-00_0	20"
VPCC061626-00_0	26"
VPCC061632-00_0	32"
VPCC061638-00_0	38"
VPCC061644-00 0	44"



Sitting Height Mobile Cabinets



Specifications

Steel Sitting Height Mobile Cabinets are 261/4" high and 22" deep. Lengths as shown. Designed to fit under 30" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

Note: Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slighty from catalog illustration due to style and hardware selection.

Cupboard Cabinets





M10C262215L0000 M10C262218L0000 M10C262221L0000 M10C262224L0000 24"L

15"L 18"L 21"L

15"L

18"L

21"L

24"L

M10C262215-0000 15"L M10C262218-0000 18"L M10C262221-0000 21"L M10C262224-0000 24"L

Drawer/Cupboard Cabinets

Drawer/Cupboard Cabinets with Pullboards









M20C262215L M20C262218L M20C262221L M20C262224L

M20C262215-15"L M20C262218-18"L M20C262221-21"L M20C262224-24"L M25C262215L 15"L M25C262218L 18"L M25C262221L 21"L 24"L M25C262224L

M25C262215-15"L M25C262218-_ 18"L M25C262221-21"L M25C262224-24"L

Drawer Cabinets

Drawer Cabinets with Pullboards





M30C262215-15"L 18"L M30C262218-___ 21"L M30C262221-M30C262224-24"L



M62C262215-M62C262218-_ 18"L M62C262221-M62C262224-24"L





M45C262215-15"L M45C262218-____ 18"L 21"L M45C262221-M45C262224-24"L



M66C262215-15"L M66C262218-18"L M66C262221-21"L M66C262224-24"L

Steel Standing Height Mobile Cabinets



Specifications


Steel Standing Height Mobile Cabinets are 33" high and 22" deep. Lengths as shown. Designed to fit under 36" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

Note: Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slighty from catalog illustration due to style and hardware selection.

Cupboard Cabinets






M10C332215L0000 15"L M10C332218L0000 18"L M10C332221L0000 21"L M10C332224L0000 24"L

M10C332215-0000 15"L M10C332218-0000 18"L M10C332221-0000 21"L M10C332224-0000 24"L

Drawer/Cupboard Cabinets





15"L


18"L

21"L


24"L







Drawer/Cupboard Cabinets with Pullboards




M20C332215L M20C332218L M20C332221L M20C332224L M20C332215-15"L M20C332218-18"L M20C332221-21"L M20C332224-24"L

M25C332215L 15"L M25C332218L 18"L M25C332221L 21"L M25C332224L 24"L


M25C332215-15"L M25C332218-18"L M25C332221-21"L M25C332224-24"L

Drawer Cabinets





M40C332215-15"L M40C332218-18"L M40C332221-21"L M40C332224-24"L




M71C332215-15"L M71C332218-18"L M71C332221-21"L M71C332224-24"L


Drawer Cabinets with Pullboards



M46C332215-	15"L
M46C332218-	18"L
M46C332221-	21"L
M46C332224-	24"L



M56C332215-15"L M56C332218-18"L M56C332221-21"L M56C332224-24"L



M68C332215-15"L M68C332218-18"L M68C332221-21"L M68C332224-24"L



Mobile Cabinet Tops & Accessories

Specifications:


1" Thick Mobile Cabinet Tops are available in three depths, 22",23", or 24" and three materials: Kemresin, Phenolic Resin or Plastic Laminate. They may be used on any Mobile

Cabinet and add 1" to the overall height of the cabinets. All edges of top are finished.

1" Thick — Mobile Cabinet Top (adds 1" to height of Mobile Cabinet)




22" De	22" Deep flush to front of case							
		Kemresin	Black	Black				
Length	Kemresin®	Colors	Phenolic Resin	Plastic Laminate				
15"	F94R2215-00	BK = Black	F94T2215-00	F94P2215-00				
18"	F94R2218-00	GR = Grey	F94T2218-00	F94P2218-00				
21"	F94R2221-00	PT = Putty	F94T2221-00	F94P2221-00				
24"	F94R2224-00	SL = Slate	F94T2224-00	F94P2224-00				
30"	F94R2230-00		F94T2230-00	F94P2230-00				
36"	F94R2236-00		F94T2236-00	F94P2236-00				
42"	F94R2242-00		F94T2242-00	F94P2242-00				
48"	F94R2248-00		F94T2248-00	F94P2248-00				



End View




23" De	ep flush to overlay	door/drawer fr	onts - 1" overhang	of case	
		Kemresin	Black	Black	
Length	Kemresin®	Colors	Phenolic Resin	Plastic Laminate	
15"	F94R2315-00	BK = Black	F94T2315-00	F94P2315-00	
18"	F94R2318-00	GR = Grey	F94T2318-00	F94P2318-00	
21"	F94R2321-00	PT = Putty	F94T2321-00	F94P2321-00	
24"	F94R2324-00	SL = Slate	F94T2324-00	F94P2324-00	<u> </u>
30"	F94R2330-00		F94T2330-00	F94P2330-00	/ + O - / - i +
36"	F94R2336-00		F94T2336-00	F94P2336-00	Inset Cabinet
42"	F94R2342-00		F94T2342-00	F94P2342-00	End View
48"	F94R2348-00		F94T2348-00	F94P2348-00	



24" Deep 1" overhang of overlay door/drawer fronts							
Kemresin Black Black							
Length	Kemresin®	Colors	Phenolic Resin	Plastic Laminate			
15"	F94R2415-00	BK = Black	F94T2415-00	F94P2415-00			
18"	F94R2418-00	GR = Grey	F94T2418-00	F94P2418-00			
041	E04D0404 00	DT D	F04T0404 00	E04D0404 00			

F94R2421-00-__ **PT** = Putty F94T2421-00 F94P2421-00 24" **F94R2424-00-**__ **SL** = Slate F94T2424-00 F94P2424-00 30" F94R2430-00-F94T2430-00 F94P2430-00 36" F94R2436-00-F94T2436-00 F94P2436-00 42" F94R2442-00-F94T2442-00 F94P2442-00 48" F94R2448-00-__ F94T2448-00 F94P2448-00



Overlay Cabinet End View

Mobile Cabinet Tops & Accessories

Specifications:


Steel Mobile Cabinet Tops are fabricated from 18 gauge steel and finished to match the cabinet. They are designed for Steel Mobile Cabinets only.


Mobile Cabinet Grabbar allows easier movement and relocation of mobile cabinets. Fabricated from powder coated steel or stainless steel, they may be used on any

Mobile Cabinet, with or without Tops. Bars must be attached in field with supplied screws and/or bolts.

22" Deep — Steel — Mobile Cabinet Top (adds less than 1/4" to height of Mobile Cabinet)










22" Deep Top for Steel Mobile Cabinets only

Lg Painted Steel Stainless Steel F94S2215-00 15" **F94M2215-00** 18" **F94M2218-00** F94S2218-00 21" **F94M2221-00 F94S2221-00** 24" F94M2224-00 F94S2224-00 30" **F94M2230-00 F94S2230-00** 36" **F94M2236-00 F94S2236-00** 42" F94M2242-00 F94S2242-00

48" **F94M2248-00 F94S2248-00**








23" Deep — Steel — Mobile Cabinet Top (adds 1/2" to height of Mobile Cabinet)







Inset Cabinet **End View**

23" Deep Top for Steel Mobile Cabinets only

Lg Painted Steel Stainless Steel 15" **F94M2315-00** F94S2315-00 18" **F94M2318-00 F94S2318-00** 21" **F94M2321-00 F94S2321-00** 24" F94M2324-00 F94S2324-00 30" **F94M2330-00 F94S2330-00** 36" **F94M2336-00 F94S2336-00** 42" **F94M2342-00 F94S2342-00** 48" **F94M2348-00 F94S2348-00**







Mobile Cabinet Grabbar




Length	Painted Steel	Stainless Steel
15"	F73M0015-00	F73S0015-00
18"	F73M0018-00	F73S0018-00
21"	F73M0021-00	F73S0021-00
24"	F73M0024-00	F73S0024-00





Components of an Apron Assembly




Specifications:


Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.

All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

22" Deep Plain Rail Box Apron






Inset Style Apron Full Overlay Style Apron

All Plain Rail Box Aprons available in either Inset or Full Overlay Style

22" Deep Apron	.	22" Deep Apron	A	22" Deep Apron	A
A00C0422150 A00C0422160 A00C0422170	Apron Length 15" 16" 17"	A00C0422360 A00C0422370 A00C0422380	Apron Length 36" 37" 38"	A00C0422570 A00C0422580 A00C0422590	Apron Length 57" 58" 59"
A00C0422180	18"	A00C0422390	39"	A00C0422600	60"
A00C0422190	19"	A00C0422400	40"	A00C0422610	61"
A00C0422200	20"	A00C0422410	41"	A00C0422620	62"
A00C0422210	21"	A00C0422420	42"	A00C0422630	63"
A00C0422220	22"	A00C0422430	43"	A00C0422640	64"
A00C0422230	23"	A00C0422440	44"	A00C0422650	65"
A00C0422240	24"	A00C0422450	45"	A00C0422660	66"
A00C0422250	25"	A00C0422460	46"	A00C0422670	67"
A00C0422260	26"	A00C0422470	47"	A00C0422680	68"
A00C0422270	27"	A00C0422480	48"	A00C0422190	69"
A00C0422280	28"	A00C0422490	49"	A00C0422700	70"
A00C0422290	29"	A00C0422500	50"	A00C0422710	71"
A00C0422300 A00C0422310 A00C0422320	30" 31" 32"	A00C0422510 A00C0422520 A00C0422530	51" 52" 53"	A00C0422720	72"
A00C0422330 A00C0422340 A00C0422350	33" 34" 35"	A00C0422540 A00C0422550 A00C0422560	54" 55" 56"		

Components of an Apron Assembly




Specifications:


Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.

All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

28" Deep Plain Rail Box Apron






All Plain Rail Box Aprons available in either Inset or Full Overlay Style

28" Deep Apron	Apron	28" Deep Apron	Apron	28" Deep Apron	Apron
A00C0428150 A00C0428160 A00C0428170	Length 15" 16" 17"	A00C0428360 A00C0428370 A00C0428380	26" 36" 37" 38"	A00C0428570 A00C0428580 A00C0428590	Length 57" 58" 59"
A00C0428180	18"	A00C0428390	39"	A00C0428600	60"
A00C0428190	19"	A00C0428400	40"	A00C0428610	61"
A00C0428200	20"	A00C0428410	41"	A00C0428620	62"
A00C0428210	21"	A00C0428420	42"	A00C0428630	63"
A00C0428220	22"	A00C0428430	43"	A00C0428640	64"
A00C0428230	23"	A00C0428440	44"	A00C0428650	65"
A00C0428240	24"	A00C0428450	45"	A00C0428660	66"
A00C0428250	25"	A00C0428460	46"	A00C0428670	67"
A00C0428260	26"	A00C0428470	47"	A00C0428680	68"
A00C0428270	27"	A00C0428480	48"	A00C0428190	69"
A00C0428280	28"	A00C0428490	49"	A00C0428700	70"
A00C0428290	29"	A00C0428500	50"	A00C0428710	71"
A00C0428300 A00C0428310 A00C0428320	30" 31" 32"	A00C0428510 A00C0428520 A00C0428530	51" 52" 53"	A00C0428720	72"
A00C0428330 A00C0428340 A00C0428350	33" 34" 35"	A00C0428540 A00C0428550 A00C0428560	54" 55" 56"		



Components of an Apron Assembly




Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.

All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

22" Deep Box Apron with Drawers




1 Drawer 2 Drawers 3 Drawers

22" Deep Apron

	Apron Length	Dwr Size
A51C042220-1	20"	18"
A51C042224-1	24"	21"
A51C042226-1	26"	24"
A51C042230-1	30"	27"
A51C042232-1	32"	30"
A51C042236-1	36"	33"

22" Deep A	pron		22" Deep Apron			
	Apron Length	Dwr Size		Apron Length	Dwr Size	
A52C042226-1 A52C042230-1 A52C042232-1 A52C042236-1 A52C042238-1 A52C042242-1 A52C042244-1	30" 32" 36" 38"	12" 12" 15" 15" 18" 18" 21"				
A52C042248-1 A52C042250-1 A52C042254-1 A52C042256-1 A52C042260-1 A52C042262-1 A52C042268-1 A52C042272-1	56" 60" 62" 66"	21" 24" 24" 27" 27" 30" 30" 33" 33"	A53C042250-1 A53C042254-1 A53C042256-1 A53C042260-1 A53C042262-1 A53C042268-1 A53C042272-1	50" 54" 56" 60" 62" 66" 68" 72"	15" 15" 17" 17" 19" 21" 21"	

Components of an Apron Assembly




Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.

All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

28" Deep Box Apron with Drawers



1 Drawer	2 Drawers	3 Drawers

28" Deep Apron

	Apron Length	Dwr Size
A51C042820-1	20"	18"
A51C042824-1	24"	21"
A51C042826-1	26"	24"
A51C042830-1	30"	27"
A51C042832-1	32"	30"
A51C042836-1	36"	33"

28" Deep A	pron		28" Deep Apron			
	Apron Length	Dwr Size		Apron Length	Dwr Size	
	00"	4.01				
A52C042826-1 A52C042830-1	26" 30"	12" 12"				
A52C042832-1 A52C042836-1	32" 36"	15" 15"				
A52C042838-1 A52C042842-1	38" 42"	18" 18"				
A52C042844-1 A52C042848-1		21" 21"				
A52C042850-1 A52C042854-1		24" 24"	A53C042850-1 A53C042854-1	50" 54"	15" 15"	
A52C042856-1 A52C042860-1	56" 60"	27" 27"	A53C042856-1 A53C042860-1	56" 60"	17" 17"	
A52C042862-1 A52C042866-1	62" 66"	30" 30"	A53C042862-1 A53C042866-1	62" 66"	19" 19"	
A52C042868-1 A52C042872-1	68" 72"	33" 33"	A53C042868-1 A53C042872-1	68" 72"	21" 21"	




Specifications:

corner work tops as well as the extra steel panels when specified with wood storage space of an apron drawer. The drawer heads. blind end may either be attached to the wall or supported by a leg or pedestal assembly.


Corner Box Aprons provide support for Corner Box Aprons are furnished with


2" High Box Aprons are designed for use when knee space height is critical. They are available 22" deep and in lengths from 15" to 72" on 1" increments. 2" High Box Aprons are not designed for use with legs.


Corner Box Apron



Furnished with steel panels when specified with wood drawer heads.







Size

	Length	Size
A05C042248-1	48"	21"
A05C042254-1	54"	21"
A05C042260-1	60"	21"
A05C042266-1	66"	21"
A05C042272-1	72"	21"

28" Deen			28" Deen		
2-1	72"	21"	A04C042272-1	72"	2
60-1 66-1	60" 66"	21" 21"	A04C042260-1 A04C042266-1	60" 66"	2
94-I	54	21	AU4CU42254-1	54	2

A04C042248-1

	Apron Length	Dwr Size
A05C042848-1	48"	21"
A05C042854-1	54"	21"
A05C042860-1	60"	21"
A05C042866-1	66"	21"
A05C042872-1	72"	21"

	Apron Length	Dwr Size
A04C042848-1	48"	21"
A04C042854-1	54"	21"
A04C042860-1	60"	21"
A04C042866-1	66"	21"
A04C042872-1	72"	21"

Specifications:


Apron Rails may be used in place of Box Aprons to provide worksurface support between two units. They are available 2" or 4" high and from 15" to 72" long in 1" increments. They may be installed vertically to provide a 2"

or 4" high rail, or turned 90 degrees to provide a 11/2" high rail. Apron rails should be ordered in pairs to provide both front and rear top support.

Work Top Support Rails and Assemblies may be used anywhere

maximum knee space is wanted but top support is required. Fabricated of heavy gauge steel channels, it is only 7/8" high. The Keyboard Tray may be added to any of the frames or may be used independently.

Inset 4" High Apron Rail




A06C0402 ____-**00_0** 4" High

Inset Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

Example: A06C0402**48**-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

Overlay 4" High Apron Rail




A06C0403 ____-**1**___**0** 4" High

Overlay Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

Example: A06C0402**48**-1610 is a 48" long apron rail with Oak overlay panel.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

4" High Adjustable Rail




A56C040236-00_0 36"L **A56C040248-00_0** 48"L **A56C040260-00_0** 60"L **A56C040272-00_0** 72"L

Same as A06 Apron Rails except furnished with a loose angle at one end to allow the rail to be cut in the field to fit available space.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

2" High Plain Rail Box Apron




A00C0222 -00 0

2" high Box Aprons are available from 15" to 72" long in 1" increments. To specify apron, insert the length required into the part number.

Example: A00C0202**48**-0000 is a 48" long 2" high cold rolled steel box apron.

Not for use with legs

2" High Apron Rail




A06C0202 ____-**00_0** 2" High

Apron rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

Example: A06C0202**48**-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

Support Rails



A81C0102 ____-0000

Support rails may be used individually or assembled to form a support framework. Rails are available from 12" to 72" long in 1" increments and are furnished with angles attached to each end. To specify rail, insert the length required in the part number.

Example: A810102**48**-0000 is a 48" long support rail.

Support rails should be ordered in pairs to provide both front and rear worksurface support.



Steel Leg Assemblies & Pedestals

Specifications:

Legs are 2" square and all except the Top Mounted Leg are joined in pairs with a 2" high stretcher and a 4" high apron rail. Dimensions are as shown.


For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a


spreader is used, order the A08 or A28 leg assembly which are designed for use with a spreader.


All legs may be ordered with 4" swivel locking casters by adding -B10 to the end of the part numbers. Example: A07C292202-0000-B10


Adjustable Height Legs are adjustable in 1" increments from 29" to 35" high.


End Leg Assemblies











For use with Spreader

223/16" Deep	
A07C292202-00_0	29"H
A07C322202-00_0 A07M352202-00 0	32"H 35"H
_	30 H
28 ³ / ₁₆ " Deep	
A07C292802-00_0	29"H 32"H
A07C322802-00_0 A07C352802-00_0	32 ⊓ 35"H

223/16" Deep A08C292202-00_0 29"H A08C322202-00_0 32"H A08C352202-00 0 35"H 28³/₁₆" Deep A08C292802-00 0 29"H A08C322802-00 0 32"H A08C352802-00 0 35"H


223/16" Deep A27C292202-00 0 29"H A27C322202-00_0 32"H 35"H A27C352202-00 0 283/16" Deep A27C292802-00 0 29"H A27C322802-00 0 32"H A27C352802-00 0 35"H

Adjustable Height Intermediate Leg Ass'y

223/16" Deep A28C292202-00_0 A28C322202-00_0 32"H A28C352202-00 0 35"H 283/16" Deep A28C292802-00 0 29"H A28C322802-00 0 32"H A28C352802-00 0 35"H


Each Leg Assembly adds 2" to the overall length of the apron assembly to which it is attached..

Adjustable Height End Leg Assemblies





283/16" Deep A57C352802-00 0 Adjustable Height 29"-35"




For use with Spreader 223/16" Deep A58C352202-00 0

Adjustable Height 29"-35


Adj.

283/16" Deep A58C352802-00 0 Adjustable Height 29"-35"



223/16" Deep A77C352202-00 0 Adj. Adjustable Height 29"-35"

283/16" Deep A77C352802-00 0 Adj. Adjustable Height 29"-35"



For use with Spreader

223/16" Deep A78C352202-00 0 Adjustable Height 29"-35

283/16" Deep A78C352802-00 0 Adjustable Height 29"-35"

Steel Apron Rails & Panels

Specifications:


Knee Space panels are designed to close off the space between the floor and the bottom of the apron or countertop. They are sized for 29" high, 32" high, and 35" high cabinets. Each Knee Space Panel is made with a seperate base rail, to allow the removal

of the panel without removing the base molding.

Knee Space Panels are furnished with two angles to attach to panel to the cabinets.

Pedestal Legs are 2" wide, depths and heights as shown. They are designed to be used with or without Box Aprons to form knee spaces at the end of assemblies.

Knee Space Panels



Knee space panels are available from 15" to 72" long in 1" increments. To specify knee space panel, insert the length required into the part number.

Example: A10C3100**48**-0000 is a 48" long knee space panel.

Sitting Height Knee Space Panels

Apron Part Number	Height	Length	for use with:
A10C250000_0	25"	15"-72"L	4" high apron
A10C280000_0	28"	15"-72"L	A80 support frame
A10C290000_0	29"	15"-72"L	no apron

ADA Height Knee Space Panels


Apron Part Number	Height	Length	for use with:
A10C280000_0	28"	15"-72"L	4" high apron
A10C310000_0	31"	15"-72"L	A80 support frame
A10C320000_0	32"	15"-72"L	no apron or frame


Standing Height Knee Space Panels


Apron Part Number	Height	Length	for use with:
A10C310000_0	31"	15"-72"L	4" high apron
A10C340000_0	34"	15"-72"L	A80 support frame
A10C350000_0	35"	15"-72"L	no apron or frame


Top Mounted Leg Wall Pedestal Leg Center Pedestal Leg Spreader

29"H









Left Side 22" Deep

Left Side 28" Deep

Right Side 22" Deep

A11C292802-00 0

A11C322802-00 0

A11C352802-00_0

A12C292202-00 0

A12C322202-00 0

A12C352802-00_0

A11C292202-00_0	29"H	A13C295202-00_0	29"H
A11C322202-00_0	32"H	A13C325202-00_0	32"H
A11C352202-00_0	35"H	A13C355202-00_0	35"H

29"H

32"H

35"H

29"H

32"H

35"H

A09C020230-00 0 30"L **A09C020236-00_0** 36"L **A09C020242-00_0** 42"L **A09C020248-00 0** 48"L **A09C020254-00_0** 54"L **A09C020260-00_0** 60"L **A09C020266-00 0** 66"L **A09C020272-00_0** 72"L

32"H A15C320202-00_0 A15C350202-00_0 35"H

A15C290202-0000

Adj. A65C350202-00 0 Adjustable Height 29"-35"


> 2" square leg welded to 4" square steel plate

Each Pedestal Leg adds 2" to the overall length of the assembly to which it is attached.

Pedestal Legs may be ordered with the topspace 7/8" deeper so that the front edge lines up with the doors and drawer fronts on adiacent Full Overlay cabinets. To specify the deeper toe space, add a **C** to the end of the part

Example: A11C352202-0000-C





Specifications:

Open Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12³/₄" and 16". Lengths are as shown. All are furnished with one, two, or three adjustable shelves as shown.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

Open Storage

	12 ³ / ₄ " Deep			123/4" Deep	
24" High	W00C241315-00_0 W00C241318-00_0 W00C241321-00_0 W00C241324-00_0	15"L 18"L 21"L 24"L	24" High	W00C241330-00_0 W00C241336-00_0 W00C241342-00_0 W00C241348-00_0	30"L 36"L 42"L 48"L
	16" Deep			16" Deep	
	W00C241615-00_0 W00C241618-00_0 W00C241621-00_0 W00C241624-00_0	15"L 18"L 21"L 24"L		W00C241630-00_0 W00C241636-00_0 W00C241642-00_0 W00C241648-00_0	30"L 36"L 42"L 48"L
30" High	12 ³ / ₄ " Deep		30" High	12 ³ / ₄ " Deep	
	W00C301315-00_0 W00C301318-00_0 W00C301321-00_0 W00C301324-00_0	15"L 18"L 21"L 24"L		W00C301330-00_0 W00C301336-00_0 W00C301342-00_0 W00C301348-00_0	30"L 36"L 42"L 48"L
	16" Deep			16" Deep	
	W00C301615-00_0 W00C301618-00_0 W00C301621-00_0 W00C301624-00_0	15"L 18"L 21"L 24"L		W00C301630-00_0 W00C301636-00_0 W00C301642-00_0 W00C301648-00_0	30"L 36"L 42"L 48"L
36" High	12 ³ / ₄ " Deep		36" High	12 ³ / ₄ " Deep	
	W00C361315-00_0 W00C361318-00_0 W00C361321-00_0 W00C361324-00_0	15"L 18"L 21"L 24"L		W00C361330-00_0 W00C361336-00_0 W00C361342-00_0 W00C361348-00_0	30"L 36"L 42"L 48"L
	16" Deep			16" Deep	
	W00C361615-00_0 W00C361618-00_0	15"L 18"L		W00C361630-00_0	30"L 36"L
	W00C361621-00_0 W00C361624-00_0	21"L 24"L		W00C361636-00_0 W00C361642-00_0 W00C361648-00_0	42"L 48"L
48" High		21"L	48" High	W00C361642-00_0	42"L
48" High	W00C361624-00_0	21"L	48" High	W00C361642-00_0 W00C361648-00_0	42"L
48" High	W00C361624-00_0 12 ³ / ₄ " Deep W00C481315-00_0 W00C481318-00_0 W00C481321-00_0	21"L 24"L 15"L 18"L 21"L	48" High	W00C361642-00_0 W00C361648-00_0 123/4" Deep W00C481330-00_0 W00C481336-00_0 W00C481342-00_0	42"L 48"L 30"L 36"L 42"L



Specifications:

Plate Glass Door Wall Mounted Storage Cabinets are available in three heights, 24", 30", 36" and 48" and two depths; 123/4" and 16". Lengths are as shown. Cabinets are furnished with one, two, or three adjustable shelves as shown.


Doors are ¹/₄" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

Sliding Plate Glass Doors — Clear Glass

Wall Mounted Cabinet Hanging Rails

24" High




. – / . – сор	
W10C241330-00_0	30"L
W10C241336-00_0	36"L
W10C241342-00 0	42"L
W10C241348-00 0	48"I

16" Deep

123/4" Deep

. с 200р	
W10C241630-00_0	30"L
W10C241636-00_0	36"L
W10C241642-00_0	42"L
W10C241648-00 0	48"L

30" High




12 3/4"	Deep
----------------	------

W10C301330-00_0	30"L
W10C301336-00_0	36"L
W10C301342-00_0	42"L
W10C301348-00_0	48"L

16" Deep

W10C301630-00_0	30"L
W10C301636-00_0	36"L
W10C301642-00_0	42"L
W10C301648-00_0	48"L

36" High



123/4" Deep

W10C361330-00_0 W10C361336-00_0 W10C361342-00_0 W10C361348-00_0	30"L 36"L 42"L 48"I
16" Deep	.0 _
W10C361630-00_0	30"L

36"L


42"L

W10C361636-00 0

W10C361642-00_0


W10C361648-00_0

48"	High	



123/4" Deep


W10C481330-00_0	30"L
W10C481336-00 0	36"L
W10C481342-00 0	42"L
W10C481348-00_0	48"L
16" Deep	
W10C481630-00_0	30"L
W10C481636-00 0	36"L
W10C481642-00 0	42"L
W10C481648-00_0	48"L



L01-020109-00	for 15" wall cabined
L01-020112-00	for 18" wall cabined
L01-020115-00	for 21" wall cabined
L01-020118-00	for 24" wall cabined
L01-020124-00	for 30" wall cabined
L01-020130-00	for 36" wall cabined
L01-020136-00	for 42" wall cabined
L01-020142-00	for 48" wall cabined

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.





Specifications:

Sliding Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12³/₄" and 16". Lengths are as shown. Cabinets are glazed with ¹/₈" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown. Doors are

mounted on friction-free roller glides.


Sliding Door Wall Mounted Storage Cabinets are available in Door Style 1 and 6 only.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

Sliding Glazed Doors

Sliding Panel Doors

24" High




/4 Doop	123/4" Deep
---------	-------------

W41C241330-00	30"L
W41C241336-0 0	36"I
W41C241342-00	42"L
W41C241348-0 0	48"L
16" Deep	

10 Беер		
W41C241630-0	0	30"L
W41C241636-0	0	36"L
W41C241642-0	0	42"L
W41C241648-0	n	48"I

24" High




123/4" Deep

W12C241330-00	30"L
W12C241336-00	36"L
W12C241342-00	42"L
W12C241348-0_0	48"L
16" Deep	

W12C241630-00	30"L
W12C241636-00	36"L
W12C241642-00	42"L
W12C241648-0 0	48"I

30" High




12 ³ /4"	Deep
---------------------	------

W41C301330-00 W41C301336-00 W41C301342-00 W41C301348-00	30"L 36"L 42"L 48"L
16" Deep	
W41C301630-00	30"L

W41C301630-0	_0	30"L
W41C301636-0	_0	36"L
W41C301642-0	_0	42"L
W41C301648-0	_0	48"L

30" High




123/4" Deep

W12C301330-00	30"
W12C301336-00	36"
W12C301342-00	42"
W12C301348-00	48"l
16" Door	

16" Deep

•		
W12C301630-0_	0	30"L
W12C301636-0_	0	36"L
W12C301642-0	0	42"L
W12C301648-0	0	48"L

36" High




123/4" Deep

W41C361330-00 W41C361336-00 W41C361342-00 W41C361348-00	30"L 36"L 42"L 48"L
16" Deep	
W41C361630-00	30"L

30"L
36"L
42"L
48"L

36" High



123/4" Deep

•	
W12C361330-00	30"L
W12C361336-00	36"L
W12C361342-00	42"L
W12C361348-00	48"L
16" Deep	
W12C361630-00	30"L

36"L

42"L

48"L

__0

48" High




123/4" Deep

W41C481330-0_	_0	30"L
W41C481336-0_	_0	36"L
W41C481342-0_	_0	42"L
W41C481348-0	_0	48"L
16" Deep		
W41C481630-0	_0	30"L
W41C481636-0	0	36"L
W41C481642-0		42"I

W41C481648-0__0 48"L

48" High



123/4" Deep

W12C361636-0

W12C361642-0

W12C361648-0

12 /4 Deep	
W12C481330-00 W12C481336-00 W12C481342-00 W12C481348-00	30"L 36"L 42"L 48"L
16" Deep	

W12C481630-0	0	30"L
W12C481636-0	0	36"L
W12C481642-0	0	42"L
W12C481648-0	_0	48"L



Specifications:


Double Sided Plate Glass Pass-Thru Cabinets are available 24", 30", and 36" high, 123/4" and 16" deep and in lengths shown. Double Sided Pass-Thru Cabinets have doors on both sides and are designed to be mounted on center bench reagent racks.

Doors are 1/4" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

Double Sided Pass-Thru Cabinet - Clear Glass

24" High




12 ³ / ₄ "	Deep
----------------------------------	------

W13C241330-00_0	30"L
W13C241336-00_0	36"L
W13C241342-00_0	42"L
W13C241348-00_0	48"L
16" Doop	

W13C241630-00_0	30"L
W13C241636-00_0	36"L
W13C241642-00_0	42"L
W13C241648-00 0	48"L

30" High




12 3/4	" D	ee	p
---------------	-----	----	---

W13C301330-00_0	30"L
W13C301336-00_0	36"L
W13C301342-00 0	42"L
W13C301348-00_0	48"L

16" Deep

•	
W13C301630-00_0	30"L
W13C301636-00_0	36"L
W13C301642-00_0	42"L
W13C301648-00 0	48"L

36" High



123/4" Deep

W13C361330-00_0	30"L
W13C361336-00_0	36"L
W13C361342-00_0	42"L
W13C361348-00_0	48"L

16" Deep

W13C361630-00_0	30"L
W13C361636-00_0	36"L
W13C361642-00_0	42"L
W13C361648-00 0	48"l





Specifications:

Swinging Glazed Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 123/4" and 16". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown.

Note: Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.

Swinging Glazed Doors

24" High	W50C241315L W50C241318L W50C241321L W50C241324L W50C241615L W50C241618L W50C241621L W50C241624L	W50C241321- W50C241324- W50C241615- W50C241618-	18"L 21"L 24"L 16" D 15"L 18"L 21"L	24" High	12³/4" Deep W50C241330
30" High	W50C301315L W50C301315L W50C301321L W50C301324L W50C301615L W50C301618L W50C301621L W50C301624L	W50C301618 W50C301621-	18"L 21"L 24"L 16" D 15"L 18"L 21"L	30" High	12³/4" Deep W50C301330 30"L W50C301336 42"L W50C301342 48"L 16" Deep W50C301630 30"L W50C301636 36"L W50C301642 42"L W50C301648 48"L
36" High	W50C361315L W50C361318L W50C361321L W50C361324L W50C361615L W50C361618L W50C361621L W50C361624L	W50C361324 W50C361615 W50C361618 W50C361621	18"L 21"L 24"L 16" D 15"L 18"L 21"L	36" High	12³/4" Deep W50C361330- 30"L W50C361336- 42"L W50C361348- 48"L 16" Deep W50C361630- 30"L W50C361636- 36"L W50C361648- 42"L W50C361648- 48"L



Specifications:


Swinging Panel Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 123/4" and 16". Lengths are as shown. Cabinets are furnished with one, two or three adjustable shelves as shown.

Note: Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.

Swinging Panel Doors

	Left Hand	Right Hand	12³/4" D		12 ³ / ₄ " Deep	
24" High	W25C241315L	W25C241315	15"L	24" High	W25C241330	30"L
	W25C241318L	W25C241318 W25C241321-	18"L 21"L		W25C241336 W25C241342-	36"L
	W25C241321L W25C241324L	W25C241321 W25C241324	21 L 24"L		W25C241342 W25C241348	42"L 48"L
	172002110212	W200211021	16" D		16" Deep	10 L
	-				<u> </u>	
	W25C241615L	W25C241615	15"L		W25C241630	30"L
	W25C241618L W25C241621L	W25C241618 W25C241621			W25C241636 W25C241642	36"L 42"L
	W25C241621L	W25C241621			W25C241648-	48"L
						10 2
30" High	Left Hand	Right Hand	12³/4" D	30" High	123/4" Deep	
g	W25C301315L	W25C301315	15"L	00 IIIg.i	W25C301330	30"L
	W25C301318L	W25C301318			W25C301336	36"L
	W25C301321L	W25C301321			W25C301342	42"L
	W25C301324L	W25C301324			W25C301348	48"L
			16" D		16" Deep	
	W25C301615L	W25C301615			W25C301630	30"L
	W25C301618L	W25C301618			W25C301636	36"L
	W25C301621L W25C301624L	W25C301621 W25C301624-	21"L 24"L		W25C301642 W25C301648-	42"L
	W25C3U1624L	W25C3U1624	24 L		W25C3U1040	48"L
36" High	Left Hand	Right Hand	12 ³ / ₄ " D	36" High	12 ³ / ₄ " Deep	
36" High	W25C361315L	W25C361315	15"L	36" High	W25C361330	30"L
36" High	W25C361315L W25C361318L	W25C361315 W25C361318	 15"L 18"L	36" High	W25C361330 W25C361336	36"L
36" High	W25C361315L W25C361318L W25C361321L	W25C361315 W25C361318 W25C361321	15"L 18"L 21"L	36" High	W25C361330 W25C361336 W25C361342	36"L 42"L
	W25C361315L W25C361318L	W25C361315 W25C361318	15"L 18"L 21"L 24"L	36" High	W25C361330 W25C361336 W25C361342 W25C361348	36"L
	W25C361315L W25C361318L W25C361321L	W25C361315 W25C361318 W25C361321	15"L 18"L 21"L	36" High	W25C361330 W25C361336 W25C361342	36"L 42"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615-	15"L 18"L 21"L 24"L 16" D	36" High	W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630-	36"L 42"L 48"L 30"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361618-	15"L 18"L 21"L 24"L 16" D	36" High	W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636-	36"L 42"L 48"L 30"L 36"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361618- W25C361621-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L	36" High	W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642-	36"L 42"L 48"L 30"L 36"L 42"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361618-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L	36" High	W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636-	36"L 42"L 48"L 30"L 36"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361618- W25C361621-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L	36" High	W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642-	36"L 42"L 48"L 30"L 36"L 42"L
 	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L W25C361624L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361618- W25C361621- W25C361624- Right Hand	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L 21"L 24"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642- W25C361648- 12 ³ / ₄ " Deep	36"L 42"L 48"L 30"L 36"L 42"L 48"L
 	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L W25C361624L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 24"L 123/4" D 15"L 18"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642- W25C361648-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L W25C361624L Left Hand W25C481315L W25C481318L W25C481321L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318- W25C481321-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L 24"L 123/4" D 15"L 18"L 21"L 24"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642- W25C361648- 12 ³ / ₄ " Deep W25C481330- W25C481336- W25C481342-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L 42"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361618L W25C361621L W25C361624L Left Hand W25C481315L W25C481318L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 21"L 24"L 123/4" D 15"L 18"L 21"L 24"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642- W25C361648- W25C361648- W25C361648- W25C361648-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361621L W25C361624L W25C361624L W25C361624L W25C481315L W25C481318L W25C481321L W25C481324L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318- W25C481321- W25C481324-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 24"L 123/4" D 15"L 18"L 24"L 16" D		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361636- W25C361642- W25C361648- 12 ³ / ₄ " Deep W25C481330- W25C481336- W25C481342- W25C481348- 16" Deep	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L 42"L 48"L
	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361621L W25C361621L W25C361624L Left Hand W25C481315L W25C481318L W25C481321L W25C481324L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318- W25C481321- W25C481324- W25C481615-	15"L 18"L 21"L 24"L 16" D 15"L 24"L 123/4" D 15"L 18"L 24"L 16" D 15"L 18"L 18"L 18"L 18"L 18"L 18"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361642- W25C361648- W25C361648- 123/4" Deep W25C481330- W25C481336- W25C481342- W25C481348- 16" Deep W25C481630-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L 42"L 48"L
48" High	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361621L W25C361624L W25C361624L Left Hand W25C481315L W25C481318L W25C481324L W25C481324L W25C481615L W25C481615L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318- W25C481321- W25C481324- W25C481615- W25C481615- W25C481618-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 24"L 123/4" D 15"L 18"L 24"L 16" D 15"L 18"L 18"L 18"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361642- W25C361648- W25C361648- 123/4" Deep W25C481330- W25C481336- W25C481342- W25C481348- 16" Deep W25C481630- W25C481636-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L 48"L
48" High	W25C361315L W25C361318L W25C361321L W25C361324L W25C361615L W25C361621L W25C361621L W25C361624L Left Hand W25C481315L W25C481318L W25C481321L W25C481324L	W25C361315- W25C361318- W25C361321- W25C361324- W25C361615- W25C361621- W25C361624- Right Hand W25C481315- W25C481318- W25C481321- W25C481324- W25C481615-	15"L 18"L 21"L 24"L 16" D 15"L 18"L 24"L 123/4" D 15"L 18"L 21"L 24"L 16" D 15"L 24"L		W25C361330- W25C361336- W25C361342- W25C361348- 16" Deep W25C361630- W25C361642- W25C361648- W25C361648- 123/4" Deep W25C481330- W25C481336- W25C481342- W25C481348- 16" Deep W25C481630-	36"L 42"L 48"L 30"L 36"L 42"L 48"L 30"L 36"L 42"L 48"L







Specifications:

Corner Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12¾" and 16" and with either glazed or panel swinging doors. Cabinets are furnished with one, two or three adjustable shelves as shown.

Note: Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.

Swinging Door Corner Cabinets






Specifications:

Full Height Storage Cabinets are 84" high and 123/4", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

CAUTION: Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

NOTE: Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.


Open Storage



12%4" Deep	
S00C841318-00_0	18"L
S00C841324-00 0	24"1


16" Deep	
S00C841618-00_0	18"L
S00C841624-00_0	24"L

22" Deep	
S00C842218-00_0	18"L
S00C842224-00_0	24"L



12¾" Deep	
\$00C841330-00_0	30"L
\$00C841336-00_0	36"L
\$00C841342-00_0	42"L
\$00C841348-00_0	48"L
16" Deep	
\$00C841630-00_0	30"L
\$00C841636-00_0	36"L
\$00C841642-00_0	42"L
\$00C841648-00_0	48"L
22" Deep	
S00C842230-00_0	30"L
S00C842236-00_0	36"L
S00C842242-00_0	42"L
S00C842248-00_0	48"L


Sliding Glazed Doors



12¾" Deep		
S41C841330-00	30"L	
S41C841336-00	36"L	
S41C841342-00	42"L	
S41C841348-00	48"L	


16" Deep		
S41C841630-00 S41C841636-0 0	30"L 36"L	
S41C841642-00 S41C841648-00	42"L 48"L	
22" Deep		
S41C842230-00	30"L	
S41C842236-00	36"L	
S41C842242-00	42"L	
S41C842248-0 0	48"l	

Sliding Panel Doors



12¾" Deep	
S12C841330-00	30"L
S12C841336-00	36"L
S12C841342-00	42"L
S12C841348-00	48"L
16" Deep	
S12C841630-00	30"L
S12C841636-00	36"L
S12C841642-00	42"L
S12C841648-00	48"L
22" Deep	
S12C842230-00	30"L
S12C842236-00	36"L
S12C842242-00	42"L
S12C842248-00	48"L





Specifications:

Full Height Storage Cabinets are 84" high and 123/4", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

CAUTION: Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

NOTE: Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

Swinging Glazed Doors


Left Hand	Right Hand	12 ³ / ₄ " D
S50C841318L		
S50C841324L	S50C841324	24 L
		16" D
S50C841618L		
S50C841624L	S50C841624	24"L
		22" D
S50C842218L	S50C8422-18_	18"L
S50C842224L	S50C8422-24_	24"L

]				1
]=	=======	ŀ	======	
J				

12¾" Deep	
S50C841330	30"L
S50C841336	36"L
S50C841342	42"L
S50C841348	48"L
16" Deep	
S50C841630	30"L
S50C841636	36"L
S50C841642	42"L
S50C841648	48"L
22" Deep	
S50C842230	30"L
S50C842236	36"L
S50C842242-	42"L
S50C842248-	48"L

Swinging Panel Doors

	Left Hand	Right Hand	12³/₄" D
:========	S25C841318L	S25C841318	18"L
	S25C841324L	S25C841324	24"L
			16" D
	S25C841618L	S25C841618	18"L
	S25C841624L	S25C841624	24"L
:========			22" D
	S25C842218L	S25C842218	18"L
	S25C842224L	S25C842224	24"L



12¾" Deep	
S25C841330	30"L
S25C841336	36"L
S25C841342	42"L
S25C841348	48"L
16" Deep	
S25C841630	30"L
S25C841636	36"L
S25C841642	42"L
S25C841648	48"L
22" Deep	
S25C842230	30"L
S25C842236	36"L
S25C842242	42"L
S25C842248	48"L



Specifications:


Microscope Storage Cabinets are 84" high and 16" or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

Full height swinging doors are are hung on three hinges per door.

CAUTION: Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

NOTE: Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

Microscope Storage Cabinets



16" Deep S51C841624L____ 24"L

22" Deep S51C842224L____ 24"L

16" Deep

S51C841624-___ 24"| **22" Deep**

S51C842224-____ 24"L

16" Deep

S51C841636-____ 36"L **22" Deep**

S51C842236-____ 36"L

16" Deep

S51C841648-____ 48"L **22" Deep**

S51C842248-____ 48"L

Provides 15 storage compartments

Provides 15 storage compartments

Provides 20 storage compartments

Provides 30 storage compartments





Specifications:

Full Height Storage Cabinets are 84" high and 12¾", 16", or 22" deep. Lengths as shown. Split door cabinets are furnished with a fixed center shelf and four adjustable shelves.

Wardrobe Storage Cabinets are 84" high and 22" deep. Lengths as shown. They are provided with a hat shelf and a horizontal bar for hanging garments.


42"L

48"L

CAUTION: Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

NOTE: Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

Split Swinging Glazed Doors




12¾" Deep	
S60C841330	30"L
S60C841336	36"L
S60C841342	42"L
S60C841348	48"L
16" Deep	
S60C841630	30"L
S60C841636	36"L

22" Deep	
S60C842230 S60C842236 S60C842242 S60C842248	30"L 36"L 42"L 48"L

S60C841642-___

S60C841648-___

Split Swinging Panel Doors



S35C841330 S35C841336 S35C841342 S35C841348 16" Deep	30"L 36"L 42"L 48"L
S35C841630 S35C841636 S35C841642 S35C841648	30"L 36"L 42"L 48"L
22" Deep S35C842230 S35C842236 S35C842242 S35C842248	30"L 36"L 42"L 48"L

123/4" Deep

Wardrobe Storage Cabinets



Left Hand — 22"	Deep
S40C842218L	18"L
S40C842224L	24"L

Right Hand — 22	' Deep
S40C842218	18"L
S40C842224-	24"L



22" Deep		
S40C842236- 36"L		
S40C842248-	48"L	

Research Collection

Steel Full Height Storage Cabinets



Specifications:

Full Height Storage Cabinets are 84" high. Lengths as shown.


Acid Storage Cabinets ar furnished with a polyethylene corrosion resistant lining and five full-width corrosion resistant shelves. Ventilation is provided by louvers at the bottom of the door and a 31/2" O.D. vent outlet at the top of the cabinet.

Full height swinging doors are are hung on three hinges per door.

CAUTION: Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.


NOTE: Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

Acid Storage Cabinets



Left Hand — 16" D	еер
S70C841624L	24"L
Left Hand — 22" D	еер
S70C842224L	24"L






16" Deep		
S70C841636	36"L	
S70C841648	48"L	
22" Deep		
S70C842236	36"L	
S70C842248	48"L	



Base Cabinet Scribes & Fillers

End Scribes



To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe and wall angle.

Available in widths from 1"- 24" in one inch increments

Sitting Height — 29" High

BESC2900 ___-00_0


ADA Height — 32" High

BESC3200 ...-00_0

Standing Height — 35" High

BESC3500 ___-00_0

End Scribe Applique Panels




End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

End Scribe must be ordered separately

Steel Square Edge Widths from 3"- 24" in one inch increments	Steel Beveled Edge Widths from 3" - 24" in one inch increments	Wood Square Edge Widths from 3" - 24" in one inch increments
Sitting Height — 29" High		
BEAC2900∐-11_0	BEAC2900∐-1700	BEAC290016_0
ADA Height — 32" High — Steel		
ADA Height — 32" High — St	eel	
ADA Height — 32" High — St BEAC3200□-11_0	eel BEAC3200∐-1700	BEAC3200□-16_0
	BEAC3200∐-1700	BEAC3200□-16_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel. Example: BESC3500**06**-0000 for 6" wide End Scribe or BEAC3500**06**-1100 for 6" wide Steel Square Edge Applique Panel

Intermediate Fillers



To close opening between two adjacent base cabinets. Includes detached toe space filler.

Available in widths from 6"- 48" in one inch increments

Sitting Height — 29" High

BIFC2900 ___-00_0


ADA Height — 32" High

BIFC3200 ...-00_0

Standing Height — 35" High

BIFC3500Ⅲ-00_0

Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.


Intermediate Filler must be ordered separately

Steel Square Edge Widths from 6"- 48" in one inch increments	Steel Beveled Edge Widths from 6" - 48" in one inch increments	Wood Square Edge Widths from 6" - 48" in one inch increments
Sitting Height — 29" Hig	h	
BIAC2900	BIAC2900 □-1700	BIAC2900∐-16_0
ADA Height — 32" High		
BIAC3200 □-11_0	BIAC3200 □ -1700	BIAC3200Ⅲ-16_0
Standing Height — 35" H	ligh	

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel. Example: BIFC3500**12**-0000 for 12" wide Intermediate Filler or BIAC3500**12**-1100 for 12" wide Steel Square Edge Applique Panel

Base Cabinet Scribes & Fillers

Corner Fillers



To close opening between two adjacent base cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1"- 12" in one inch increments

Sitting Height — 29" High

BCFC2900 □ -00_0


ADA Height — 32" High

BCFC3200 ...-00_0

Standing Height — 35" High

BCFC3500 □ -00_0

Corner Applique Panels




Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

Corner Filler must be ordered separately

Steel Square Edge Widths from 3"- 12" in one inch increments	Steel Beveled Edge Widths from 3"- 12" in one inch increments	Wood Square Edge Widths from 3" - 12" in one inch incrementss
Sitting Height — 29" High	h	
BCAC2900	BCAC2900	BCAC2900 □ -16_0
ADA Height — 32" High		
ADA Height — 32" High BCAC3200□-11_0	BCAC3200⊞-1700	BCAC3200∐-16_0
		BCAC3200⊞-16_0

Insert width required before dash in part number to specify size of Corner Filler or Applique Panel. Example: BCFC3500**04**-0000 for 4" Corner Filler or BCAC3500**04**-1100 for 4" Steel Square Edge Applique Panel

Rear Scribes



To close opening between wall and rear of adjacent base cabinet. Includes wall mounting angle.

Widths from 1"- 24" in one inch increments

Sitting Height — 29" High


BRSC2900 ...-00_0

ADA Height — 32" High BRSC3200□-00_0

Standing Height — 35" High

BRSC3500 -00_0

Pipe Chase Fillers



To close opening between two adjacent base cabinets installed back-to-back or a cabinet and a finished back.

Widths from 1"- 24" in one inch increments

Sitting Height — 29" High

BPCC2900 ...-00_0

ADA Height — 32" High


BPCC3200 ...-00_0

Standing Height — 35" High

BPCC3500 ...-00_0

Insert width required before dash in part number to specify width of Rear Scribe or Pipe Chase Filler. Example: BRSC3500**8**-0000 for 8" wide Rear Scribe or BPCC3500**8**-0000 for 8" wide Pipe Chase Filler.

Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of base cabinets to cabinet faces.

Wood Square Edge

Sitting Height — 29" High		
BPAC290023-16_0 BPAC290029-16_0	23" W 29" W	
ADA Height — 32" High		
BPAC320023-16_0	23" W	
BPAC320029-16_0	29" W	
Standing Height — 35" High		


BPAC350023-16_0 23" W **BPAC350029-16_0** 29" W



Wall Cabinet Scribes & Fillers

End Scribes

End Scribe Applique Panels



To close opening between wall and front of adjacent wall cabinet. Includes wall mounting angles.

12¾" Deep	16" Deep
Available in widths fro	m 1"- 24" in one inch increments


24" High

WESC2413 - 00_0 WESC2416 - 00_0
30" High

WESC3013 - 00_0 WESC3016 - 00_0
36" High

WESC3613 - 00_0 WESC3616 - 00_0
48" High

WESC4813 - 00_0 WESC4816 - 00_0



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

End Scribe must be ordered separately

Steel Square Edge Widths from 3"- 24" in one inch increments	Steel Beveled Edge	Wood Square Edge Widths from 3"- 24" in one inch increments
24" High		
WEAC2400□-11_0	not available	WEAC2400∐-16_0
30" High		
WEAC3000□-11_0	not available	WEAC3000□-16_0
36" High		
WEAC3600□-11_0	not available	WEAC3600□-16_0
48" High		
WEAC4800 □ -11_0	not available	WEAC4800□-16_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel. Example: WESC3613**06**-0000 for 6" wide End Scribe or WEAC3600**06**-1100 for 6" wide Steel Square Edge Applique Panel


Intermediate Fillers

D. W

To close opening between two adjacent wall cabinets.

16" Deep

Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

Intermediate Filler must be ordered separately

Available in widths from 6	- 48" in one inch increments
24" High	
WIFC2413Ⅲ-00_0	WIFC2416Ⅲ-00_0
30" High	
WIFC3013Ⅲ-00_0	WIFC3016Ⅲ-00_0
36" High	
WIFC3613Ⅲ-00_0	WIFC3616Ⅲ-00_0
48" High	
WIFC4813Ⅲ-00_0	WIFC4816□-00_0

Steel Square Edge Widths from 6"- 48" in one inch increments	Steel Beveled Edge	Wood Square Edge Widths from 6"- 48" in one inch increments
24" High		
WIAC2400 □ -11_0	not available	WIAC2400 □ -16_0
30" High		
WIAC3000 □-11_0	not available	WIAC3000 □-16_0
36" High		
WIAC3600 □-11_0	not available	WIAC3600 □-16_0
48" High		
WIAC4800 □ -11_0	not available	WIAC4800 □ -16_0


Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel. Example: WIFC361312-0000 for 12" wide Intermediate Filler or WIAC360012-1100 for 12" wide Steel Square Edge Applique Panel

123/4" Deep


Wall Cabinet Scribes & Fillers

Corner Fillers

Corner Applique Panels



To close opening between two adjacent wall cabinets joined at right angles.



Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.


Corner Filler must be ordered separately

12¾" Deep Available in depth and widincrements.	16" Deep ths from 1"- 12" in one inch
24" High	
WCFC2413Ⅲ-00_0	WCFC2416Ⅲ-00_0
30" High	
WCFC3013Ⅲ-00_0	WCFC3016Ⅲ-00_0
36" High	
WCFC3613Ⅲ-00_0	WCFC3616Ⅲ-00_0
48" High	
WCFC4813∐-00_0	WCFC4816□-00_0

Steel Square Edge Widths from 3"- 12" in one inch increments	Steel Beveled Edge	Wood Square Edge Widths from 3"- 12" in one inch increments
24" High		
WCAC2400∐-11_0	not available	WCAC2400Ⅲ-16_0
30" High		
WCAC3000	not available	WCAC3000∐-16_0
36" High		
WCAC3600⊞-11_0	not available	WCAC3600Ⅲ-16_0
48" High		
WCAC4800⊞-11_0	not available	WCAC4800 □-16_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel. Example: WCFC3613**06**-0000 for 6" wide Corner Filler or WCAC3600**06**-1100 for 6" wide Steel Square Edge Applique Panel

Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of wall cabinets to cabinet faces.


Wood Square Edge

24" High	
WPAC240013-16_0	12 ³ /4" W
WPAC240016-16_0	16" W
30" High	
WPAC300013-16_0	12 ³ /4" W
WPAC300016-16_0	16" W
36" High	
WPAC360013-16_0	12 ³ / ₄ " W
WPAC360016-16_0	16" W
48" High	
WPAC480013-16_0	12 ³ / ₄ " W
WPAC480016-16_0	16" W



Full Height Cabinet Scribes & Fillers

Full Height End Scribes



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe and wall mounting angles.

Available in widths from 1"- 24" in one inch increments

123/4" Deep

FESC8413 ___-00_0

16" Deep

FESC8416Ⅲ-00_0

22" Deep

FESC8422-00_0

Full Height End Applique Panel



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.


End Scribe must be ordered separately

Steel Square Edge	Steel Beveled Edge	Wood Square Edge
Widths from 3"- 24"		Widths from 3"- 24"
in one inch increments		in one inch increments

FEAC8400 ☐ -11_0 not available FEAC8400 ☐ -16_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel. Example: FESC8413**06**-0000 for 6" wide End Scribe or FEAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

Full Height Intermediate Fillers



To close opening between two adjacent full height cabinets. Includes detached toe space filler.

Available in widths from 6"- 48" in one inch increments

123/4" Deep

FIFC8413 ___-00_0


16" Deep

FIFC8416**□**-00_0

22" Deep

FIFC8422 -00 0

Full Height Intermediate Applique Panel



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

Intermediate Fillers must be ordered separately

Steel Square EdgeSteel Beveled EdgeWood Square EdgeWidths from 6"- 48"Widths from 6"- 48"in one inch incrementsin one inch increments


FIAC8400 ☐ -11_0 not available FIAC8400 ☐ -16_0

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel. Example: FISC8413**06**-0000 for 6" wide Intermediate Filler or FIAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

Research Collection

Full Height Cabinet Scribes & Fillers

Corner Filler



To close opening between two adjacent full height cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1"- 12" in one inch increments.

123/4" Deep

FCFC8413 ___-00_0

16" Deep

FCFC8416-00_0

22" Deep

FCFC8422-00_0

Corner Fillers Applique Panel




Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

Corner Filler must be ordered separately

Steel Square Edge Widths from 3"- 12" in one inch increments	Steel Beveled Edge	Wood Square Edge Widths from 3'- 12' in one inch increments
FCAC8400∐-11_0	not available	FCAC8400 □-16_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel. Example: FCFC8413**06**-0000 for 6" wide Corner Filler or FCAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

Full Height Rear Scribe




To close opening between wall and rear of adjacent full height cabinets.

Available in widths from 1"- 24" in one inch increments

56" High 84" High FRSC5600 □-00_0 FRSC8400 □-00_0

Insert width required before dash in part number to specify width of Rear Scribel. Example: FRSC8400**6**-0000 for 6" wide Rear Scribe.

Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of full height cabinets to cabinet faces.

Wood Square Edge

84" High	
FPAC840013-16_0	12 ³ / ₄ " W
FPAC840016-16_0	16" W
FPAC840022-16_0	22" W



Cabinet Finished Backs

Base Cabinet Finished Back

Wall Cabinet Finished Back

Full Height Finished Back





To cover backs of base cabinets exposed to view. Back adds 1" to the depth of the base cabinet.

42"W

48"W

60"W

Sitting Height — 29" High		
BFBC290112-00_0	12"W	
BFBC290115-00_0	18"W	
BFBC290118-00_0	18"W	
BFBC290121-00_0	21"W	
BFBC290124-00_0	24"W	
BFBC290130-00_0	30"W	
BFBC290136-00_0	36"W	

ADA Height — 32" High

BFBC290142-00_0

BFBC290148-00_0


BFBC290160-00_0

BFBC320112-00_0	12"W
BFBC320115-00_0	15"W
BFBC320118-00_0	18"W
BFBC320121-00_0	21"W
BFBC320124-00_0	24"W
BFBC320130-00_0	30"W
BFBC320136-00_0	36"W
BFBC320142-00_0	42"W
BFBC320148-00_0	48"W
BFBC320160-00_0	60"W

Standing Height — 35" High

Standing neight — 35	підп
BFBC350112-00_0	12"W
BFBC350115-00_0	15"W
BFBC350118-00_0	18"W
BFBC350121-00_0	21"W
BFBC350124-00_0	24"W
BFBC350130-00_0	30"W
BFBC350136-00_0	36"W
BFBC350142-00_0	42"W
BFBC350148-00_0	48"W
BFBC350152-00_0U3*	52"W
BFBC350158-00_0U3*	58"W
BFBC350160-00_0U3*	60"W
BFBC350160-00_0	60"W

^{*} Includes toespace notches for End Sink Cabinets



To cover backs of wall cabinets exposed to view. Adds 1" to depth of cabinet.

24" High

24 nigii	
WFBC240115-00_0	15"W
WFBC240118-00_0	18"W
WFBC240121-00_0	21"W
WFBC240124-00_0	24"W
WFBC240130-00_0	30"W
WFBC240136-00_0	36"W
WFBC240142-00_0	42"W
WFBC240148-00_0	48"W

30" High

15"W
18"W
21"W
24"W
30"W
36"W
42"W
48"W
15"W
18"W

21"W

24"W 30"W

36"W 42"W

48"W

48" High


WFBC360121-00_0

WFBC360124-00_0

WFBC360130-00_0 WFBC360136-00_0

WFBC360142-00 0 WFBC360148-00_0

WFBC480115-00_0 15"W WFBC480118-00_0 18"W WFBC480121-00_0 21"W WFBC480124-00_0 24"W WFBC480130-00_0 30"W WFBC480136-00_0 36"W WFBC480142-00_0 42"W WFBC480148-00_0 48"W		
WFBC480121-00_0 21"W WFBC480124-00_0 24"W WFBC480130-00_0 30"W WFBC480136-00_0 36"W WFBC480142-00_0 42"W	WFBC480115-00_0	15"W
WFBC480124-00_0 24"W WFBC480130-00_0 30"W WFBC480136-00_0 36"W WFBC480142-00_0 42"W	WFBC480118-00_0	18"W
WFBC480130-00_0 30"W WFBC480136-00_0 36"W WFBC480142-00_0 42"W	WFBC480121-00_0	21"W
WFBC480136-00_0 36"W WFBC480142-00_0 42"W	WFBC480124-00_0	24"W
WFBC480142-00_0 42"W	WFBC480130-00_0	30"W
	WFBC480136-00_0	36"W
WFBC480148-00_0 48"W	WFBC480142-00_0	42"W
	WFBC480148-00_0	48"W




To cover backs of full height cabinets exposed to view. Adds 1" to depth of cabinet.

84" High

FFBC840118-00_0	18"W
FFBC840124-00_0	24"W
FFBC840130-00_0	30"W
FFBC840136-00_0	36"W
FFBC840142-00_0	42"W
FFBC840148-00_0	48"W

Steel Adjustable Slip Scribes



Specifications:


Steel Adjustable Slip Scribes are designed the same as standard scribes with the addition of a Slip Scribe Angle at the wall edge to dress-up the edge of the scribe.

The Steel Adjustable Slip Scribe is designed to be cut to width in the field to more easily fit the space between the wall and the end of the furniture assembly.

Adjustable Base Cabinet Scribe


Adjustable Wall Cabinet Scribe

Adjustable Full Height Scribe




To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe, end scribe and slip scribe angle.

Sitting Height — 29" High			
ABSC290012-00_0	12"W		
ADA Height — 32" High			
ABSC320012-00_0	12"W		
Standing Height — 35" High			
ABSC350012-00_0	12"W		



To close opening between wall and front of adjacent wall cabinet. Includes end scribe and slip scribe angle.

24" High		
AWSC241312-00_0	12¾" D	12"W
AWSC241612-00_0	16" D	12"W
30" High		
AWSC301312-00_0	12¾" D	12"W
AWSC301612-00_0	16" D	12"W
36" High		
AWSC361312-00_0	12¾" D	12"W
AWSC361612-00_0	16" D	12"W
48" High		
AWSC481312-00_0	12¾" D	12"W
AWSC481612-00 0	116" D	12"W



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe, end scribe, and slip scribe angle.

84" High		
AFSC841312-00_0	12¾" D	12"W
AFSC841612-00_0	16" D	12"W
AFSC842212-00_0	22" D	12"W

Slip Scribe Angle




May be used with any standard scribe or filler to make it adjustable. Available in lengths from 3" to 96" in one inch increments.

SCAC0101 -00_0

Insert length required before dash in part number to specify length of angle.

Spring Clips



18 guage spring steel clip, $^{3}/_{4}$ " x $^{1}/_{2}$ ". May be used with any scribe or filler to make them easily removable.

F-4995-PK

package of 10

13³/₄" Hiah



Steel Sloping Tops

Sloping Tops


Sloping Top Fillers

83/8" High



May be used on wall cabinets and full height cabinets to keep the tops dust free.

83/8" High	10 ³ / ₄ " High	13 ³ / ₄ " High	
12 ³ / ₄ " Deep	16" Deep	22" Deep	
SLTC081315-00_0	SLTC111615-00_0	SLTC142215-00_0	15"L
SLTC081318-00_0	SLTC111618-00_0	SLTC142218-00_0	18"L
SLTC081321-00_0	SLTC111621-00_0	SLTC142221-00_0	21"L
SLTC081324-00_0	SLTC111624-00_0	SLTC142224-00_0	24"L
SLTC081330-00_0	SLTC111630-00_0	SLTC142230-00_0	30"L
SLTC081336-00_0	SLTC111636-00_0	SLTC142236-00_0	36"L
SLTC081342-00_0	SLTC111642-00_0	SLTC142242-00_0	42"L
SLTC081348-00_0	SLTC111648-00_0	SLTC142248-00_0	48"L
SLTC081360-00_0	SLTC111660-00_0	SLTC142260-00_0	60"L
SLTC081372-00_0	SLTC111672-00_0	SLTC142272-00_0	72"L
SLTC081384-00_0	SLTC111684-00_0	SLTC142284-00_0	84"L
SLTC081396-00_0	SLTC111696-00_0	SLTC142296-00_0	96"L




To close opening between wall and front of adjacent sloping top. Includes wall mounting angle.

10³/₄" Hiah

12 ³ / ₄ " Deep	16" Deep	22" Deep	
STFC081301-00_0	STFC111601-00_0	STFC142201-00_0	1"W
STFC081306-00_0	STFC111606-00_0	STFC142206-00_0	6"W
STFC081312-00_0	STFC111612-00_0	STFC142212-00_0	12"W

Corner Cabinet Sloping Tops




To be used on corner wall cabinets.

83/8" High	10 ³ / ₄ " High	
12 ³ / ₄ " Deep	16" Deep	

CSTC081327-00_0 CSTC111627-00_0 267/8"L

Sloping Top Corner Fillers



May be used on wall cabinets and full height cabinets when cabinets are joined at right angles.


8%" High	10¾" High	13¾" High
12¾" Deep	16" Deep	22" Deep


SCFC081301-00_0 SCFC111601-00_0 SCFC142201-00_0 1"W SCFC081306-00_0 SCFC111606-00_0 SCFC142206-00_0 6"W

Steel Ceiling & Pipe Enclosures

Ceiling Enclosures

Corner Ceiling Enclosures





To close opening between ceiling and top of wall or full height cabinets.

To close opening between ceiling and top of corner wall cabinets.

Front Ceiling Enclosure

12" High	24" High	
FCEC120018-00_0	FCEC240018-00_0	18"L
FCEC120024-00_0	FCEC240024-00_0	24"L
FCEC120030-00_0	FCEC240030-00_0	30"L
FCEC120036-00_0	FCEC240036-00_0	36"L
FCEC120042-00_0	FCEC240042-00_0	42"L
FCEC120048-00_0	FCEC240048-00_0	48"L
FCEC120060-00_0	FCEC240060-00_0	60"L
FCEC120072-00_0	FCEC240072-00_0	72"L
FCEC120084-00_0	FCEC240084-00_0	84"L
FCEC120096-00_0	FCEC240096-00_0	96"L

End Ceiling Enclosure

U01C


96=96"

24" High	
ECEC240013-00_0	12 ³ / ₄ "D
ECEC240016-00_0	16"D
ECEC240022-00_0	22"D
	ECEC240013-00_0 ECEC240016-00_0

12" High	24" High	Depth	Length
CCEC121320-00_0	CCEC241320-00_0	12 ³ / ₄ "	20"
CCEC121615-00 0	CCEC241615-00 0	16"	15 ³ / ₈ "

Pipe Enclosures





Pipe Enclosures conceal vertical service lines coming down from the ceiling to the worksurface. One side is removable for easy access with no exposed fasteners. Designed to sit on top of curb on worksurface and penetrate drop ceiling. Available in powder coated steel or stainless steel.


-00_0 Wall Bench Enclosure **U02C** -00_0 Center Bench Enclosure Available Available Available Available Materials Depths Heights Lengths **60**=60" **05**=5" 12 = 12" 0 = Powder Coated Steel**66**=66" **06**=6" **18**= 18" **3** = Stainless Steel **72**=72" **07**=7" **24**=24" **78**=78" **08**=8" **30**=30" **84**=84" **09**=9"

36=36"

48=48"

90=90" **10**=10" **42**=42"

Center Bench Enclosure





Fume Hoods & Ventilated Enclosures

Ducted Ventilated Products

Kewaunee offers a full line of Laboratory Fume Hoods to meet the needs of every lab application. In addition to our newest line of **Supreme Air Venturi Fume Hoods**, we also offer uniquely designed

hoods for education and specialty hoods for custom applications.

Supreme Air – General Purpose Fume Hoods

The design is based on rigid frame construction that assures a solid installation and low vibration and sound levels. Access panels are easily removable to access service lines and fittings. Radius corner posts and airfoils, plus, easy operator control of interior baffle settings assure a high level of comfort, safety and efficiency.


Supreme Air Venturi – Next Generation Fume Hoods

The Supreme Air Venturi utilizes proven containment enhancing techniques coupled with innovative laboratory convenience features, including: adjustable LED lights, electromechanical sash stop, ultra-smooth sash system, and push-button sash management. The Supreme Air Venturi provides the safest user friendly fume hood manufactured today.

TruView – Teaching Fume Hoods

Designed for today's classroom allowing clear line-of-sight through multiple hoods to allow easier supervision of sudents by the instructor. May be used as a standalone hood or in assemblies, both single-sided for use against the wall, or double-sided for use in the center of the laboratory.







For more information on Kewaunee's **Supreme Air Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirHoods/

For more information on Kewaunee's **Supreme Air Venturi Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirVenturi/

For more information on Kewaunee's **TruView Teaching Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/FumeHoods/TruViewHood/

Research Collection

Fume Hoods & Ventilated Enclosures

Ductless Ventilated Products

Kewaunee offers a full line of filtered products to meet every protection requirement. **DetectAir** ductless fume hoods provide a safe alternative to conventional chemical fume hoods

DetectAir – Ductless/Filtered Fume Hoods

With the introduction of DetectAir, the filtered fume hood has come of age. DetectAir combines Kewaunee's traditional fume hood expertise with the filtering technology of AirClean Systems.

when total exhaust is not an option. The Interceptor Biological Safety Cabinets are offered in both class II A2 and B2 models, (NSF Certified) for both recirculating and direct exhaust


the Interceptor – Biological Safety Cabinet

The Interceptor Class II, A2 and Class II, B2 Biological Safety Cabinets are Kewaunee's latest design. Developed with comprehensive features created to enhance both safety and ergonomics... the Interceptor leads its class.


requirements. The **Interceptor VLF** provides ISO 14644-1 Standard Class 5 conditions, providing a sterile environment for product protection.

the Interceptor – Vertical Laminar Flow Workstation

Interceptor Vertical Laminar Flow Workstations feature a front loading HEPA filter, hinged front panel with self-support brackets, 10° angled safety glass, Minihelic gauge and T5 florescent lamp.



For more information on Kewaunee's **DetectAir - Ductless/Filtered Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/DuctlessHoods/



For more information on the Interceptor Biological Safety Cabinet from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/BioSafetyCabinets/



For more information on the Interceptor - Vertical Laminar Flow Workstation from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: http://www.kewaunee.com/LaboratoryHoods/LaminarFlowHoods/






Work Tops

Work Top material selections are among the most critical of the laboratory furniture decisions. The intended applications for your work tops dictate material choice. Surfaces could be exposed to reagents and caustic substances. Heavy equipment may need to be supported or only light duty tasks may be performed. Kewaunee can recommend and furnish materials fabricated to your individual needs, complete with appropriate sinks, drains and accessories.

Kemresin® Work Tops

manufactured exclusively by Kewaunee



Appearance: Ultra-smooth non-porous surface.

Available Colors: Black, Grey, Putty, and Slate

Composition: 1" thick molded epoxy resin.

Physical Properties: Excellent resistance to practically all acids,

alkalies and solvents. Highly resistant to

heat and normal physical abuse.

Applications: Use wherever optimum physical and

chemical resisting properties are required.

Other Work Top Materials

			OUENION	
MATERIAL	APPEARANCE	COMPOSITION	CHEMICAL RESISTANCE	APPLICATIONS
STAINLESS STEEL	No.4 Finish	Type 304 furnished as standard, but other types are available when specified. Type 316 supplied for Perchloric Acid fume hoods.	Good chemical resistance to a wide range of reagents.	Recommended for sterile and radiological laboratories.
HARDWOOD, OAK OR MAPLE	Natural varnish, or penetrating oil finish	Solid edge grain Maple or Oak strips, glued together electron- ically, and finished with baked- on clear synthetic varnish or penetrating oil.	Fair chemical resistance, but not recommended for use in chemical laboratories.	Excellent top material for Physics tables, Electronic benches, and other tables where mechanical or physical work levels are prevalent.
COMPOSITE RESIN	Hard, smooth non-glare surface	Cellulose fiber reinforced phenolic resin core with a polyurethane copolymer surface. Available in multiple colors and patterns.	Good chemical resistance to wide range of acids, alkalies and solvents.	Use wherever a high degree of chemical resistance is required. (Not recommended for fume hood worksurfaces.)
HIGH-PRESSURE LAMINATES	Hard, smooth surface	Plastic laminated surface attached with water resistant adhesive to solid particle board. Available in choice of colors and patterns.	Fair chemical resistance to wide range of acids, alkalies and solvents.	Recommended for light duty purposes.

Specifications:


The configurations shown are for Kewaunee base cabinets . Standard top thickness is 1", standard curb height is 4".


To order, select configuration and top material desired. Specify number of linear or square feet per product number.

Dimensions are nominal. Not all con-


figurations are available in all materials. Indented tops available on special

Wall Bench Tops






31"



Standard	18"	
Sizes:	23"	
	24"	
	30"	
	31"	
-	— в ——∥	C

Top & Ledge with Applied curb Style			
Standard	A 30"	B 23"	C 7"
Sizes:	31"	23"	8"



Corner Top			
no curbs .		Style	
	Α	В	С
Standard	33"	24"	9"
Sizes:	39"	30"	9"
	40"	31"	9"

Corner Top with Applied Curbs Style L				
	Α	В	С	
Standard	33"	24"	9"	
Sizes:	39"	30"	9"	
	40"	31"	9"	


Center Bench Tops

Standard Sizes:




23"

23"





Top with Center Joint			Style I	
	Α	В	С	
Standard	54"	27"	27"	
Sizes:	54"	23"	31"	
	60"	30"	30"	
	62"	31"	31"	



Top with Center Ledge			Style J	
	Α	В	С	D
Standard	54"	23"	8"	23"
Sizes:	60"	23"	14"	23"
	62"	23"	16"	23"

Ledges and Shelves



Ledge Unfinishe	ed Edges	Style A
Standard	8"	
Sizes:	10"	
	12"	
	14"	
	16"	



Shelf Finished	Edges	Style N
Standard	8"	
Sizes:	10"	
	12"	
	14"	
	16"	



Center Bench Shelving Assemblies

Specifications:


Shelving Support Assemblies are designed to add shelving to center benches. Each assembly consists of four uprights with slots on 1" centers and a top frame. They are 12" deep in lengths as shown.

Each assembly is furnished with mounting brackets for attaching to the counter top and cabinets below but no fasteners. Fasteners must be furnished in field depending on top material and cabinet configuration.

Shelving Assembly Support Struts are especially designed to provide added support to the Shelving Support Assemblies while allowing space for connection of plumbing or electrical connections.

Shelving Assembly Support Strut

12" Deep Shelving Support Assembly




Length	Part Number
24"	F17-4824-00
30"	F17-4830-00
36"	F17-4836-00
42"	F17-4842-00
48"	F17-4848-00
54"	F17-4854-00
60"	F17-4860-00
66"	F17-4866-00
72"	F17-4872-00



Width	Part Number
12"	K17-3612-00
1//"	K17-3614-00

Shelving Uprights



Powder Coated Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	F13M000224-00	F12M000224-00	each
30"	F13M000230-00	F12M000230-00	each
36"	F13M000236-00	F12M000236-00	each
42"	F13M000242-00	F12M000242-00	each
48"	F13M000248-00	F12M000248-00	each
72"	F13M000272-00	F12M000272-00	each
79"	F13M000279-00	F12M000279-00	each
84"	F13M000284-00	F12M000284-00	each

Stainless Steel Uprights

	. •		
Length	Single Slot Upright	Double Slot Uprigh	nt
24"	F13S000224-00	F12S000224-00	each
30"	F13S000230-00	F12S000230-00	each
36"	F13S000236-00	F12S000236-00	each
42"	F13S000242-00	F12S000242-00	each
48"	F13S000248-00	F12S000248-00	each
72"	F13S000272-00	F12S000272-00	each
79"	F13S000279-00	F12S000279-00	each
84"	F13S000284-00	F12S000284-00	each

Designed to be used with Adjustable Module Shelving Assemblies to create wall mounted shelving.

Center Bench Shelving Assemblies

Specifications:

Pass-Thru Shelf Assemblies for 12" Island Modules are six inches deeper than the Shelf Bracket to allow the Shelf to pass through the module core and form a continuous surface. These shelves may also be moved forward and mounted to form


a 2" cord gap between the shelves. Shelf Brackets are 11 gauge steel, adjustable on 1" centers.

Painted Steel Shelves – 1" thick, formed from 16 gauge cold rolled steel


Veneer Plywood Shelves – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

Phenolic Resin Shelves – 1" thick phenolic resin.

Pass-Thru Shelf Assemblies



Pass-Thru Shelves



Pass-Thru Shelves with Cord Gap

L	anath	Daintad Staal	Vancar Physicad	Phonolia Pooin
	ength	Painted Steel	Veneer Plywood	Phenolic Resin
4	24"	F21M0924-60	F21W0924-6_	F21T0924-60
# #	30"	F21M0930-60	F21W0930-6_	F21T0930-60
SHELF acket	36"	F21M0936-60	F21W0936-6_	F21T0936-60
_	42"	F21M0942-60	F21W0942-6	F21T0942-60
	48"	F21M0948-60	F21W0948-6	F21T0948-60
DEEP Deep B	54"	F21M0954-60	_	
	60"	F21M0960-60		
ص ال	66"			
-	72"	F21M0972-60		
	24"	F21M1224-60	F21W1224-6_	F21T1224-60
SHELF tracket	30"	F21M1230-60	F21W1230-6_	F21T1230-60
SHEI Bracket	36"		F21W1236-6_	
S	42"	F21M1242-60	F21W1242-6_	F21T1242-60
	48"	F21M1248-60	F21W1248-6_	
" DEEP 2" Deep E	54"		1211112400_	12111240 00
	60"	F21M1260-60		
2 12				
8	66"			
	72"	F21M1272-60		
ш	24"	F21M1824-60	F21W1824-6_	
et 🔚	30"	F21M1830-60	F21W1830-6_	F21T1830-60
SHELF Bracket	36"	F21M1836-60	F21W1836-6_	F21T1836-60
	42"	F21M1842-60	F21W1842-6_	F21T1842-60
6 8	48"	F21M1848-60	F21W1848-6	
Jeep E	54"	F21M1854-60	_	


To specify Veener Plywood Species replace blank at end of part number with number that indicates material required.

Veener Plywood Material Options

-1 = Red Oak	-5 = FSC® Red Oak
-2 = White Maple	-6 = FSC® White Maple
-7 = Steamed Beech	-8 = FSC® St'md Beech

Shelf Retaining Lips

60"



F21M1860-60

Length	Shelf Lip	
24"	F63M0224-00	each
30"	F63M0230-00	each
36"	F63M0236-00	each
42"	F63M0242-00	each
48"	F63M0248-00	each
54"	F63M0254-00	each
60"	F63M0260-00	each
66"	F63M0266-00	each
72"	F63M0272-00	each

Designed to be added to Pass-thru and Wall Shelf Assemblies. Steel Shelf Lips may be mounted to steel, veneer plywood, or phenolic resin shelves with supplied screws.



Wall Shelving Assemblies

Specifications:

Vertical Uprights are 11/4" deep and either 79" or 42" high. They are manufactured with a double-slotted face with slots on 1" centers. They are designed to be mounted between two Support Rails to hold accurate spacing.


Support Rails are available in 8 foot lengths and may be cut to the required length. Support Rails are notched on 6" centers to accept the Vertical Uprights.

Support Rail Covers are required to finish the assembly and must be

ordered in pairs to match the Vertical Upright spacing.


Optional Center Rails may be ordered to provide a horizontal member at counter top height. They should be sized to match the Vertical Upright spacing.

Full Height Wall Suspension System



Vertical Uprigi	ht	
F10M0279-00	79"H	pair
Support Rail		
F11M0196-00	96"L	each
Support Rail C	overs	
F11M0324-00	24"L	pair
F11M0330-00	30"L	pair
F11M0336-00	36"L	pair
F11M0342-00	42"L	pair
F11M0348-00	48"L	pair
F11M0354-00	54"L	pair
F11M0360-00	60"L	pair
F11M0366-00	66"L	pair
F11M0372-00	72"L	pair
0.11		
Optional Cente		
F11M0224-00	24"L	each
F11M0230-00	30"L	each
F11M0236-00	36"L	each
F11M0242-00	42"L	each
F11M0248-00	48"L	each
F11M0254-00	54"L	each
F11M0260-00	60"L	each
F11M0266-00	66"L	each
F11M0272-00	72"L	each

Above Countertop Wall Suspension System



Vertical Uprigl	nt			
F10M0242-00	42"H	pair		
Support Rail				
F11M0196-00	96"L	each		
Support Rail C	overe			
F11M0324-00	24"L	pair		
F11M0330-00	30"L	pair		
F11M0336-00	36"L	pair		
F11M0342-00	42"L	pair		
F11M0348-00	48"L	pair		
F11M0354-00	54"L	pair		
F11M0360-00	60"L	pair		
F11M0366-00	66"L	pair		
F11M0372-00	72"L	pair		

Wall Shelving Assemblies

Specifications:

Adjustable Wall Shelf Assemblies consist of a pair of 11 gauge steel Shelf Brackets and a Shelf. The Shelf Bracket engages the inner slot of the double slotted Wall Vertical Upright and is adjustable on 1" centers. Shelves are available in Painted

Steel, Wood veneer plywood, or Phenolic Resin.

Painted Steel Shelves – 1" thick, formed from 16 gauge cold rolled steel

Daintad Ctaal

1 0000110


Veneer Plywood Shelves – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

Phenolic Resin Shelves – 1" thick phenolic resin.


Dhanalia Daain

Vanaar Dharraad

Adjustable Wall Shelf Assemblies



\geqslant			
	7 /	7	



Optional Shelf Locking Tab

Sliding tab mounted on Shelf Bracket locks Shelf into upright to keep it from being accidently removed.

To specify, add an **K** to the end of the Shelf Assembly part number.

Example: F21M1236-00,K

Length	Painted Steel	Veneer Plywood	Phenolic Resin
24" 30" 36" 9 DEEP SHELF 9 26" 9 26" 9 26" 60" 60" 66" 72"	F21M0624-00 F21M0630-00 F21M0636-00 F21M0642-00 F21M0648-00 F21M0654-00 F21M0660-00 F21M0666-00 F21M0672-00	F21W0624-0_ F21W0630-0_ F21W0636-0_ F21W0642-0_ F21W0648-0_	F21T0624-00 F21T0630-00 F21T0636-00 F21T0642-00 F21T0648-00
24" 30" 36" 42" 48" 54" 60" 66" 72"	F21M0924-00 F21M0930-00 F21M0936-00 F21M0942-00 F21M0948-00 F21M0954-00 F21M0960-00 F21M0966-00 F21M0972-00	F21W0924-0_ F21W0930-0_ F21W0936-0_ F21W0942-0_ F21W0948-0_	F21T0924-00 F21T0930-00 F21T0936-00 F21T0942-00 F21T0948-00
12" DEEP SHELF 12" Deep Bracket 12" Deep Bracket 12" Deep Bracket 12" 24" 26" 26" 26" 26" 26" 26" 26" 26" 26" 26	F21M1224-00 F21M1230-00 F21M1236-00 F21M1242-00 F21M1248-00 F21M1254-00 F21M1260-00 F21M1266-00 F21M1272-00	F21W1224-0_ F21W1230-0_ F21W1236-0_ F21W1242-0_ F21W1248-0_	F21T1224-00 F21T1230-00 F21T1236-00 F21T1242-00 F21T1248-00
18" DEEP SHELF 18" Deep Bracket 18" Deep Bracket 24" 26" 26" 26" 26" 26" 26" 26" 26" 26" 26	F21M1824-00 F21M1830-00 F21M1836-00 F21M1842-00 F21M1848-00 F21M1854-00 F21M1860-00	F21W1824-0_ F21W1830-0_ F21W1836-0_ F21W1842-0_ F21W1848-0_	F21T1824-00 F21T1830-00 F21T1836-00 F21T1842-00 F21T1848-00
24" DEEP SHELF 24" Deep Bracket 30" 36" 48" 54" 60" 60" 60" 60" 60" 60" 60" 60" 60" 60	F21M2424-00 F21M2430-00 F21M2436-00 F21M2442-00 F21M2448-00 F21M2454-00 F21M2460-00		F21T2424-00 F21T2430-00 F21T2436-00 F21T2442-00 F21T2448-00 wood Species replace umber with number that uired.

Veener Plywood Material Options

-7 = Steamed Beech -8 = FSC® St'md Beech

-1 = Red Oak

-2 = White Maple

-5 = FSC® Red Oak

-6 = FSC® White Maple



Reagent Rack & Shelving Systems

Reagent Rack Uprights are

avialable in two heights and two depths. All are 11/2" thick and available in either powder coated steel to match or complement the casework color or in type 304 stainless steel. 12" deep uprights are designed for


use on Center Benches with a 14" deep shelf, 10" deep uprights are designed for use on Wall Benches with a 12" deep shelf, but both can be used on either bench.

Installation hardware is included with each upright.

Spacing between upright should not exceed 48".

Plastic Laminate Shelves are 1" thick particle board faced on top and all four edges with black plastic laminate. A backer sheet covers the bottom.


Typical Reagent Rack Assemblies



Single Tier Center Bench Reagent Rack using 18" high by 12" deep uprights and 14"deep shelving.




Double Tier Center Bench Reagent Rack using 18" high by 12" deep and 12" by 12" deep uprights and 14"deep shelving.




Single Tier Wall Bench Reagent Rack using 18" high by 10" deep uprights and 12"deep shelving.

Reagent Rack Uprights



18" H x 12" D x 1½" T



12" High x 12" Deep

12" H x 12" D x 1½" T



18" H x 10" D x 1½" T



12" High x 10" Deep

12" H x 10" D x 1½"T


RRU-1812-0000 P.C. Steel RRU-1812-0030 St. Steel

RRU-1212-0000 P.C. Steel RRU-1212-0030 St. Steel

RRU-1810-0000 P.C. Steel RRU-1810-0030 St. Steel

RRU-1210-0000 P.C. Steel RRU-1210-0030 St. Steel

Plastic Laminate Shelves



PL-4814-00 PL-7214-00

14" Deep x 48" Long 14" Deep x 72" Long

PL-4812-00 PL-7212-00

12" Deep x 48" Long 12" Deep x 72" Long

12" Deep Shelf

Reagent Rack & Shelving Systems

Fixed Wall Mounted Shelf Brackets are designed to mount directly to the wall and furnished without mounting hardware.

Adjustable Wall Mounted Shelf Brackets are designed to mount to
KF-0170-04 Shelving Uprights and

include mounting bolts and nuts.


KF-0170-04 Shelving Uprights are designed to hold a maximum of three shelves.

Brackets and Uprights are finished in Kewaunee's unique corrosion

resistant powder coat finish to match or complement the casework color.

Typical Reagent Rack Assemblies






Fixed Wall Shelving



Adjustable Wall Shelving




Adjustable Wall Shelving

Flat Top Shelf Bracket

Lipped Shelf Bracket

Book End Bracket

Shelving Upright



Fixed Shelf Brackets

bracket size

M-0031-0A 11" × 8" 12"

M-0030-0A 15" × 10" 16"

M-0032-0A 17" × 12" 18"

part number is for single backet



Fixed Shelf Brackets bracket size Size M-0038-0A $12^1/8^{\circ} \times 6^{\circ}$ 12° M-0044-0A $16^1/8^{\circ} \times 6^{\circ}$ 16° M-0039-0A $18^1/8^{\circ} \times 6^{\circ}$ 18° part number is for single backet



 Fixed Shelf Brackets


 bracket size
 shelf size

 M-0060-0A
 12" x 6"
 12"

 M-0061-0A
 16" x 6"
 16"

 M-0062-0A
 18" x 6"
 18"

 part number is for pair of backets



KF-0170-04 48" L


High-strength, cold-rolled steel shelf support upright allows adjustable shelf bracket to be positioned along 48" length for maximum use of space. A pair of uprights holds a maximum of three shelf brackets.



Adjustable Shelf Brackets		
•	bracket size	
M-0031-AD	11" x 8"	12"
M-0030-AD	15" x 10"	16"
M-0032-AD	17" x 12"	18"
part number is	for single b	oacket



Adjustable Shelf Brackets			
•	bracket		
	size	size	
M-0038-AD		12"	
M-0044-AD		16"	
M-0039-AD	18 ¹ / ₈ " x 6"	18"	
part number is	for single I	oacket	




Adjustable Shelf Brackets				
•	bracket size	shelf size		
M-0060-AD	12" x 6"	12"		
M-0061-AD	16" x 6"	16"		
M-0062-AD	18" x 6"	18"		
part number is fo	r pair of b	ackets		



Laboratory Pegboards

Kemresin Laboratory Pegboards




Pegboard without Drip Trough

29" High x 1" Thick x 20" Long

•		_
	Pegboard Color	Peg Color
X-020014-BK	Black	Black
X-020014-GR	Grey	White
X-020014-PT	Putty	White
X-020014-SL	Slate	White
4.4	614	

Peg Size: ½" dia. x 61/8" long Peg Material: Polypropylene

Number of Pegs: 32



Pegboard without Drip Trough

 29" High x 1" Thick x 32" Long

 Pegboard Color
 Peg Color

 X-020015-BK
 Black
 Black


 X-020015-GR
 Grey
 White

 X-020015-PT
 Putty
 White

 X-020015-SL
 Slate
 White

Peg Size: ½" dia. x 61/8" long Peg Material: Polypropylene

Number of Pegs: 53



Pegboard with Drip Trough


30" High x 1" Thick x 20" Long

•		_
	Pegboard Color	Peg Color
X-020016-BK	Black	Black
X-020016-GR	Grey	White
X-020016-PT	Putty	White
X-020016-SL	Slate	White
4.4	II I' 01/	11.1

Peg Size: ½" dia. x 61/8" long Peg Material: Polypropylene

Number of Pegs: 32

Trough Size: 1" high x 21/4" deep



Pegboard with Drip Trough

30" High x 1" Thick x 32" Long

oo riigir x r riilok x oz zorig							
	Pegboard Color	Peg Color					
X-020017-BK	Black	Black					
X-020017-GR	Grey	White					
X-020017-PT	Putty	White					
X-020017-SL	Slate	White					

Peg Size: ½" dia. x 61/8" long Peg Material: Polypropylene

Number of Pegs: 53

F-2940-00

Trough Size: 1" high x $2\frac{1}{4}$ " deep



Stainless Steel Drip Trough

F-2939-00 X-020014 Pegboard

Size: 1½"H x 2½"W x 20"L

Material: Type 304 Stainless Steel

Includes: ½" OD drain outlet




Stainless Steel Drip Trough

for use on X-020015 Pegboard

Size: 1½"H x 2½"W x 32"L Material: Type 304 Stainless Steel Includes: ½" OD drain outlet

Type 304 Stainless Steel Laboratory Pegboards




Pegboard without Drip Trough

30" High x 3/4" Thick x 24" Long

X-020008 with Open Back X-020010 with Finished Back

Peg Size: ½" diameter x 6" long Peg Material: Polypropylene

Peg Color: White Number of Pegs: 32




Pegboard without Drip Trough

30" High x $\frac{3}{4}$ " Thick x 30" Long

X-020009 with Open Back
X-020011 with Finished Back

Peg Size: ½" diameter x 6" long Peg Material: Polypropylene Peg Color: White

Number of Pegs: 50




Pegboard with Drip Trough

30" High x 3/4" Thick x 24" Long

X-020004 with Open Back
X-020006 with Finished Back

Peg Size: ½" diameter x 6" long Peg Material: Polypropylene

Peg Color: White Number of Pegs: 32



Pegboard with Drip Trough

30" High x 3/4" Thick x 30" Long


X-020005 with Open Back X-020007 with Finished Back

> Peg Size: ½" diameter x 6" long Peg Material: Polypropylene

Peg Color: White Number of Pegs: 50

Miscellaneous Accessories


Writing Boards & Tackboards



Black Chalkboard

F-5336-00 18¹/₄"H × 24¹/₄"L **F-5337-00** 24¹/₄"H × 36¹/₄"L **F-5338-00** 36¹/₄"H × 48¹/₄"L


Surface: washable Duracite Frame Material: anodized aluminum Trough Material: anodized aluminum



White Magnetic Markerboard

F-5342-00 18¹/₄"H × 24¹/₄"L **F-5343-00** 24¹/₄"H × 36¹/₄"L **F-5344-00** 36¹/₄"H × 48¹/₄"L


Surface: 1/4" white magnetic Frame Material: grey mylar vinyl



Green Chalkboard

F-5333-00 181/4"H x 241/4"L **F-5340-00** 241/4"H x 361/4"L **F-5341-00** 361/4"H x 481/4"L

Surface: washable Duracite Frame Material: anodized aluminum Trough Material: anodized aluminum




Cork Tackboard

F-5418-24 18"H × 24"L **F-5424-36** 24"H × 36"L **F-5436-48** 36"H × 48"L

Surface: Resilient natural cork Frame Material: anodized aluminum

Kemresin Balance Table



30" High x 24" Deep x 35" Long

K7-302435-BK Black
K7-302435-GR Grey
K7-302435-PT Putty
K7-302435-SL Slate

Work Top: 2" thick Kemresin
Work Top Size: 2"T x 24"D x 35"L
Pedestal Legs: 2" thick Kemresin
Total Weight: 300 pounds

L-110



Marvel Refrigerators & Freezers

L-110 Refrigerator/Freezer

Size: 34¹/₂" H x 24" D x 23⁷/₈" L 6.1 cubic feet

Door Swing: Right Hand

(User changeable with swing reversing kit)

Shelves: (3) Adjustable
Service: Single-phase,
60 hertz, 120 volt AC

L-105 Refrigerator

Same as above except without freezer compartment.

without freezer compartment.

Same as above except entire unit is a freezer.

Gillis Stainless Steel Cart




F-5349-00

Overall Size: 333/8"H x 19"D x 31"L Shelf Size: (3) 18"D x 27"L

(1" high lip on 3 sides of top & middle shelf)

Capacity: 400 pounds
Casters: (4) Black rubber

4" Diameter - Swivel



Lund Key Cases

F-5332-00

Height: 11"
Depth: 21/2"
Length: 175/16"
Key Capacity: 90 Keys

F-5332-01

Height: 243/4"
Depth: 21/2"
Length: 175/16"
Key Capacity: 200 Keys

Freezer




Miscellaneous Accessories


Pipe Space Kemstruts



Floor mounted support to provide stability to rear overhang of fume hood tops or worksurfaces, and for mounting piping and drainlines. Adjustable to match variances in height.




Corner Kemstruts



Floor and wall mounted support system for adding stability to corner tops in L-shaped furniture assemblies. Adjustable to match variances in height. Includes two support angles to attach to adjoining base cabinets.


NOTE: Kemstrut Corner Top Support Assemblies are not required when a corner base cabinet is used.

Kemstrut Pipe Hangers



K-0005-0A Drainline Support **M-0051-0A** Pipe Clamp
for ³/₈"-1/₂" IPS **M-0052-0A** Pipe Clamp
for ³/₄"-1" IPS


Support Angles



L-2034-06 6" L L-2034-12 12" L L-2034-18 18" L L-2034-24 24" L L-2034-48 48" L L-2034-96 96" L

Formed from 16 gauge cold-rolled steel, with $1^{1}/2^{"}$ sides.

LED Task Lighting



Part Number Length

X-030001-00 300mm - 11½"

X-030002-00 500mm - 19½"

X-030003-00 1000mm - 39¾"

LED Task Lighting consists of a thin-profile, white, LED strip, 9.5mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving and integral tabs for mounting to wood and plastic laminate.

Miscellaneous Accessories

Rod Assembly Support Rods Cross Rods Rod Clamp



Consists of two 0553-00 upright rods, one 0557-00 cross rod and two 0550-00 rod clamps.

0549-18 1/2" Dia. x 18" H 0549-24 1/2" Dia. x 24" H 0549-36 1/2" Dia. x 36" H 0553-00 3/4" Dia. x 36" H

Duraluminum support rods are available in ½" and ¾" diameter and have a tapered end to fit 0548-00 rod sockets.

0554-00 1/2" Dia. x 39" L **0556-00** 1/2" Dia. x 18" L **0557-00** 3/4" Dia. x 39" L

Duraluminum cross rods are available in 1/2" and 3/4" diameter and have rounded ends

0550-00 Rod Clamp


Made of cast aluminum.
Designed to clamp horizontal cross rods on vertical upright rods. Can be used for either 1/2" or 3/4" diameter rods

Rod Sockets

0548-00 Rod Socket

Made of chrome-plated brass. Use with ½" and ¾" diameter upright rods. 1" diameter x 17/8" long threaded shank. Complete with locknut and washer.

Label Holder



F-0260-00 Label Holder

Aluminum with adhesive back for easy installation. Holds 11/8 H x 21/2" L identification card.

Base Molding



F-0862-04 48" L

4" high black base molding. Made of pliable material to match variances in floor surface. Comes with high-strength adhesive applied to back.

Note: Specify number of 4' sections required.

Corner Guard



F-1368-00 Corner Guard

4" high black corner guard to finish off corner joints of base molding.

Wall Cabinet Hanger Rail



L01-020109-00 for 15" wall cabinet for 21" wall cabinet for 24" wall cabinet for 30" wall cabinet for 36" wall cabinet for 36" wall cabinet for 42" wall cabinet for 42" wall cabinet for 42" wall cabinet for 48" wall cabinet

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.



General Plumbing Information

Recommended Supply Line Pressures

Service & Rough-in Abbreviations A – Air

Service Fitting Index Color CodingWhen ordering service fittings, specify

Gas	0.216 psi
Air (Ordinary)	15 psi
Air (Air Powered Stirrer Op	peration)40 psi
Vacuum (Ordinary)	18"-22" of Mercury
Vacuum (High)	28.5" of Mercury
Hot and Cold Water	Local Pressure
Steam (Low Pressure)	0-15 psi as required
Steam (High Pressure)	15-30 psi as required

W – Waste
CW – Cold Water
DW – Distilled Water
HV – Hood Vent
HW – Hot Water

WV - Waste Vent

service required, so that the proper index button(s) will be supplied.

Air – Orange

Air - Orange
Cold Water - Dark Green
Distilled Water - White
Gas - Dark Blue
Hot Water - Red
Hydrogen - Pink
Nitrogen - Brown
Oxygen - Light Green
Special Gases - Light Blue
Steam - Black
Vacuum - Yellow

Recommended Pipe & Drain Line Sizes


Service	Number of Outlets	Pipe Size	Pipe Material
	1-8	3/8"	
Charm			
Steam	9-16	1/2"	In accordance with governing codes
	17-24	3/4"	
	1-10	1/2"	
*Water	11-18	3/4"	In accordance with governing codes
	19-24	1"	
	1-12	3/8"	
Gas	13-18	1/2"	
	19-24	3/4"	In accordance with governing codes
	25-30	1"	
	1-8	3/8"	
Air or Vacuum	9-16	1/2"	In accordance with governing codes
	17-24	3/4"	
	Fittings or	Pipe	Pipe Material
	Filter Pumps	Size	Pipe Material
	1-6 Fittings	11/2"	
Drain Lines	or 4 Pumps	1 72	
	7-10 Fittings		In accordance with governing codes
	or 8 Pumps	2"	
* If table is equipped v	vith a series of filte	r pumps (aspirat	tion), next larger pipe size should be used.

NOTE: For information not supplied, or if additional details are required, contact your local Kewaunee Representative.

Open Reagent Rack Horizontal Pipe Spacing


Water pipe lines should always be the lowest and a minimum of 9" above the working surface. Electric conduit should always be the uppermost service line. (Insulated pipe will require a greater spacing than indicated below.)

Spacing Pipe Center	m Vertical g Between erlines of 2 Pipe Lines
3/8" and 3/8"	2"
3/8" and 1/2"	2"
1/2" and 1/2"	2"
1/2" and 3/4"	$\dots 2^1\!/\!2^{\!"}$
3/4" and 3/4"	$\dots 2^1\!/\!2^{\!"}$
3/4" and 1"	3"
1" and 1"	3"



Minimum Spacing Between Fittings

These drawings show the allowable minimum spacing between service and electrical fittings on laboratory bench tops and ledges.





Laboratory Service Fittings

Laboratory Ball Valves

(for use with air, vacuum or gases – not to be used for liquids, steam or oxygen)








W-0242-00 Gas, Air or Vacuum Valve is chrome plated brass. Valve provides on/ off control of gas, air or vacuum at pressures up to 75 PSI. 3/8" I.P.S. male inlet.

Specify service desired.


W-0286-02 Double Gas, Air or Vacuum Fitting (curb or panel mounted) is chrome plated brass. Consists of two-way 70° outlet panel flange and two laboratory ball valves. Complete with 3/8" I.P.S. tank nipple, locknut and washer. Specify service desired.

W-0287-00 Gas, Air or Vacuum Fitting (curb or panel mounted) is chrome plated brass. Consists of a single outlet panel flange and a laboratory ball valve. Complete with a ³/₈" I.P.S. tank nipple, locknut and washer.

Specify service desired.




W-0260-00 Gas, Air or Vacuum Deck Mounted Fitting is chrome plated brass. Consists of a single outlet turret with a laboratory ball valve. Complete with a ³/₆" I.P.S. tank nipple, locknut and washer. Specify service desired.



W-0263-00 Double Gas, Air or Vacuum Deck Mounted Fitting Same as No. W-0260-00 except consists of a 2-way 90° outlet turret and two laboratory ball valves. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.

Specify service desired.



W-0266-00 Double Gas, Air or Vacuum Deck Mounted Fitting Same as No. W-0260-00, except consists of a 2-way 180° outlet turret and two laboratory ball valves.Complete with a 3/8" I.P.S. tank nipple, locknut and washer.


Specify service desired.





(for use with air, vacuum or gases – not to be used for liquids or steam)







W-0243-00 Gas, Air or Vacuum Valve is chrome plated brass. Single needle point valve with integral straight serrated hose ends. Brass floating cone and integral brass seat. Valve provides metering control of laboratory gases at pressures up to 60 PSI. 3/s" I.P.S. male inlet.

Specify service desired.

NOTE: If valve is to be used for oxygen service, order W-0243-01.

W-0400-00 Double Gas, Air or Vacuum Fitting (curb or panel mounted) with two needle point valves with integral serrated hose ends. Brass floating cone and integral brass seat. Complete with a ³/₈" I.P.S. tank nipple, locknut and washer. Specify service desired.

W-0401-00 Gas, Air or Vacuum Fitting (curb or panel mounted) Same as No. W-0400-00 except single needle point valve. Complete with a ³/8" I.P.S. tank nipple, locknut and washer.

Specify service desired.





Laboratory Service Fittings

Needle Point Cock Valves

(for use with air, vacuum or gases – not to be used for liquids or steam)







W-0259-00 Single Needle

Cock Deck Mounted Fitting Turret Base is chrome plated brass. Single needle point valve with integral straight serrated hose end. Complete with a ³/s" I.P.S. tank nipple, locknut and washer.

Specify service desired.

W-0402-00 Double Needle

Cock Deck Mounted Fitting Turret Base is chrome plated brass. Two needle point valves with integral straight serrated hose ends at 90°. Complete with a ³/_e" I.P.S. tank nipple, locknut and washer.

Specify service desired.

W-0403-00 Double Needle

Cock Deck Mounted Fitting Turret Base is chrome plated brass. Same as No. W-0402-00 except integral straight serrated hose ends are at 180°.Complete with a ³/6" I.P.S. tank nipple, locknut and washer.


Specify service desired.

W-0404-00 Quadruple

Needle Cock Deck Mounted Fitting Turret Base is chrome plated brass. Four needle point valves with integral straight serrated hose ends. Complete with a ³/₈" I.P.S. tank nipple, locknut and washer.

Specify service desired.


Water Valve



W-0306-00 Water Fitting (deck or panel mounted)


Made of chrome plated brass. Consists of a single outlet turret and a straight water valve with removable straight serrated hose connector. Complete with a 3/s" I.P.S. tank nipple, locknut and washer.

Steam Valves



W-0300-00 Steam Fitting (deck or panel mounted)


Made of chrome plated brass. Consists of a single outlet turret and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a 3/s" I.P.S. tank nipple, locknut and washer.



W-0301-00 Steam Fitting (curb or panel mounted)

Made of chrome plated brass. Consists of a single outlet panel flange and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a %e" I.P.S. tank nipple, locknut and washer.

Multiple Service Fitting




W-0273-02 Deck Mounted Combination Water/Gas

Fitting Made of chrome plated brass. Water faucet has a 6" spread gooseneck with removable aerator outlet and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two ³/₈" I.P.S. tank nipples, locknuts and washers.

(Specify service desired for ball valves.)

W-0273-VB Deck Mounted Combination Water/Gas Fitting

Same as W-0273-02 except furnished with vacuum breaker type gooseneck and removable straight serrated hose connector. (Specify service desired for ball valves.)



W-0274-VB Deck Mounted Combination Water/Gas

Fitting Made of chrome plated brass. Water faucet has a 8" spread vacuum breaker type gooseneck with removable straight serrated hose connector and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two 3/8" I.P.S. tank nipples, locknuts and washers.

(Specify service desired for valves.)

W-0274-02 Deck Mounted Combination Water/Gas Fitting

Same as W-0274-VB except furnished with standard gooseneck and removable aerator outlet.

(Specify service desired for valves.)


Laboratory Service Fittings

Single Water Faucet




W-0337-00 Deck Mounted **Gooseneck Water Faucet**

Chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valve is disc type with a renewable internal assembly. Faucet body is furnished with a 3/8" I.P.S. threaded mounting shank with locknut and washer.



W-0337-0V Deck Mounted **Vacuum Breaker Gooseneck Water Faucet**

Same as W-0337-00 except furnished with 8" spread rigid vacuum breaker type gooseneck.




W-0337-06 Deck Mounted **Gooseneck Water Faucet**

Same as W-0337-00 except furnished with 6" spread rigid gooseneck




W-0337-6V Deck Mounted Vacuum Breaker **Gooseneck Water Faucet**

Same as W-0337-0V except furnished with 6" spread rigid vacuum breaker type gooseneck.




Mounted Gooseneck Water

Faucet Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valve is disc type with renewable internal assembly. Flanged inlet is complete with 1/2" I.P.S. tank nipple, locknut and washer.




W-0350-00 Right Hand Curb W-0350-0V Right Hand Curb **Mounted Vacuum Breaker Gooseneck Water Faucet**

Same as W-0350-00 except supplied with vacuum breaker type gooseneck.




W-0351-00 Left Hand Curb **Mounted Gooseneck Water**

Faucet Same as W-0350-00 except supplied with valve handle on left.



W-0351-0V Left Hand Curb **Mounted Vacuum Breaker Gooseneck Water Faucet**

Same as W-0350-0V except supplied with valve handle on




W-0405-00 Portable Pump Gooseneck Faucet for labo-

ratory demonstration table. Made of chrome plated brass with 75/8" spread, rigid, 1/4" I.P.S. gooseneck. Faucet/pump body is complete with an integral threaded mounting shank and flanged locknut. 5/8" O.D. hose connector inlet.




Laboratory Service Fittings

Mixing Faucet




W-0340-00 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet Made

of chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.




W-0340-0V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet Same as W-0340-00 except supplied with


vacuum breaker type gooseneck.



W-0340-06 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet Same as W-0340-00 except supplied with 6" spread rigid gooseneck.




W-0340-6V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet Same as W-0340-06 except supplied with vacuum breaker type gooseneck.




W-0333-00 Deck Mounted Hot and Cold Water Mixing Swing Spout Faucet Made

of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.




W-0354-00 Curb Mounted Hot and Cold Water Mixing Gooseneck Faucet Made

of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with ³/₈" I.P.S. tank nipples, locknuts and washers on 8" centers.



W-0354-0V Curb Mounted Hot and Cold Water Mixing Vacuum Breaker

Gooseneck Faucet Same as No. W-0354-00 except supplied with vacuum breaker type gooseneck.




W-0355-00 Curb Mounted Hot and Cold Water Mixing Swing Spout Faucet


Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with ³/₈" I.P.S. tank nipples, locknuts and washers on 8" centers.

Laboratory Service Fittings


Mixing Faucet



W-0361-00 Panel Mounted Hot and Cold Water Mixing Gooseneck Faucet Chrome plated brass. Designed for mounting on facia type reagent rack plumbing box uprights. It has an 8" spread rigid gooseneck equipped with a removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets on 3" centers are complete with 3/s" I.P.S. tank nipples, locknuts and washers.



W-0361-0V Panel Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet Same as W-0361-00 except supplied with vacuum breaker type gooseneck.



W-0825-00 Floor Mounted Foot Control Double Pedal Mixing Valve is chrome plated brass. Has twin pedals with self-closing valves for mixing hot and cold water. 1/2" I.P.S. female thread inlets and outlet at the rear. Inlets are on 2" centers, 7/8" above floor. Outlet is on center 11/2" above floor.




W-0833-00 Deck Mounted Water Gooseneck Fitting is chrome plated brass. Consisting of a 8" gooseneck with mounting flange and W-0532-00 aerator. ³/₈" I.P.S. female thread inlet complete with tank nipple, locknut and washer. Intended for use as an outlet for W-0825-00.

Faucet Accessories



W-0532-00 Aerator Antisplash aerator can be used in lieu of the removable straight serrated hose connectors supplied with the gooseneck water faucet. 3/8" I.P.S. male thread inlet.



W-0510-00 Filter Pump (Aspirator) 51/2" high overall. Has serrated hose connector side and bottom outlets. 3/8" I.P.S. male thread inlet.




W-0535-00 In-Line Vacuum Breaker Chrome plated. Mounts between the outlet end of a fitting (gooseneck, flange, etc.) and a removable serrated hose connector. ³/s" I.P.S. male thread inlet and female thread outlet. (Check plumbing codes for required location of vacuum breakers before ordering.)



Laboratory Service Fittings

Distilled Water Faucet





W-0367-00 Right Hand Distilled Water (Curb or Panel Mounted) Faucet Chrome

plated with a hot tin lined interior. It has a self-closing handle and a disc type valve with a renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. 1/2" I.P.S. flanged female thread inlet complete with tank nipple, locknut and washer.


W-0368-00 Deck Mounted Distilled Water Gooseneck

Faucet Made of chrome plated brass. Has a hot tin lined interior. It has a self-closing handle and a disc type valve with renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlet with coupling nut and 1/4" I.P.S. tailpiece.

Ultra-Pure Deionized Water Faucet







W-0407-00 Distilled, Deionized and Ultra-Pure Water Faucet (deck mounted)

made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

W-0408-00 Distilled, Deionized and Ultra-Pure Water Faucet (panel mounted)


made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

W-0409-00 Distilled Deionized and Ultra-Pure Water Faucet (panel mounted)


made of chrome plated brass. All polypropylene internal construction with 8" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

Laboratory Service Fittings

Vandal-Resistant Fittings




W-0410-VR Vandal Resistant W-0411-VR Vandal Resistant W-0412-VR Vandal Resistant **Deck Mounted Gas, Air or** Vacuum Fitting with single vandal-resistant laboratory ball valve. Complete with a 3/8" I.P.S. tank nipple, locknut and washer. Specify service desired.




Deck Mounted Double Gas, Air or Vacuum Fitting with two vandal resistant laboratory ball valves with integral straight serrated hose ends at 90°

Complete with a 3/8" I.P.S. tank nipple, locknut and washer.


Specify service desired.




Deck Mounted Double Gas, Air or Vacuum Fitting Same as W-0411-VR except at 180° Complete with a 3/8" I.P.S. tank nipple, locknut and washer. Specify service desired.




W-0413-VR Vandal Resistant W-0414-VR Vandal Resistant W-0415-VR Vandal **Deck Mounted Hot and** Cold Mixing Faucet with 8" rigid vacuum breaker gooseneck. Valves are disc type with renewable and mounting shanks. internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlets with coupling nuts for slip joint connections.



Panel Mounted Hot and Cold Mixing Faucet with 8" rigid vacuum breaker gooseneck




Resistant Deck Mounted Water Faucet with 8" rigid vacuum breaker gooseneck. Faucet body is furnished with a %" I.P.S. threaded mounting shank with locknut and washer.



Hand-Held Drench Hose with single spray type outlet heads, squeeze-lever operated valve and flexible stainless steel hose. 3/8" I.P.S male inlet.

Resistant Deck Mounted


W-0417-VR Vandal



W-0416-VB Vandal Resistant **Deck Mounted Combination** Vacuum Breaker Water/Gas

Fitting Made of chrome plated base. Rigid gooseneck with disc type valve and in-line vacuum breaker. Includes two-way turret with two laboratory ball valves. Complete with two 3/8" I.P.S. tank nipples, locknuts and washers.

Specify service desired.




W-0416-VR Vandal Resistant **Deck Mounted Combination** Water/Gas Fitting Same as W-416-VB except furnished without in-line vacuum breakers. Specify service desired.




Laboratory Service Fittings


Emergency Eye Wash/Shower Equipment




W-0927-00 Panel Mounted Eye Wash Fixture Consists of stainless steel bowl, two Gentle Spray outlet heads, a "PUSH TO OPERATE" ball valve and mounting bracket. 1/2" I.P.S. female thread inlet and 11/4" I.P.S. female thread outlet.




W-0933-00 Eye/Face Wash (wall mounted) the same as No. W-0927-00, except with four Gentle Spray outlet heads for complete coverage of eyes and face.




W-0930-00 BARRIER-FREE (ADA) **Eye/Face Wash (wall mounted)** with four Gentle Spray outlet heads, stainless steel skirt and bowl and stayopen ball valve. 1/2" I.P.S. female thread inlet and 11/4" I.P.S. female thread outlet.




W-0931-00 Eye Wash (deck mounted) with two Gentle Spray outlet heads, stay-open ball valve and 90° swivel feature (right hand mounting).




W-0932-00 Eye Wash (deck mounted) the same as No. W-0931-00 but (left hand mounting).




W-0923-00 Deck Mounted Hose Spray Fitting Consists of 8 ft. reinforced PVC hose with a squeeze-lever operated valve and spray-type outlet head with a protective rubber collar. Flange, shank and locknut for attaching to table top are included. ³/₈" I.P.S. male thread swivel inlet.



W-0934-00 Deck Mounted Dual Purpose Eye Wash/Drench Hose with two Gentle Spray outlet heads angled at 45° and squeeze valve with locking clip. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




W-0935-00 Wall Mounted Dual Purpose Eye Wash/ Drench Hose the same as No. W-0934-00 except wall mounted.




W-0374-01 Emergency Shower Head Assembly (horizontally mounted).

Includes 10" diameter orange ABS plastic shower head. Operated by a rigid pull-rod and easy to grip handle. Complete with an elbow, connecting pipe and ball valve. 1" I.P.S. female thread inlet.



W-0928-00 Emergency Shower Head Assembly (vertically mounted). The same as No. W-0374-01 except vertically mounted.



W-0926-00 Drench Shower and Eye Wash Combination

A necessity for chemical stock-rooms and other areas handling a quantity of corrosive chemicals. Free-standing 10" dia. drench shower made of orange ABS plastic. Has a 1" ball valve operated by a rigid pull-rod and easy to grip handle. Eye wash fixture with stainless steel bowl, two Gentle Spray heads and a "PUSH TO OPERATE" ball valve. Eye wash and shower operate independently. Overall height is 80". 11/4" I.P.S. female thread inlet and outlet.

W-0936-00 BARRIER-FREE (ADA)
Drench Shower and Eye Wash
Combination Similar to No. W-0926-00
except this unit is barrier-free and eye wash
furnished with four Gentle Spray heads.

Cupsinks & Drain Accessories

Cupsinks



0490-BE Ore Department of the strainer. 6"x 3" inside dimension. 11/2" I.P.S. male straight thread outlet. (Overall height is 71/2")



0491-BP O491-GPO491-GPBlack Poly Grey PolyComplete with removable

strainer. 6"x 3" inside dimension. 11/2" I.P.S. male straight thread outlet. (Overall height is 81/2")



0492-BP 0492-GP

Complete with removable strainer. 9"x 3" inside dimension. 11/2" I.P.S. male straight thread outlet. (Overall height is 7")

Black Poly

Grey Poly



0493-BE 0493-GE 0493-PE 0493-SE

Black Epoxy Grey Epoxy Putty Epoxy Slate Epoxy

Complete with removable type 316 stainless steel wire screen. 13³/s"x 4¹/s" x 5¹/4" I.D. 1¹/2" I.P.S. male straight thread outlet. (Overall height is 8")

Stainless Steel Cupsinks




0950-0F Stainless Steel Made of type 304 stainless steel and has integral cross bars. 5" I.D. 11/2" I.P.S. male straight thread outlet.




0975-0F Stainless Steel Made of type 316 stainless steel and has integral cross bars. 6"x 3" I.D. 11/2" I.P.S. male straight thread outlet.

Drain Accessories



0464-00 Polyethylene Trap and Sink Outlet Adjustable riser (tailpiece) can be cut to the desired length. Sink outlet is complete with plug strainer, gasket and locknut. Has a 11/2" I.P.S. male straight thread outlet. Running type trap can be used for either "P" or "S" type drain connections. 11/2" I.P.S. male straight thread outlet.



0469-00 Polyethylene Trap and Sink Outlet Adjustable riser (tailpiece) with 41/2" adjustment that does not require cutting. The sink outlet is complete with plug, strainer, gasket and locknut. Has a 11/2" I.P.S. male straight thread outlet. Trap is the bottle "P" type, with removable bottom for cleanout. 11/2" I.P.S. male straight thread outlet.

Sink Outlet

Sink Overflow



0482-BP Black Poly **0482-GP** Grey Poly

Complete with locknut. 1¹/2" I.P.S. male straight thread outlet. For use with epoxy resin tub sinks and troughs.





 0494BP-04
 4" High

 0494BP-06
 6" High

 0494BP-08
 8" High

 0494BP-10
 10" High

 Sink Overflows are designed

Sink Overflows are designed to be used with 0482-BP sink outlets. 11/2" diameter.

St. Steel Sink Drain



0995-00 Stainless Steel Crumb Cup Strainer Complete with gaskets, locknut and basket strainer. 11/2" male straight thread outlet.



Kemresin Tub Sinks & Drain Troughs

Kemresin® Tub Sinks and Drain Troughs are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting material. They are formed in metal molds to provide one-piece units with coved corners and bottoms pitched to the drain outlet. After being removed

from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1¹/₂" sink outlets unless otherwise specified.

Kemresin Drain Troughs are molded with no longitudinal joints. A slope of 1/8" per foot assures complete drainage, and a 7/8" cove radius

in the bottom corners is provided for ease of cleaning. Five trough sections, or any part of them, are available to be cemented together to provide Open End Troughs up to 21' long, Closed End Troughs with an end drain outlet up to 21'-8" long and Closed End Troughs with a center drain outlet up to 42'-8" long. Section TR-4 is provided with a center drain recess for a standard 11/2" sink outlet.

Epoxy Resin Drain Troughs Epoxy Resin Tub Sinks Closed End Trough Double Tub Sink Open End Trough Section TR-01-00-Single Tub Sink 53/8" Section TR-02-00-**Center Bench End Sink** with Trough Cutout Section TR-03-00-NOT LESS THAN 11/8 31/81 Section TR-05-00-BOTTOM PITCHED TO DRAIN OUTLET \neg 23/4 31/81 **OVÉRALL HEIGHT** THAN 11/8 Section TR-04-00-NOT MORE THAN 5°

TR-02-00-

TR-01-00-

TR-04-00-

TR-01-00

TR-05-00-

TR-03-00-

Kemresin Tub Sinks & Drain Troughs

Kemresin Colors Color Codes

Kemresin® Tub Sinks and Drain Troughs are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number. **Example: 1005-00-BK**

BK = Black **PT** = Putty **GR** = Grey **SL** = Slate

Kemresin Tub Sinks

		Inside Dimensions (at Top of Sink)			Drain Location		
Catalog Number	Length	Width	Height	Α	В	Height	
1000-00	14"	10"	6"	4"	5"	67/8"	
1001-00	16"	8"	7"	4"	4"	8"	
1003-00	18"	15"	11 ¹³ /16"	4"	4"	12 ⁷ /8"	
(1) 1003-DT-	37"	15"	11"	41/2"	41/2"	12"	
1005-00	25"	15"	10"	4"	4"	11½"	
1006-00	18"	15"	8"	4"	4"	9"	
(2) 1006-DT-	36¾"	15"	81/4"	41/2"	41/2"	9"	
1007-00	24"	16"	8"	4"	8"	91/8"	
1008-00	16"	12"	8"	4"	6"	91/8"	
1011-00	18"	61/2"	6"	9"	31/4"	7"	
(3) 1014-00	25"	15"	10"	4"	4"	111/4"	
(4) 1015-00	233/8"	133/4"	131/4"	11 ¹¹ /16"	91/2"	173/4"	
(<u>5)</u> 1016-00	23%"	13¾"	131/4"	11 ¹ / ₁ 6"	91/2"	17¾"	
1021-00	16"	16"	71/2"	4"	4"	8¾"	
1022-00	21"	18"	10"	41/2"	9"	11"	
1023-00	21½"	15½"	7"	4"	4"	8%"	
1024-00	21½"	151/2"	11"	4"	4"	121/8"	
(6) 1024-DT-	44"	15½"	11"	4"	4"	12 ¹ /2"	
1025-00	28"	20"	91/8"	14"	6"	101/8"	
1026-00	351/2"	19½"	10"	173/4"	93/4"	11½"	
1027-00	42"	20"	10"	21"	10"	11 ¹ / ₄ "	
1030-00	16"	12"	8"	8"	6"	9"	
1031-00	24"	16"	8"	12"	8"	9"	
1032-00	293/4"	143/4"	6"	14 ⁷ /8"	73/8"	7"	
1034-00	21½"	191/4"	8"	103/4"	95/8"	91/8"	
1035-00	45¾"	173/8"	121/2"	22 ⁷ /8"	83/4"	13 ³ / ₈ "	
1036-00	485/8"	241/2"	5 ³ / ₈ "	33/4"	121/4"	6"	
(7) 1037-00-	231/4"	141/4"	123/4"	115/8"	91/4"	223/8"	
(8) 1038-00-		ameter	75/8"	cer	8%"		
(9) 1039-00		ameter	81/8"		nter	8 ⁷ /8"	
SM-01-00	. 12"	6"	6"	6"	3"	7"	
SM-02-00	19½"	161/2"	8"	93/4"	81/4"	9"	
SM-03-00	. 24"	21"	8"	6"	41/2"	91/8"	
SM-04-00	251/2"	161/2"	19"	41/2"	41/2"	201/8"	
SM-05-00	_ 28"	18"	10"	4"	4"	111/8"	
SM-06-00-	28"	18"	10"	24"	4"	11½"	
SM-07-00	30"	16"	10"	4"	4"	111/8"	
SM-08-00	30"	16"	10"	26"	4"	111/8"	
SM-09-00	30"	18"	18"	41/2"	43/4"	18 ³ / ₄ "	
SM-10-00	36"	21"	43/4"	4"	4"	53/4"	
SM-11-00	48"	13"	121/4"	24"	61/2"	131/4"	
SM-12-00-	8"	8"	8"	4"	4"	8 ⁷ ⁄8"	

(1) Double tub sink. Distance between centerlines of drain outlets is 193	/8".
---	------

- (2) Double tub sink. Distance between centerlines of drain outlets is 19".
- (3) Sink has 71/8" long by 55/8" high trough cutout.
- (4) Center table end sink has 7¹/s" long by 5⁵/s" high trough cutout and 24¹/₂" long by 3³/₄" high back splash. Table side of sink is not tapered.
- (5) Center table end sink has 241/2" long by 33/4" high back splash without trough cutout. Table side of sink is not tapered.
- (6) Double tub sink. Distance between centerlines of drain outlets is 221/2".
- (7) Assembly of tub sink, back splash, ledge and end fillers. Used for T69W4525-00 Scrub Sink Assembly.
- (8) Cylindrical sink -12" inside diameter at top of sink with center drain outlet.
- (9) Cylindrical sink -18" inside diameter at top of sink with center drain outlet.
- (10) Sides have 11/2" taper per foot in height.
- (11) Exterior of one 21" long side is not tapered.

		de Dimens Top of Si	Dra Loca	Overall Sink		
Catalog Number	Length	Width	Height	Α	В	Height
SM-13-00	14"	10"	8"	4"	5"	9"
SM-14-00	231/4"	14 ¹ / ₄ "	123/4"	115/8"	91/4"	135%"
SM-15-00	19 ³ / ₄ "	18"	24"	97/8"	9"	25"
SM-16-00	40"	18"	181/4"	20"	9"	191/4"
SM-17-00	29¾"	231/2"	6"	14 ⁷ /8"	113/4"	7"
SM-18-00	18"	15"	15 ¹ / ₄ "	43/4"	43/4"	16½"
SM-19-00	351/2"	19½"	15 ³ / ₈ "	31/2"	41/2"	161/2"
SM-20-00	21"	18"	151/4"	41/2"	41/2"	16½"
SM-21-00	601/2"	22"	71/4"	6"	6"	81/8"
SM-22-00	24"	16"	15 ³ / ₈ "	12"	8"	161/%"
SM-23-00	28"	18"	151/2"	51/2"	4"	16½"
SM-24-00	48"	8"	8"	24"	4"	9"
SM-25-00	21"	18"	151/8"	101/2"	9"	16½"
SM-26-00	21"	11%"	151/8"	101/2"	5 ¹ 1/16"	161/8"
SM-27-00	36"	20"	16 ¹ / ₄ "	18"	10"	17 ¹ / ₄ "
SM-28-00	30"	18"	41/8"	15"	9"	51/4"
SM-29-00	36"	24"	121/4"	18"	12"	131/4"
SM-30-00	40"	30"	81/8"	20"	15"	9"
SM-31-00	24"	15¾""	121/8"	43/4"	8"	13"
SM-32-00	59"	21"	125/8"	291/2"	101/2"	133/8"
SM-33-00	24"	18"	11"	4"	9"	11 ⁷ /8"
(10) SM-34-00-	15"	14"	101/2"	71/2"	7"	1111/4"
SM-35-00	51"	17"	121/2"	251/2"	81/2"	131/4"
SM-36-00	14"	10"	6"	7"	5"	6 ⁷ /8"
(11) SM-37-00-	24"	20"	105/8"	63/4"	10"	111/4"
SM-39-00	35	65/8"	43/4"	17½"	35/16"	51/2"
SM-40-00	333/8"	291/8"	153/4"	6"	14%16"	16 ¹ /2"
SM-41-00	18"	15"	81/4"	9"	71/2"	9"
SM-42-00	25"	15"	101/2"	121/2"	71/2"	111/4"
SM-43-00	21"	18"	103/8"	101/2"	9"	111/8"
SM-44-00	36"	16"	11"	12"	4"	12"
SM-45-00-	273/4"	143/4"	81/2"	101/2"	71/2"	93/8"
SM-47-00	43"	15½"	111/8"	21½"	73/4"	12"
SM-48-00	56" 36"	24" 16"	8 ¹ / ₄ " 8 ³ / ₈ "	28" 6"	12" 8"	9" 9"
SM-49-00	24"	24"	12 ¹ /2"	12"	12"	
SM-50-00- SM-51-00-	30"	15 ⁷ /8"	12½" 16 ³ /8"	6"	7 ¹⁵ / ₁₆ "	13 ¹ / ₄ " 17 ¹ / ₈ "
SM-51-00	45½"	111/2"	8 ¹ / ₈ "	229/16"	53/4"	9"
SM-53-00-	31½"	16"	16 ¹ /2"	71/8"	8"	9 17½"
SM-53-00	16"	16"	10 /2	4 ¹ / ₄ "	4 ¹ / ₄ "	17 74
SM-55-00-	421/4"	18 ¹ / ₈ "	12 ⁵ / ₈ "	211/8"	91/16"	13 ³ /8"
3141-30-00	42 74	10 78	1278	2178	3 /16	1378

General Notes:

- Sinks with center drain outlets, except 1011-00, can be furnished with a drain recess sized for a crumb cup strainer or garbage disposal when specified.
- Any sink may be furnished with trough cutout. Specify location and size. (Trough location must be at least 11/2" from inside corner and 2" from bottom.)
- Most sinks can be joined together to form multiple tub sinks provided adjoining sides of sinks are the same size. Sinks longer than 30", wider than 24" or deeper than 12" cannot be factory joined, but they can be factory prepared for field assembly. Assemblies of three or more tub sinks should be no longer than 60" overall. For longer assemblies, one factory joint can be made with a second joint prepared for field assembly.



Kemresin Drop-in Sinks


Kemresin Drop-in Sinks are formed in metal molds to provide a lipped one-piece tub with coved corners and bottoms pitched to the drain outlet. They are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting


material. After being removed from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1½" sink outlets unless otherwise specified.

They are designed to be mounted in Kemresin tops in a specially rabbetted cutout so that no support from below is required.


Drop-in Kemresin Sinks









	I.D.	Depth	O.D.	Drain Location
Part Number	A B	С	D E	
1000-DI	13 ⁷ /8" 9 ⁷ /8"	6"		½" 7" 5"
1001-DI	8" 16"	7"	10¼" 18	1/4" 4" 4"
1003-DI	17 ¹⁵ /16"14 ¹⁵ /10			3½" 3½"
1005-DI	25" 15"	10"	27¼" 17	
1006-DI	17 ⁷ /8" 14 ⁷ /8'		201/8" 17	
1008-DI	16" 12"	8"		1/4" 4" 6"
1011-DI	18" 6½"	6"	201/4" 83	
1021-DI	16" 16"	7½"	181/4" 18	
1022-DI	21" 18"	10"	231/4" 20	
(1) 1038-DI-	12" diamete		14½" diame	
(2) 1039-DI-	18" diamete	r 8½"	201/4" diame	
1040-DI	15%" 11%"			.6" 4" 4"
1041-DI	23 ³ / ₈ " 15 ³ / ₈ "		25.6" 17	.6" 4" 4"
1042-DI	28" 15"	11.8"	301/4" 17	
1043-DI	30" 16"	17.8"	321/4" 18	1/4" 5" 8"
1044-DI	12" 8" 21" 17"	6"	14 ¹ / ₄ " 10	1/4" 6" 4"
1045-DI		10"		1/4" 41/2" 81/2"
1046-DI	21½" 15½"		23 ³ / ₄ " 17 26 ¹ / ₄ " 20	³ / ₄ "
1047-DI	24" 18" 17½" 12¾"	11" ' 6½"		5/8" 4 4 8 4 6 5/8" 8 5/4" 6 6 7 1 6 1 6 1
(3) 1048-DI- (4) 1049-DI-	18" 15"	11"		78
1070-DI-	24" 16"	15½"		¹ / ₄ "
	24 10	13/2	20 /4 10	/4 12 0
ADA Compliant				
1000-ADA-DI	13 ⁷ / ₈ " 9 ⁷ / ₈ ""			1/8" 4" 4"
1005-ADA-DI	25" 15"	4.8"	27¼" 17	
1006-ADA-DI	17 ⁷ /8" 14 ⁷ /8'	' 5"	20½" 17	1/8" 31/2" 31/2"
(1) Cylindrical sink, 12" inside o			iangular sink used	in X92 and X94 assemblies







SECTION VIEW 1049-DI-

Kemresin Colors

(2) Cylindrical sink, 18" inside diameter

Kemresin® Drop-in Sinks are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.

(4) Stepped sink, 43/4" deep at front

Example: 1005-DI-BK

Color Codes

BK = Black **PT** = Putty **GR** = Grey **SL** = Slate

Stainless Steel Drop-in Sinks


Seamless die-drawn from nickel bearing type 304, 18-8, 18 gauge stainless steel. Allows easy cleaning and stain resistance. Interior and top surfaces polished to a gleaming finish with a highlighted bowl rim. Corners and bottom turns are radiused 13/4". Sinks are thoroughly sound deadened with environmentally

safe compounds, to ensure proper sound insulation and to reduce condensation.

All sinks are self-rimming with stainless steel mounting channels for easy, top mount, snap-in clip installation. Sinks are compatible with any type worksurface material from 3/4" to 21/2" thick.

Each sink is center punched to receive a standard cup strainer. ADA compliant sinks are punched in the center, 4-1/2" from the rear of the bowl for extra clearance.


Stainless Steel Sinks




	_				
D	Front to	Side	Front	Side	Sink
Part Number	Back	Side	Back	Side	Depth
F-5353-08	15"	15"	12"	12"	71/2"
F-5354-08	17"	19"	14"	16"	7½"
F-5355-08	19"	22"	16"	19"	7 1/2"
F-5356-08	19"	25"	16"	22"	7½"
F-5357-08	19"	31"	16"	28"	71/2"
F-5353-09	15"	15"	12"	12"	9"
F-5354-11	17"	19"	14"	16"	10½"
F-5355-11	19"	22"	16"	19"	10½"
F-5356-11	19"	25"	16"	22"	10½"
F-5357-11	19"	31"	16"	28"	10½"
F-5354-12	17"	19"	14"	16"	12"
F-5355-12	19"	22"	16"	19"	12"
F-5356-12	19"	25"	16"	22"	12"
F-5357-12	19"	31"	16"	28"	12"

ADA Compliant

	0	O.D.		D.	
Part Number	Front to Back	Side to Side	Front to Back	Side to Side	Sink Depth
F-5353-06	15"	15"	12"	12"	51/2"
F-5354-06	17"	19"	14"	16"	5½"
F-5355-06	19"	22"	16"	19"	5½"
F-5356-06	19"	25"	16"	22"	5½"
F-5357-06	19"	31"	16"	28"	5½"




Stainless Steel Sinks with Fixture Ledge




	U	.ט.	1.	I.D.		Г
	Front	Side	Front to	Side to	Sink	Fixture
Part Number	to Back	to Side	Back	Side	Depth	Centers
F-5358-08	18"	15"	12"	12"	71/2"	4"
F-5359-08	20"	19"	14"	16"	71/2"	6"
F-5360-08	22"	22"	16"	19"	71/2"	6"
F-5361-08	22"	25"	16"	22"	8"	8"
F-5362-08	22"	31"	16"	28"	71/2"	8"
F-5358-09	18"	15"	12"	12"	9"	4"
F-5359-11	20"	19"	14"	16"	10½"	6"
F-5360-11	22"	22"	16"	19"	101/2"	6"
F-5361-11	22"	25"	16"	22"	101/2"	8"
F-5362-11	22"	31"	16"	28"	10½"	8"
F-5360-12	22"	22"	16"	19"	12"	6"
F-5361-12	22"	25"	16"	22"	12"	8"
F-5362-12	22"	31"	16"	28"	12"	8"

ADA Compliant with Fixture Ledge

	0	O.D.		I.D.		F
Part Number	Front to Back	Side to Side	Front to Back	Side to Side	Sink Depth	Fixture Centers
F-5358-06	18"	15"	12"	12"	5½"	4"
F-5359-06	20"	19"	14"	16"	51/2"	6"
F-5360-06	22"	22"	16"	19"	5½"	6"
F-5361-06	22"	25"	16"	22"	51/2"	8"
F-5362-06	22"	31"	16"	28"	51/2"	8"







Laboratory Electrical Fixtures

Electrical Receptacles

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. Pedestals are furnished with device, wrinkled black epoxy powder coated pedestal box, stainless steel face plate, tank nipple and locknut. All devices and boxes are UL listed. (Ship loose)









			_			_	
Type 2 and Type 3 pedestal boxes furnished witank nipple, Type 7 and Type 8 pedestal boxes with 1" IPS tank nipple.		Flush mount single gang box	Type 2 – Si	le Pedestal ingle Faced ouble Faced	Flush mount double gang box	Type 7 – Si	de Pedestal ngle Faced ouble Faced
Electrical Device	Color	Type 1	Type 2	Type 3	Type 6	Type 7	Type 8
Specification Grade 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20	Black Ivory White Grey Red	S581-1K S581-1V S581-1W S581-1G S581-1R	S581-2K S581-2V S581-2W S581-2G S581-2R	S581-3K S581-3V S581-3W S581-3G S581-3R	S581-6K S581-6V S581-6W S581-6G S581-6R	S581-7K S581-7V S581-7W S581-7G S581-7R	S581-8K S581-8V S581-8W S581-8G S581-8R
Specification Grade 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20	Black Ivory White Grey Red	S681-1K S681-1V S681-1W S681-1G S681-1R	S681-2K S681-2V S681-2W S681-2G S681-2R	S681-3K S681-3V S681-3W S681-3G S681-3R	S681-6K S681-6V S681-6W S681-6G S681-6R	S681-7K S681-7V S681-7W S681-7G S681-7R	S681-8K S681-8V S681-8W S681-8G S681-8R
Hospital Grade 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20	Black Ivory White Grey Red	H581-1K H581-1V H581-1W H581-1G H581-1R	H581-2K H581-2V H581-2W H581-2G H581-2R	H581-3K H581-3V H581-3W H581-3G H581-3R	H581-6K H581-6V H581-6W H581-6G H581-6R	H581-7K H581-7V H581-7W H581-7G H581-7R	H581-8K H581-8V H581-8W H581-8G H581-8R
Specification Grade 240 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 6-20	Black Ivory White Grey Red	S582-1K S582-1V S582-1W S582-1G S582-1R	\$582-2K \$582-2V \$582-2W \$582-2G \$582-2R	\$582-3K \$582-3V \$582-3W \$582-3G \$582-3R	S582-6K S582-6V S582-6W S582-6G S582-6R	S582-7K S582-7V S582-7W S582-7G S582-7R	S582-8K S582-8V S582-8W S582-8G S582-8R
Ground Fault Protected 120 Volt AC, 20 amp, specification grade, duplex receptacle NEMA 5-20	Black Ivory White Grey Red	S656-1K S656-1V S656-1W S656-1G S656-1R	S656-2K S656-2V S656-2W S656-2G S656-2R	S656-3K S656-3V S656-3W S656-3G S656-3R			
Ground Fault Protected 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20	Black Ivory White Grey Red				S657-6K S657-6V S657-6W S657-6G S657-6R	S657-7K S657-7V S657-7W S657-7G S657-7R	S657-8K S657-8V S657-8W S657-8G S657-8R
Ground Fault Protected 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20	Black Ivory White Grey Red				S658-6K S658-6V S658-6W S658-6G S658-6R	S658-7K S658-7V S658-7W S658-7G S658-7R	S658-8K S658-8V S658-8W S658-8G S658-8R

Laboratory Electrical Fixtures

Electrical Switches and Data Ports

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. All devices and boxes are UL listed. (Ship loose)





			Flush mount single gang bo	υX		n mount gang box	
Electrical Device	Color		Cat-6	Cat-6a	Cat-6	Cat-6a	
Specification Grade Single pole, 120/240 Volt AC, 20 amp toggle switch.	Black Ivory White Grey Red	S655-1K S655-1V S655-1W S655-1G S655-1R					
RJ-45 Data Port Duplex Cat-6 Data Port	Black Ivory White Grey Red		0720-1K-6 0720-1V-6 0720-1W-6 0720-1G-6 0720-1R-6	0720-1K-6a 0720-1V-6a 0720-1W-6a 0720-1G-6a 0720-1R-6a			
RJ-45 Data Port Duplex Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, duplex receptacle. NEMA 5-20	- I write				\$721-6K-6 \$721-6V-6 \$721-6W-6 \$721-6G-6 \$721-6R-6	\$721-6K-6a \$721-6V-6a \$721-6W-6a \$721-6G-6a \$721-6R-6a	
RJ-45 Data Port Duplex Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, GFI duplex receptacle. NEMA 5-20	Black Ivory White Grey Red				\$751-6K-6 \$751-6V-6 \$751-6W-6 \$751-6G-6 \$751-6R-6	\$751-6K-6a \$751-6V-6a \$751-6W-6a \$751-6G-6a \$751-6R-6a	



Laboratory Seating

Specifications:

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- Upholstered seat/backrest 5-legged base with 18" fixed footring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
F-4940-00	163/4"-22"	81/2"
F-4941-00	221/4"-271/2"	81/2"
F-4942-00	261/4"-311/2"	12 ¹ / ₂ "



- Upholstered seat/backrest 5-legged base with 18" fixed foot ring
- Tubular base
- Backrest with lumbar support

	Seat Height	Footring
Model No.	Adjustment	Height
F-4943-00	16 ³ / ₄ "–22"	81/2"
F-4944-00	221/4"-271/2"	8 ¹ /2"
F-4945-00	261/4"-311/2"	12 ¹ / ₂ "

Laboratory Seating

Specifications:

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- 13" round chrome-plated seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
F-4946-00	18 ¹ / ₂ "–22 ³ / ₄ "	8 ¹ /2"
F-4947-00	22 ¹ / ₂ "–27 ³ / ₄ "	12 ¹ /2"



- 13" round wood seat
- Tubular 4-legged base with 16" fixed foot ring
- Black tubular steel base

Model No.	Seat Height
F-5350-00 F-5351-00	18" 24"
F-5352-00	30"




- 15" round upholstered seat
- 5-legged base with 18" fixed foot ring
- Tubular base


Model No.	Seat Height Adjustment	Footring Height
F-4948-00	21 ¹ / ₂ "-24 ³ / ₄ "	8 ¹ /2"
F-4949-00	25 ¹ / ₂ "-20 ³ / ₄ "	12 ¹ /2"



Inset Style Wall Bench End Views




Inset Style Center Bench End Views



Center Bench End/Section View

Center Bench End Sink Cabinet

Inset Style Full Height Storage Cabinet End Views




12³/₄" Deep Full Height Cabinet

16" Deep Full Height Cabinet

22" Deep Full Height Cabinet


Inset Style Center Bench Plan View




Plan View
Center Bench End Sink Cabinet



Full Overlay Style Wall Bench End Views




Full Overlay Style Center Bench End Views



Center Bench End/Section View

Center Bench End Sink Cabinet

Full Overlay Style Full Height Storage Cabinet End Views




12³/₄" Deep Full Height Cabinet

16" Deep Full Height Cabinet

22" Deep Full Height Cabinet

Full Overlay Style Center Bench Plan View



Plan View Center Bench End Sink Cabinet



Lock Schedule Sample and Examples

Lock Schedule Examples

Keyed Different & Masterkeved

An example of a simple lock schedule with disk tumbler locks. All locks on the project are keyed different with a single master key that can open all locks. Either disk tumbler or pin tumbler locks could be used for this arrangement.

LOCK SCHEDULE

Project Na	Project Name: SAMPLE PROJECT							Project/Quote Number Date:			or Order Dept. use
Lock Type	Disk Tu Pin Tur		Number of Keys per Lock	2	2 Number of Master Keys per change 4 Number of Grand-Master Keys per change						
Room Number	Item Number	Cabinet (all cabin	nets with locks mus	t be listed)		Number of Locks	Keyed Different	Keyed Key Alike Group	Masterkey Group	Grand- Masterke Group	у
100	1	D30C	3518-0101-J1			4	V		А		
100	1	D36C	3536-0101-J	1		6	V		А		
100	2	E40C	3524-0101-J1			2	V		А		
105	6	D30C	3524-0101-J1			4	V		А		
105	6	E41C3	5536-0101-J1			3	V		А		
105	6	W20C	2936-0101-J	1		1	V		А		

Project/Quote Number

Keyed Alike in Sets

An example of another simple lock schedule with disk tumbler locks. On this project all the locks in each room are keyed alike, but each room is different from the other. A master key could be used in this type of keying arrangement, but in this example it is not. Either disk tumbler or pin tumbler locks could be used for this arrangement.

LOCK SCHEDULE

SAMPLE PROJECT

Sheet of

for Order Dept. use

Lock Type	Disk Tu Pin Tur		Number of Keys per Lock	2	r of Master er change	r		er of Gra er chan	and-Master ge		
Room Number	Item Number	Cabinet (all cabir	nets with locks must	t be listed	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand- Masterkey Group	
115	5	D30C	3518-0101-K1		4		1	Α			
115	5	E56C	3536-0101-K1		4		V	Α			
115	5	E41C3	5536-0101-K1		3		V	Α			
126	16	E40C	3524-0101-K1		2		1	В			
126	16	E42C	3536-0101-K1		5		✓	В			
126	17	CO3C	3536-0101-K1		1		✓	В			
126	17	D30C	3524-0101-K1		4		V	В			
							1				1

Keyed Different - Masterkeyed & Grand-Masterkeyed

An example of a more complex lock schedule with pin tumbler locks. All locks are keyed different, with a separate master key for each room, and a grand-master key for the entire project. This keying arrangement could not be used with disk tumbler locks.

LOCK SCHEDULE

SAMPLE PROJECT

Sheet of

for Order Dept. use

Lock ☐ Type ▼	Disk Tu Pin Tur		Number of Keys per Lock	2		er of Maste oer change		Number of Gra Keys per chan		5	
Room Number	Item Number	Cabinet (all cabir	nets with locks mus	t be listed))	Number of Locks	Keyed Different	Keyed Key Alike Group	Masterkey Group	Grand- Masterkey Group	
102	1	E43C	3536-0109-J	12		6	V		Α	Α	
102	1	E40C	3518-0109-J	2		2	V		Α	Α	
102	1	D30C	3524-0109-J	12		4	V		Α	Α	
105	3	D40C	3524-0109-J	2		2	V		В	Α	
105	3	D30C	3524-0109-J	12		4	V		В	Α	
1113	9	D32C	3536-0109-J	12		8	V		С	Α	
1113	9	CO3C	3536-0109-J	J2		1	V		С	Α	
1113	10	D30C	3524-0109-J	12		4	V		С	Α	
								1			

Project/Quote Number

LOCK SCHEDULE

Sheet	of

Project Na	me:					Project/Q	uote Number:		Date:		for Order Dept. use
Lock Type	Disk Tu Pin Tu	umbler mbler	Number of Keys per Lock		Numbe Keys pe	I r of Maste er change	r	Number of G Keys per ch	irand-Master ange		
Room Number	Item Number	Cabinet (all cabin	nets with locks mus	st be listed		Number of Locks	Keyed Different	Keyed Key Alike Grou	Masterkey Group	Grand Masterk Group	ey



Alphabetical Index

A		Drench Shower
	Acid Storage Cabinet-ADA Height	E
	Acid Storage Cabinet-Standing Height 27 & 32	Electrical Fixtures
	Acid Storage Fume Hood Cabinet	Emergency Eye Wash/Shower90
	Acid Storage Cabinet - Full Height	End Leg Assembly
	ADA Height Cabinets	End View Dimensions
	ADA Compliant Sink	Epoxy Resin Cupsink91
	Adjustable Slip Scribes	Epoxy Resin Drain Trough92
	Aerator87	Epoxy Resin Top 70-71
	Air Valves83-84	Epoxy Resin Tub Sink 92-93
	Applique Panels	Eye Wash90
	Apron Rails	
	Aprons	F
	Aspirator	Face Wash90
		Faucet Accessories
В		Filler Panel
	Ball Valves	Filter Pump87
	Blind Corner Apron	Finished Back
	Blind Corner Cabinet	Fixtures-Electrical
	Balance Table	Foot Valve87
	Base Molding	Freezer
	2400	Full Height Storage Cabinets 53-57
C		Full Overlay Applique Panels
	Cabinet Dimensions 5 & 100	Fume Hood Acid Storage Cabinet–ADA Ht 20
	Cabinet Styles	Fume Hood Acid Storage Cabinet–Standing Ht31
	Cart - Stainless Steel	Fume Hood Base Cabinet-ADA Height 20
	Caster Options6	Fume Hood Base Cabinet–Standing Height 30
	Catalog Number Explanation	Fume Hoods
	Center Bench Shelving	1 diffe 1100d3 00-00
	Ceiling Enclosure	G
	Chalkboard79	Gas Valves
	Cold Water Faucet	Gooseneck
	Combination Faucet84	Gooseneck Coupling
		Gooseneck Coupling
	Construction Details	••
	Corner Box Apron	H
	Corner Cabinets–ADA Height	Hoods
	Corner Cabinets–Sitting Height	Hose Spray
	Corner Cabinets–Standing Height29	Hot & Cold Water Mixing Faucet 86-87
	Corner Cabinets–Wall Mounted	
	Corner Fillers	1
	Corner Guard81	Intermediate Leg Assembly44
	Crumb Cup Strainer91	
	Cupboard/Drawer Cabinets-ADA Height 16-17	K
	Cupboard/Drawer Cabinets-Sitting Height 12-13	Kemresin Work Tops 70-71
	Cupboard/Drawer Cabinets-Standing Height. 24-25	Kemresin Trough
	Cupboard Cabinets-ADA Height14	Kemresin Tub Šink92-94
	Cupboard Cabinets-Sitting Height	Kemstruts80
	Cupboard Cabinets-Standing Height	Key Case
	Cupsink	Keying Options
	·	Knee Space Panel
D		·
	Deionized Water Faucet	L
	Diagonal Corner Cabinets-ADA Height 18	Label Holder
	Diagonal Corner Cabinets–Sitting Height 13	Laboratory Cart
	Diagonal Corner Cabinets—Standing Height 29	Laboratory Service Fittings 83-90
	Distilled Water Faucet	LED Lighting
	Double-Sided Wall Cabinet	Leg Assemblies
	Drain Outlet	Lock Keying
		Lock Panels
	Drain Trough	Lock Schedule Examples
	Drawer Cabinets—Sitting Height	Lock Schedule Sample
	Drawer Cabinets-Standing Height 22-23	Locks8

Alphabetical Index (continued)

M		Stainless Steel Sink
	Markerboard79	Stainless Steel Sink Outlet
	Mobile Cabinets	Standing Height Base Cabinets
	Mixing Faucet	Steam Valves
	Multiple Service Fittings84	Stools
		Styles
N		Support Angle
	Needle Point Cock Valves 83-84	Support Strut
	Needle Point Cock valves	Support Strut
_		Swinging Glazed Door Full Height Cabinet 54
0		Swinging Glazed Door Wall Cabinet 50
	Open Cabinets-Sitting Height	Swinging Panel Door Full Height Cabinet 54
	Open Cabinets-ADA Height	Swinging Panel Door Wall Cabinet 51
	Open Cabinets-Standing Height22	Switch
	Open Cabinets–Full Height	
	Open Cabinets-Wall Mounted46	т
	Open Cabinets-Wall Mounted40	Tackboard
_		
P		Task Lighting
	Part Number Explanation4	Tops
	Pass-Thru Wall Cabinet	Trap
	Pedestal Legs	Trough
	Pegboard78	Tub Šink92-93
	Pipe Enclosure67	
	Plastic Laminate Work Tops 70-71	U
	Plastic Laminate Lodges & Chalves 70	
	Plastic Laminate Ledges & Shelves	Ultra-Pure Water Faucet
	Plumbing Information82	Under Counter Refrigerator
	Portable Pump Faucet	Upright Rod Assembly81
R		V
	Reagent Rack Assemblies	Vacuum Breaker
	Reagent Rack Ledges76	Vacumm Pump Storage Fume Hood Cabinet 20&33
	Rear Scribes & Filler	Vacumm Pump Cart
	Receptacle - Electric	Vacuum Valves
	Defricantes	
	Refrigerator	Vandal-Resistant Fittings
	Rods & Accessories81	Ventilated Enclosures 68-69
_		
S		W
	Safety Shower	Wall Cabinets
	Scribe	Wall Shelving
	Seating	Wardrobe Storage Cabinet
	Service Fittings 82-90	Waste Storage Cabinet
	Shelf Brackets	Water Faucet
	Shelf Lips	Water Mixing Faucet
	Shelving Assemblies	Water Valves
	Sink	Water/Gas Fitting
	Sink Cabinets-ADA Height	Work Tops
	Sink Cabinets-Sitting Height	Work Top Support Frame43
	Sink Cabinets-Standing Height26	Wrist Blade Handle87
	Sink Outlet	
	Sink Overflow	
	Sink Trap	
	Sitting Height Base Cabinets	
	Sloping Tops	
	Sliding Glazed Door Full Height Cabinet53	
	Sliding Glazed Door Wall Cabinet 48	
	Sliding Panel Door Full Height Cabinet 53	
	Sliding Panel Door Wall Cabinet 48	
	Sliding Plate Glass Door Wall Cabinet	
	Solvent Storage Cabinet 19-20, 28 & 32	
	Spreader	
	Spring Clips	
	Stainless Steel Cart	
	Stainless Steel Cupsink	



Numerical Index

Part Number	Page	Part Number	Page
019623 -	019627	C03C292215L -	C03C292248 11
010020	010027	C03C322215L	
0464-00 -	91	C03C352212L	
0482-BP -	0482-GP 91	CCEC121320-00_0 -	CCEC241615-00_067
0490-BE -	0490-SE91	CSTC081327-00 0 -	_
0491-BP -	0491-GP 91		00.00=: 00_00
0492-BP -	0492-GP 91	D00C292215	D00C29224811
0493-BE -	0493-SE91	D00C322215	D00C322248 15
0494BP-04 -	0494BP-10 91	D00C352215	D00C35224822
0548-00 -	0557-00	D02C322230	
		D02C352230	D02C35224822
0720-1K-6 -	0720-1R-6a 97	D08C352215	D08C352248 22
		D26C292218	
0950-0F -	0975-0F	D30C292215	D30C292224 11
0995-00 -		D30C322215	D30C322248 15
		D30C352212	D30C35224823
1000-00	1039-00 93	D31C292230	D31C292248 11
1000-DI	1070-DI 94	D31C322230	D31C322248 15
1000-ADA-DI	1006-ADA-DI 94	D31C352230	D31C35224823
		D32C292230	D32C29224811
4740-00	81	D32C322230	D32C322248 15
		D32C352230	
A00C022215-00_0 -	7.000022272 00_0 111111 10	D36C322230	
A00C0422150 -	7.0000 IZE7E000	D36C352230	D36C35224823
A00C0428150 -	7.0000 12072000	D50C322215	D50C322224 15
A04C042248-1	, 10000 1207 E 1 12	D50C352215	
A06C020215-00_0 -	7.0000 10272 00_0 : : : : 10	D51C322230	D51C322248 15
A06C040315-10 -	A06C040372-10 43	D51C352230	20.00022.0
A07C292202-00_0 -	A08C352802-00_0 44	D52C322230	
A09C020230-00_0 -	7.000020272 00_0 10	D52C352230	D52C35224823
A15C290202-00_0 -	A15C350202-00_0 45	D81C292218	D81C292224 12
A10C250015-00_0 -	A10C350072-00_0 45	D81C352218	
A11C292202-00_0 -	A13C355202-00_0 45	D83C322215	D83C322224 15
A27C292202-00_0 -	A28C352802-00_0 44	D84C292218	
A51C042220-1	A53C042272-1 40	D84C322215	
A51C042820-1	A53C042872-1 41	D84C352218	D01000222111111120
A56C040236-00_0 -	A56C040272-00_0 43	D85C292218	D85C292224 12
A57C352202-00_0 -	7.000002002 00_0 : : : : : : : :	F0F0000001	F05000010
A65C350202-00_0		E05C292230L	
A77C352202-00_0 -	A78C352802-00_0 44	E05C322230L	
A DCC 000010 00 0	AD00050040.00.0	E05C352230L F40C292215I -	
	ABSC350012-00_0 65	L 1002022 10L	L 10020222 1 12
	AWSC481612-00_065	E40C322215L	E40C32224816
AFSC841312-00_0 -	AFSC842212-00_0 65	E40C352212L	
B00C00001E 00 0	B00C000048 00 0 11	E41C292230 E41C322230	E41C292248 12 E41C322248 16
	B00C292248-00_0 11		
B00C322215-00_0 - B00C352215-00_0 -		E41C352230 E42C322230L	
B00C332213-00_0 -	D00C332240-00_0 22	E42C352230L	
BCAC200002 11 0	BCAC350012-17_059	E43C322230L	
BCFC290003-11_0 -	——————————————————————————————————————	E43C352230L	E43C35224824
BEAC290003-11_0 -			E45C292248 13
BESC290001-00_0 -			E45C32224816
BFBC290112-00 0 -		E45C352230L	
BIAC290006-11_0 -			E56C352248 16
BIFC290006-11_0 -		E56C352236	
BPCC290001-00- 0 -		E60C352212L	E60C352248 25
_	BRSC350024-00_0 59	E61C352230	
BPAC290023-16 0 -		E62C322230L	
	222222 12_3 00		E62C352248

Numerical Index (continued)

Part Number	Page	Part Number	Page
E63C322230L	- E63C32224817	G35C322024L	G35C322248 20
E63C352230L	- E63C35224825	G35C352024L	
E65C322230L	- E65C32224817		
E65C352230L	- E65C35224825	G63C322018L_100 -	G63C322248100 19
		G63C352024L_100 -	
ECEC120013-00_0	- ECEC240022-00_0 67	G68C322024L_100 -	G68C322248100 20
		G68C352024L_100 -	
F-0260-00		G80C322018L	
F-0862-04		G80C352018L	G80C352260 31
F-1368-00			
F-2939-00	- F-2940-00 78	H581-1K -	H581-8R 96
F-4940-00	- F-4945-00 98		
F-4946-00	- F-5352-00 99	K-0005-0A -	K-0020-1A
F-4995-PK		K7-302435-BK -	K7-302435-SL79
F-5332-00	- F-5332-01 79	K12-2905-0A -	K12-3520-0A80
F-5333-00	- F-5348-00	K17-3612-00 -	K17-3614-00
F-5349-00		KF-0170-04	
F-5353-06	- F-5362-12 95	K-MSS-12 -	K-MSS-60 13, 19 & 26
F-8440-15	- F-8441-48 27		
	F. (1) (1) (1)	L-105 -	L-115
F10M0242-00	- F10M0279-00	L-2034-06 -	L-2034-96 80
F11M0196-00	- F11M0372-00	L01-020109-00 -	L01-020142-00 47 & 81
F12M000224-00	- F13S000284-0072		
F17-4824-00	- F17-4872-00	M-0051-0A -	M0052-0A 80
F21M0624-00	- F21W2448-0	M-0030-0A -	M-0044-AD
F21M0924-60	- F21W1848-6	M-0060-0A -	M-0062-AD 77
F63M0224-00	- F63M0272-00	M400000451	N40000004
F73M0015-00	- F73S0024-0037	M10C262215L	
F94M2215-00	- F94M2348-00	M10C332215L	
F94P2215-00	- F9412448-00	M20C262215L	
FCAC840003-11 0	- FCAC840012-16_0 63	M20C332215L M25C262215L	
FCEC120018-00 0	- FCEC240096-00_0 67	M25C332215L	M25C332224 34
FCFC841301-00_0	- FRSC840024-00_0 63	M30C262215	
FEAC840003-11 0	- FEAC840024-16_0 62		
FESC841301-11_0	- FESC842224-00_0 62	M35C262215 M40C332215	M40C332224 35
FFBC840118-00 0	- FFBC840148-00_0 64	M45C262215	
FIAC840006-11 0	- FIAC840024-16 0 62	M46C332215	M46C332224 35
FIFC841301-00_0	- FIFC842224-00_0 62	M56C332215	M56C332224 35
FPAC840013-16_0	- FPAC840022-16_0 63	M62C262215	
FRSC560001-00 0	- FRSC840024-00_0 63	M66C262215	
		M68C332215	M68C332224 35
G00C292215L	- G00C292248 13	M71C332215	
G00C322212L	- G00C322248 19		
G00C352212L	- G00C352260 26	PL-4812-00 -	PL-7214-00 76
G03C322236	- G03C322248 21		
G03C352230	- G03C352248 26	RRU-1210-0000 -	RRU-1812-0030 76
G04C352252	- G04C352260 26		
G07C322236	- G07C322248 21	S00C841318-00-0 -	S00C842248-00_0 53
G08C322012L	- G08C322248 20	S12C8413300 -	S12C842248053
G08C352012	- G08C352260 30	S25C841318L	S25C84224854
G12C322248	- G13C322248 18	S35C841330	00000122101111100
G12C352248	- G13C352248 29	S40C842218L	
G14C292214		S41C8413300 -	0 1 1 0 0 1 1 2 1 0 0 1 1 1 1 1 1 0 0
	- G15C352214 18	S50C841318L	
G14C352214		S51C841624L	
G20C352215	_ G20C352224 29	S60C841330	
G23C352218		S70C841624L	S70C84224857
G21C322018L	- G21C322248 18		
G21C352015L	- G21C352248 27		
G23C322218	19		





Numerical Index (continued)

Part Number \$581-1K \$582-1K \$655-1K \$656-1K \$657-6K \$658-6K \$681-1K \$721-6K-6 \$751-6K-6		Page \$581-8R 96 \$582-8R 96 \$655-1R 97 \$656-3R 96 \$657-8R 96 \$658-8R 96 \$681-8R 96 \$721-6R-6a 97 \$751-6R-6a 97
SCAC010103-00_0 SCFC081301-00_0 SLTC081315-00_0 SM-01-00 STFC081301-00_0	- - - -	SCAC010196-00_0 65 SCFC142206-00_0 66 SLTC142296-00_0 66 SM-56-00- 93 STFC142212-00_0 66
TR-01-00	-	TR-05-00 92
U01C600512-00_0 VPCC061414-00_0	-	U02C961048-00_0 67 VPCC061644-00_0 33
W-0242-00 W-0259-00 W-0260-00 W-0273-02 W-0286-02 W-0300-00 W-0337-00 W-0340-00 W-0350-00 W-0351-00 W-0354-00 W-0354-00 W-0361-00 W-0367-00 W-0374-01 W-0400-00 W-0402-00 W-0402-00 W-0405-00 W-0410-VR W-0510-00 W-0520-00 W-0532-00 W-0535-00 W-0535-00 W-0833-00 W-0923-00 W-0928-00 W-0928-00 W-0931-00		W-0243-00. 83
W-0932-00 W-0933-00 W-0934-00 W-0935-00 W-0936-00	- - -	

Part Number			Page
W00C241315-00_0 -	-	W00C481648-00_0	46
W10C241330-00_0 -	-	W10C481648-00_0	47
W12C2413300 -	-	W12C4816480	
W13C241330-00_0 -	-	W13C361648-00-0	49
W25C241315	-	W25C481648	
W26C241320	-	W26C481615	52
W41C2413300 -	-	W41C4816480	
W50C241315	-	W50C481648	50
W51C241320	-	W51C481615	
WIAC240006-11_0	- - - -	WCAC480012-16_0 WCFC481612-00_0 WEAC470024-16_0 WESC481624-00_0 WFBC480148-00_0 WIAC480048-16_0 WIFC481648-00_0 WPAC480016-16_0	61 60 60 64 60
X-020004 - X-030001-00 -	-	X-020017 X-030003-00	78 80



Typical Installations









Trademarks

Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation.

Part numbers, illustrations, descriptions, specifications and other information contained in this catalog are based on the most current design and production information available

Patents:

Products in this catalog are protected by one or more of the following patents:

4,856,420 5,241,788 5,302,030 5,407,389 5,582,472 5,728,626 D388,986 DES. 346,437 Additional patents are pending.

at the time of publication. Kewaunee Scientific Corporation reserves the right to make changes at any time without notice.



Typical Installations









We are pleased to have been selected to supply laboratory furniture for the facilities of some of the best companies, organizations and institutions in the world. We believe the repeat opportunities, over time, to participate in their construction projects is good evidence of our product quality and experience in successfully completing their laboratory projects.

Whether your laboratory projects are large or small, new construction or remodeling, Kewaunee has the qualified people and quality products to transform your visions into successful reality.

From pre-planning to final completion, experienced Kewaunee representatives are available to help you through the intricacies of laboratory design and layout, product selection, ordering, scheduling and installation. Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer Total Laboratory Solutions, Innovative Designs, and Technological Expertise with a worldwide distribution network to provide laboratories that are truly World Class.

Laboratory Products from Kewaunee®:

ALPHA® – Adaptable Casework System

ALPHA Columns – Expandable Workstations

ALPHA Overhead Carrier Systems

BasikBench® – Welded C-Leg Bench

DetectAir – Filtered/Ductless Fume Hoods

EarthLine® – Environmentally Friendly Casework

Enterprise® – Moveable Workstations

Evolution® – Column Based Workstations

Interceptor® – Biological Safety Cabinet

Interceptor® – Laminar Flow Cabinet

Kemresin® – Epoxy Resin Counters & Sinks

Kewaunee Matrix® – Educational Stations

Research Collection® - Steel Casework

Signature® Series – Wood Casework

Spektrum – Laminate Architectural Millwork

Sturdilite® – Flexible Pedestal System

Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to www.kewaunee.com, or contact your local sales representative directly.











Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation



P.O. Box 1842 • Statesville, NC 28687-1842 Phone: (704) 873-7202 • Fax: (704) 873-5160

Kewaunee Labway India Pvt. Ltd.

No. CA-9A, 2nd Floor, Jigani Link Road Jigani, Bangalore 562 106, India Phone: 9180-27826725 • Fax: 9180-27826724 Web Site: www.kewaunee.in

Kewaunee Labway Asia Pte. Ltd.

194 Pandan Loop, #06-22 Pantech Business Hub, Singapore 128 383 Phone: 65-6773-0288 • Fax: 65-6773-2322

Kewaunee Scientific (Suzhou) Co Ltd

333 Xingpu Road, Modern Industrial Square Building 4, Unit #2-A • Suzhou Industrial Park • China Phone: (0512) 6262 5090 • Fax: (0512) 6262 5089 Web Site: china.kewaunee.com

Kewaunee Steel Cabinet Construction	2
	3
Kewaunee Steel Cabinet Styles	
Steel Cabinet Information	4
Steel Cabinet Options	6
Lock Keying Options	9
, , ,	
Steel Sitting Height Base Cabinets	11
Steel ADA Height Base Cabinets	14
Steel Standing Height Base Cabinets	22
Fume Hood Base Cabinets	30
Steel Sitting Height Mobile Cabinets	34
Steel Standing Height Mobile Cabinets	35
Mobile Cabinet Tops & Accessories	36
Steel Apron Assemblies	38
22" Deep Apron	38
22" Deep Apron	38
Steel Box Apron	38
22" Deep Apron	38
• •	
Applique Panel	38
Steel Box Apron	38
28" Deep Apron	39
	39
28" Deep Apron	
Steel Box Apron	39
28" Deep Apron	39
Applique Panel	39
	39
Steel Box Apron	
Steel Leg Assemblies & Pedestals	44
Steel Apron Rails & Panels	45
Steel Wall Mounted Storage Cabinets	46
_	53
Steel Full Height Storage Cabinets	
Base Cabinet Scribes & Fillers	58
Wall Cabinet Scribes & Fillers	60
Full Height Cabinet Scribes & Fillers	62
•	
Steel Cabinet Finished Backs	64
Steel Adjustable Slip Scribes	65
Steel Sloping Tops	66
Steel Ceiling & Pipe Enclosures	67
• .	
Fume Hoods & Ventilated Enclosures	68
Work Tops	70
Center Bench Shelving Assemblies	72
Wall Shelving Assemblies	74
Reagent Rack & Shelving Systems	76
Pegboards	78
Miscellaneous Accessories	79
Miscellaneous Accessories	80
K12-	80
-0A	80
General Plumbing Information	82
•	
Laboratory Service Fittings	83
Cupsinks & Drain Accessories	91
Kemresin Tub Sinks & Drain Troughs	92
Kemresin Drop-in Sinks	94
·	95
Stainless Steel Drop-in Sinks	
Laboratory Electrical Fixtures	96
Laboratory Seating	98
Steel Cabinet Technical Information	100
Lock Schedule Sample and Examples	104
Alphabetical Index	106
Numerical Index	108
Typical Installations	112
Typical Indianations	114